

Instruction for Use

NOVASTREAK Microbial Contamination Monitoring Device (Cat. No. BD-508) Baird Parker Agar and Violet Red Bile Agar (VRBA)



INTENDED USE

NOVASTREAK MCMD BD-508 is a convenient semi-quantitative screening culture device for sampling and assessing microbial contamination of food and dairy products, industrial fluids and surfaces of sanitary importance. A unique streaking mechanism permits the isolation of single colonies even when the original bacterial population of the sample was as high as 10^7 organisms per milliliter. NOVASTREAK MCMD BD-508 is intended for use in the food industry.

SUMMARY AND EXPLANATION

NOVASTREAK MCMD BD-508 comprises a plastic paddle with two types of agar attached back-to-back, housed in a closed transparent plastic tube. A ring with elongated prongs is attached to the end of the paddle so that there are prongs on each side of the slide. The ends of the prongs are dipped into the liquid sample. Upon re-insertion into the plastic tube, the prongs are prevented from moving and the agar surfaces are inoculated with the sample as the paddle passes over the prongs. The result is a series of streaks of decreasing bacterial concentration, which permits isolation of single colonies even when the original bacterial population of the sample was as high as 10^7 organisms per milliliter. NOVASTREAK MCMD BD-508 can be used to monitor microbial growth wherever the potential may exceed 10^2 microorganisms in ml of sample. NOVASTREAK MCMD BD-508 unit consists of two different agar modifications: Baird Parker Agar and Violet Red Bile Agar (VRBA), attached back-to-back on a plastic sampling paddle, which is permanently fastened to the

cap for comfort of handling during use.

SAMPLING

| SAMPLE TYPE | MATERIAL TO BE TESTED | PROCEDURE | READING |
|-----------------------------|--|--|--|
| Liquid samples | Milk (raw and pasteurized) | Dip sampling procedure | Compare with Colony Density Chart No. 2 |
| | Industrial water (waste, recycled, cooling or process water), dairy products (starter cultures, sour cream, yogurt and other fermented products) | General streaking procedure | Compare with Colony Density Chart No. 1 |
| Viscous and friable samples | Syrups, pasts and dehydrated products (vegetables, fruit, egg powder, milk powder, powdered soups, instant desserts, cocoa, etc.) | Dilute 1:1 or 1:10 in sterile water with following General streaking procedure | Compare with Colony Density Chart No. 1 and multiply the result by 2 or 10 |
| Solid samples | Raw material, frozen and chilled products (meat, fish and sea food products) | Homogenize and suspended 1:1 or 1:10 in sterile water with following General streaking procedure | Compare with Colony Density Chart No. 1 and multiply the result by 2 or 10 |
| Surfaces | Utensils, work surfaces | Touch surface with paddle for several seconds | Compare with Colony Density Chart No. 3 |

GENERAL PROCEDURE

A. STREAKING SAMPLING

1. Unscrew the NOVASTREAK MCMD cap. Pull the paddle out. Do not touch any part but the cap.
2. Hold the paddle vertically and dip the white prongs into the sample up to about half of their length (see below).
3. Return the paddle to its container in a quick, continuous and vertical motion and tighten cap.
4. Transport the tube to laboratory for incubation and examination
5. Before incubation, loosen cap one-half turn.
6. Incubate the entire container at $(35^{\circ}\text{C} \pm 2^{\circ}\text{C})$ for 18-24 hours in a vertical position.
7. Interpret the results by simple visual comparison of bacterial growth on the agar surface with Colony Density Chart No. 1 provided. **No actual colony counting is necessary.**

B. DIP SAMPLING (DIPSLIDE TECHNIQUE)

1. Unscrew the NOVASTREAK MCMD cap. Pull the paddle out. Do not touch any part but the cap.
2. Immediately return the paddle to the tube (in order to move the prongs out of the way) and then pull the paddle out again.
3. Dip the culture paddle into a diluted/undiluted sample or pour the sample over agar surfaces, if the volume of sample is not adequate to fully immerse the agar surfaces.

4. Replace inoculated culture paddle in its protective NOVASTREAK MCMD vial and close cap.
5. Transport NOVASTREAK MCMD vial to laboratory for incubation and examination.
6. Place inoculated NOVASTREAK MCMD vial upright in incubator (35°C±2°C) for 18-24 hours. Before incubation, loosen cap one-half turn.
7. Interpret the results by simple visual comparison of bacterial growth on the agar surface with Colony Density Chart No. 2 provided. **No actual colony counting is necessary.**

C. SURFACE CONTACT SAMPLING

1. Unscrew the NOVASTREAK MCMD cap. Pull the paddle out. Do not touch any part but the cap.
2. Immediately return the paddle to the tube (in order to move the prongs out of the way) and then pull the paddle out again.
3. Gently touch the agar faces onto the surface to be tested. The agar should remain in contact for about 20 seconds.
4. Replace inoculated culture paddle in its protective NOVASTREAK MCMD vial and close cap.
5. Transport NOVASTREAK MCMD vial to laboratory for incubation and examination.
6. Place inoculated NOVASTREAK MCMD vial upright in incubator (35°C±2°C) for 18-24 hours. Before incubation, loosen cap one-half turn.
7. Interpret the results by simple visual comparison of bacterial growth on the agar surface with Colony Density Chart No. 3 provided. **No actual colony counting is necessary.**

MATERIALS PROVIDED

| CAT. No | REAGENTS | COLOR OF REAGENTS | EXPECTED TYPES OF MICROORGANISMS |
|---------|--|--|--|
| BD-507 | Side 1: Baird Parker Agar Side 2: Violet Red Bile Agar (VRBA) | Agar Color: Yellow, Opaque Agar Color: Purple | Staphylococci Total Coliforms Count |

CLASSICAL COMPOSITION (g/liter)

1. **Baird Parker Agar:** Pancreatic digest of Casein 10; Beef extract 5; Yeast extract 1; Lithium chloride 5; Agar 17; Glycine 12; Sodium pyruvate 10.
2. **Violet Red Bile Agar (VRBA):** Pancreatic digest of Gelatin 7; Yeast extract 3; Bile salts mixture 1.5; Lactose 10; Sodium chloride 5; Agar 15; Neutral red 0.03; Crystal violet 0.002; MUG (4-methylumbelliferyl-beta-D-glucuronide) 0.1.
- 3.

APPLICATION FIELDS

| MATERIAL TO BE TESTED | TYPE OF MICROORGANISMS | | |
|--|------------------------|-----------------------|---------------|
| | Total Coliforms Count | Total Bacterial Count | Yeast & Molds |
| Water (waste, recycled, cooling or process water) | • | • | • |
| Raw Milk | • | • | |
| Dairy Products (pasteurized milk, starter cultures, sour cream, yogurts and other fermented dairy products) | • | • | • |
| Meat, fish, sea food (raw material, frozen and chilled products) | • | • | |
| Surfaces | • | • | • |
| Syrups, pasts and dehydrated products (vegetables, fruit, egg powder, milk powder, powdered soups, instant desserts, cocoa, etc.) | • | • | • |

MATERIAL REQUIRED BUT NOT PROVIDED

Incubator (35 ± 2°C)

Incubation Stand

WARNING AND PRECAUTIONS

1. For *In Vitro* Diagnostic Use.
2. Use aseptic technique and established laboratory procedure in handling and disposing of infectious material.

STORAGE

1. Store NOVASTREAK MCMD at 2-8°C up to 6 months, refer to product label.
2. Protect contents from direct light to ensure product stability through the expiration date, shown on the tube cap or packaging label.

KIT CONTENTS

| REAGENTS | EXPECTED RESULTS |
|---|---|
| Baird Parker Agar is used for the selective isolation and enumeration of coagulase-positive staphylococci from food, skin, soil, air and other materials. It may also be used to identify staphylococci on the basis of their ability to clear egg yolk (lecithinase production). Sodium pyruvate is incorporated to impart <i>S. aureus</i> their black colony appearance. Glycine and lithium chloride have inhibitory action on organisms other than <i>S. aureus</i> . | Typical colonies of <i>S. aureus</i> are black, shiny, convex and surrounded by clear halos of approximately 2 to 5 mm. Coagulase-negative staphylococci do not grow well; if at all, the typical clear halos are absent. Other organisms are inhibited or may grow sparsely, producing white to brown colonies with no clearing of the egg yolk. |
| Violet Red Bile Agar (VRBA) is a selective medium for the detection of coliform organisms in water, milk, and other materials of sanitary importance. The medium is selective due to the presence of inhibitors, bile salts and crystal violet. Differentiation of enteric microorganisms is achieved by the combination of lactose and the neutral red indicator. Colorless or pink-to-red colonies are produced depending on the ability to ferment lactose. | Lactose fermenters are rose-red in color and generally surrounded by a halo of precipitated bile. <i>E. coli</i> colonies are entire edged, 1mm or more in diameter. <i>E. aerogenes</i> are larger, often mucoid and pinkish. Lactose non-fermenters produce colorless colonies. Enterococci occasionally grow to produce rose colonies pinpoint in size |

DISPOSAL

The used *NOVASTREAK* MCMD is disposed by standard methods of biohazard disposal.

EXPIRATION DATE

1. The expiration date applies to the product in its intact container when stored as directed.
2. Do not use *NOVASTREAK* MCMD exhibiting any of the following characteristics: discoloration, dehydration, wrinkling or shrinkage of an agar surface; microbial growth prior to inoculation; or an atypical cultural response in Quality Control procedures.


INCUBATION CONDITIONS

| REAGENTS | OBJECT | TEMPERATURE (°C) | PRELIMINARY RESULTS (Hrs) | FINAL RESULTS (Hrs) |
|--------------------------|-----------------------|------------------|---------------------------|---------------------|
| Baird Parker Agar | Staphylococci | 35-37 | 24 | 48 |
| VRBA | Total Coliforms Count | 35-37 | 18 | 24 |


INTERPRETATION OF RESULTS

1. BACTERIAL COUNT


If more than 200 colonies grow on the *NOVASTREAK* MCMD, the growth may become semi-confluent and the presence of more than 100,000 bacteria per ml is indicated. If fewer than 20 colonies are counted, less than 10,000 bacteria per ml is indicated.




COLONY DENSITY CHART No. 1 FOR SAMPLING BY STREAKING




COLONY DENSITY CHART No. 2 FOR TOTAL GROWTH BY DIP SAMPLING




COLONY DENSITY CHART No. 3 FOR TOTAL GROWTH BY SURFACE CONTACT SAMPLING



COLONY DENSITY CHART OF YEAST GROWTH BY SURFACE CONTACT SAMPLING



COLONY DENSITY CHART OF MOLD GROWTH BY SURFACE CONTACT SAMPLING



LIGHT CONTAMINATED (up to 1 colony)
 MODERATELY CONTAMINATED (up to 10 CFU)
 HEAVILY CONTAMINATED (confluent growth)

2. COLONIES MORPHOLOGY

Preliminary identification of the microorganisms made on the base of type and color of the colonies.

| ORGANISMS | Baird Parker (Agar Color: Yellow, Opaque) | VRBA (Agar Color: Purple) |
|-----------------------------------|---|---------------------------|
| <i>E. coli</i> , <i>Coliforms</i> | No growth | Red or pink |
| <i>S. typhimurium</i> | No growth | Colorless |
| <i>S. epidermidis</i> | Grey-Black colonies | No growth |
| <i>C. albicans</i> | No growth | No growth |

QUALITY CONTROL

Baird Parker Agar:

- Identity Specifications:** (1) sterility of media: as per sterility test of Standard Operating Procedures; (2) physical appearance: Yellow, opalescent; (3) pH of media: pH 7.0±0.2; (4) weight of media: 3.15±0.1g;
- Cultural Response:** (5) inoculate challenged media with the following microorganisms (as per inoculation procedure of Standard Operating Procedure):

| Microorganisms | ATCC no. | Growth | Appearance |
|-----------------------------------|----------|------------|--|
| <i>Escherichia coli</i> | 25922 | Inhibition | None |
| <i>Staphylococcus aureus</i> | 25923 | Growth | Black colonies surrounded by clear zones |
| <i>Staphylococcus epidermidis</i> | 12218 | Growth | Black colonies |

Violet Red Bile Agar (VRBA):

- Identity Specifications:** (1) sterility of media: as per sterility test of Standard Operating Procedures; (2) physical appearance: Reddish purple, slightly opalescent; (3) pH of media: pH 7.4±0.2; (4) weight of media: 3.15±0.1g;
- Cultural Response:** (5) inoculate challenged media with the following microorganisms (as per inoculation procedure of Standard Operating Procedure):

| Microorganisms | ATCC no. | Growth | Appearance |
|-------------------------------|----------|------------|----------------------|
| <i>Salmonella typhimurium</i> | 14028 | Growth | Colorless colonies |
| <i>Escherichia coli</i> | 25922 | Growth | Rose to red colonies |
| <i>Staphylococcus aureus</i> | 25923 | Inhibition | None |
| <i>Shigella flexneri</i> | 12022 | Growth | Colorless colonies |