



# FALL MANAGEMENT

## SYSTEMS AND ACCESSORIES



The PM100 Personal Monitor features a dome and string mechanism that is attached to the patient's garment. When an unassisted exit takes place, the dome is pulled from its magnetic socket and an alarm sounds. The alarm is silenced when the dome is placed back into the socket. The Personal Monitor can be applied quickly, and moved easily from one location to another.

## PM100 PERSONAL MONITOR

### EASE OF USE

- "Alarm", "Low Battery", and "In Use" indicator LEDs show status at a glance

### FLEXIBILITY

- Nurse call jack allows direct notification to nurse call
- Bed, chair, and wall mount options are available

### SAFETY

- Magnetic dome allows activation when pulled from any direction
- Volume control located within the battery compartment deters patients from changing the volume level

### DURABILITY

- Durable thick plastic casing and impact-resistant rubberized edging reduces replacement/repair costs
- Two Year Warranty



# MOUNTING OPTIONS AND ACCESSORIES




All Curbell monitors feature a clip on the rear of the monitor that allows them to be mounted onto a wheel chair, or onto any of the optional mounting brackets listed below.







The clip slides easily into the bed mounting bracket (pictured above), Magnetic Mounting Bracket, or Wall Mounting Bracket, so you can mount the monitor virtually anywhere.




A Velcro loop is also supplied with every monitor. Just slip the loop through the bracket and hang the monitor on a wheelchair or bed rail.

Model	Description	Model Number
 Bed Mounting Bracket	Adjustable to fit 7/8" to 2 5/8" width headboards and footboards	CSM-BMB

# REPLACEMENT PARTS AND ACCESSORIES

Model	Description	Model Number
 Velcro Strap	For all monitors	CSM-VS
 Pull String/Magnetic Dome Replacement with Linen Clip	For the CSM-PM100 Personal Monitor	CSM-PS-L40
 Y Adaptor for Nurse Call Bell System	The doubling adaptor allows two call cords or fall management monitors to function off a single 1/4" patient station receptacle.	JU96-195
 Jumper Cable	This 8' jumper cable connects the monitor to the nurse call system	JU96-295

# PERSONAL MONITOR

Model	Description	Model Number
 PM100 Personal Monitor	Personal monitor - Pull string	CSM-PM100

