

URETHROtech®

i n n o v a t i o n f o r s a f e t y

Urethral Catheterisation Device (UCD)®

Integrated Guide Wire



“Empowering Non-Specialist Healthcare Professionals
to provide safe urethral catheterization”

SAFE URETHRAL CATHETERIZATION


Difficult male urethral catheterization is not an infrequent problem and injuries are expensive:

- 4% occur in cardiovascular surgery.
- 6% of ALL urological referrals over 1 year in one teaching hospital
- 0.7/1000 adult male hospital admissions
- 32% of men sustained significant urethral injury
- additional cost of \$371,790.00 in just two UK hospitals over a 6 month study period. [1-5]

Problem Solution

The UCD® is a simple yet effective catheterization device to solve difficult urethral catheterization. A non-traumatic hydrophilic Nitinol guide wire is integrated into a standard urethral 3-way catheter design providing a ready-to-use medical device.

The result is successful catheterization and uninterrupted patient care, avoiding the medical emergency of 'failed urethral catheterization'.



Safe Catheter Insertion

The non-traumatic guide wire forms one unit with the catheter and is inserted into the bladder first, bypassing enlarged prostate lobes.


Universally Suitable

The UCD® has a 5-year shelf life and is suitable for any clinical environment and for every healthcare professional who normally performs urethral catheterization. Avoidance of more invasive specialist procedures and hospital admissions.


1. Özkan B et al. J Neph Urol Res. 2014; 2: 32-36 2. Bacsu C et al. Can Urol Assoc J. 2013; 7:343-7 3. Thomas AZ et al. BJU Int. 2009; 104: 1109-1112 4. Davis NF et al. J Urol 2016; 196: 1473-77 5. Kashefi C et al. J Urol. 2008; 179: 2254-2257

SAFE

The UCD® device enables an unassisted 'One-man' bed side procedure. Insertion is easy and safe with the soft integrated guide wire. No trauma to urethra and prostate.



EASY



The guide wire passes into the bladder first and automatically drags the catheter behind as the guide wire curls up in the bladder.

PASS

Information

The UCD® is currently not marketed for sale in the US (awaiting FDA approval), but is used successfully in hospitals in the UK, Europe and Middle East.

We are passionate about the goal of eradicating Urethral Catheterization Injury (UCI). We would like to keep in touch and keep you up-to-date with the availability of the UCD® and the progress of FDA approval. We would welcome your involvement and evaluation.

Please get in touch, visit www.urethrotech.com

Product Code	Description	Package
UCD-0035-090-0016	Urethral Catheterisation Device (UCD®)	1/box

Distributed in Australia by:



Endotherapeutics

T: 1300 00 3636 F: 02 9869 7723
info@endotherapeutics.com.au
www.endotherapeutics.com.au

Manufactured by:

URETHROtech®



NICE accredited

www.nice.org.uk/accreditation