VITA Easyshade® Compact

After Shade Measurement, Prescribe VITA Materials for Superior Esthetic Results!





VITAVM®13 Restoration on tooth 9 Dentistry and ceramics by Dr. Ed McLaren





VITA**VM®9** restoration on tooth 8 Dentistry and ceramics by Dr. Ed McLaren





VITAVM®7 restoration on tooth 8 Dentistry by Dr. Todd Snyder, ceramics by Mark Baker, CDT





VITA**PM®9** restoration on tooth 8 Dentistry by Dr. Xavier Saab, ceramics by Victor Castro, CDT



Ordering Information

Part Number	Product Description
E75-201	VITA Easyshade® Compact





© Vident 2008. All rights reserved. L-9404V Rev 0908. Easyshade, 3D-Master, VITABLOCS, Bleachedguide, and Linearguide are trademarks or registered trademarks of VITA Zahnfabrik.


VITA Easyshade® Compact







Total Shade Control in the Palm of Your Hand.




Cordless, Sleek and Incredibly Fast.

The pressure of getting a perfect shade match has never been higher due to patients' increased attention to esthetic appearance. The VITA Easyshade Compact quickly and precisely measures 16 VITA Classical Shades, 29 VITA System 3D-Master® shades and even those hard-to-see intermediate shades. High-esthetics with VITA materials has never been easier to achieve.



Powerful Plug-and-Play Features.

The amazing design of Easyshade Compact hints at the powerful technology working underneath its sleek shell. Spectrophotometer technology, the latest OLED display and white-LED color-corrected lighting just for starters.





OLED display with easy-to-read icons

Versatility and Flexibility.

Measure a single central, check crown accuracy before seating, get a more detailed shade analysis with a detailed map — it's all here with Easyshade Compact. And with built-in memory that stores the last 25 measurements, you never have to worry about forgetting the result.


Select the appropriate measurement setting easily using the navigation buttons and easy-to-follow icons:

- Tooth single mode
- Tooth areas mode: maps the cervical, middle and incisal shades
- Verify restoration mode: check a restoration for shade accuracy
- Shade tab mode: practice proper placement of probe on VITA shade tab



Quality Control at the Push of a Button.

For the dental laboratory, Easyshade Compact verifies shade accuracy during the buildup process to determine where adjustment may be required.



By displaying which color parameters need attention (value, chroma or hue), Easyshade Compact makes it a simple matter to know precisely which VITA 3D-Master or VITA Classical shades should be used to make shade corrections.









VITA 3D-Master® Toothguide



VITA Bleachedguide®

Welcome to the Evolution.

Today, millions of satisfied patients around the world enjoy beautiful smiles from their dentists, aided by the superior esthetics of VITA shades and materials. Explore the entire system of VITA materials, including CAD/CAM blocks, all-ceramic and metal porcelains, denture teeth and composites.





VITABLOCS® Guide 3D-Master



VITA Linearguide™