

Welcome to the 2011 Edition of the Vident, a VITA Company, Product Catalog

Dear Dental Professional:

You've come to rely on VITA as the industry benchmark for quality and consistency. VITA materials offer seamless compatibility and flawless fabrication, regadless of material type, and now VITA Made restorations enable you to maintain esthetic integrity without sacrificing simplicity and productivity. We invite you to explore our product offerings in this catalog and to embrace the future of efficiency with VITA Made restorations and quality VITA products: from the Easyshade Compact shade-taking device and Vacumat 6000M furnace to the next generation of VMK porcelain – VMK Master. And remember, if it's a VITA Shade, it has to be VITA Made!

For those who are not familiar with Vident... for 26 years we have prided ourselves on offering unique products of the highest quality, including the world renowned VITA line of shade determination and restorative products. We are different from broad-line dental dealers in that we focus on innovative, specialty products, sourced from around the world or developed here at Vident. Our products are designed for dental labs and dentists that want to save time, enhance productivity, and provide the most natural looking, reliable outcome for your clients. These products are supported by Vident's extensive technical support and professional education resources.

For those who are already Vident customers... thank you for your continuing business. We know it is important to earn your loyalty with each order you place. Vident continues to invest heavily to improve and increase availability of technical support and professional education courses.



"We strive to be the best in all we do with one single ojective in mind: helping you succeed."

VIDENT

Fred Ketcho President

Vident Culture Book

At Vident, we have a new culture book that contains our core beliefs and values as well as provides the foundation of how we make our decisions, all with the goal of satisfying our customers. It contains our main ideals like Honesty and Integrity, Trust, Anticipating Customer Needs, Educational Leadership and more.

We take our culture very seriously and look forward to demonstrating our devotion to you and our core beliefs every day. You can access a downloadable PDF at www.vident.com Or call us at 800-828-3839 USA, 800-263-4778 Canada and will send you a free copy.

Where to find...

SHADE MANAGEMENT	3-4
METAL CERAMICS	5-17
ALL-CERAMICS AND PRESSABLES	18-29
STAINS	30-31
INDIRECT COMPOSITE	32-34
MACHINABLE CERAMICS	35-41
CHAIRSIDE PRODUCTS	42-44
DENTURE TEETH	45-47
FURNACES	48-50
ACCESSORIES	50-58
BOOKS	59

The prices shown in this catalog are effective January 1, 2010 through December 31, 2010

Who to call...

800-828-3839 in USA 800-263-4778 in Canada	
Place an order	Call 800
Check shipment status	number
Payment/Credit status	above
Return goods	
Course registration	Ext. 238
Technical help desk	(888) 249-1640
CEREC support	Ext. 240
Payments/Credits A-K	Ext. 252
Payments/Credits L-Z	Ext. 215

help@vident.com

Vident reserves the right to change prices or modify/discontinue products without prior notice.





VITA Easyshade® Compact

progress with the push of a button.

SHADE MANAGEMENT

The all-new VITA Easyshade Compact delivers the fastest, easiest way to measure tooth shades under any lighting condition. Easyshade Compact gives you the ability to conquer the trickiest shade challenges all day, every day, instantly and accurately. Now you can measure natural teeth, check restoration shades and track bleaching

The pressure of getting a perfect shade match has never been higher due to patients' increased attention to esthetic appearance. The VITA Easyshade Compact quickly and precisely

measures 16 VITA Classical Shades, 29 VITA System 3D-Master shades and even those hard-to-see intermediate shades. High esthetics with VITA materials has never been easier to achieve.

Measure a single central, check crown accuracy before seating, get a more detailed shade analysis with a detailed map — it's all here with Easyshade Compact. And with built-in memory that stores the last 25 measurements, you never have to worry about forgetting the result. Plus, Easyshade Compact seamlessly integrates with VITA ShadeAssist software for a complete shade communication solution.

Select the appropriate measurement setting easily using the navigation buttons and easy-to-follow icons:

- Tooth single mode
- Tooth areas mode: maps the cervical, middle and incisal shades
- Verify restoration mode: check a restoration for shade accuracy
- Shade tab mode: practice proper placement of probe on VITA shade tab

VITA Easyshade Compact Unit Order No. E75-201. . \$2,295.00

Reorder Items

Infection Control Sleeves, 160 pk Order No. E75-012. . \$66.00

VITA ShadeAssist Software

VITA ShadeAssist software allows the dental professional to enter VITA shades manually or import them from the VITA Easyshade Compact. ShadeAssist allows the user to create a detailed prescription and case planning document that can be printed or emailed.

Case-related, importable photos provide the ability to adjust images with cropping and contrast tools for optimal display of critical characterization details. ShadeAssist also features the ability to draw a detailed characterization sketch and add notes.

ShadeAssist can be used as stand-alone software for case planning and patient data storage or in conjunction with the VITA Easyshade Compact for

a complete digital shade solution.

VITA ShadeAssist Software Order No. GX10-003... \$99.00

Easyshade® Quality Manager Software

An efficient way to analyze the shade of ceramic restorations prior to staining and glazing. If shade correction is needed, Quality Manager recommends the necessary adjustment and can also provide a printed quality certificate. Easily connects to the VITA Easyshade for fast, accurate results.

Easyshade Quality Manager Software Order No. GX10-002... \$795.00

Easymatch™ Software

Easymatch Software saves time and money by predicting the optimal VITA 3D-Master shaded VITABLOCS before placing a milling block into your CEREC or Sirona inLab device.

Easymatch Software Order No. GX10-001. \$795.00 Patterson Order No. 706-8992

Orders: 800-828-3839 USA | 800-263-4778 Canada | www.vident.com



VITA Linearguide 3D-Master®

VITA Linearquide 3D-Master®

The all-new VITA Linearguide 3D-Master enables the quick determination of precise tooth shades and uses the same scientific principles found in the popular VITA 3D-Master shade guide. The Linearguide features a sleeker, linear design that makes the process of precise shade determination even faster and easier. In two simple steps the final shade is achieved, first by selecting from five value tabs, then by choosing the proper mix of chroma and hue within the selected value range.

The VITA Linearguide 3D-Master contains the same 29 shade samples as the VITA Toothquide 3D-Master.

Advantages:

- Classic linear handling
- Self-explanatory, simple use
- Quick, reliable and precise shade determination
- Modern design which reflects perfect esthetics and clear arrangement

VITA Linearquide 3D-Master Order No. GX69-913.. \$229.00



VITA Bleachedquide 3D-Master

VITA Bleachedguide 3D-Master®

The new VITA Bleachedquide 3D-Master offers a true alternative to the traditional VITAPAN Classical shade guide, and sets a new standard for the planning and monitoring of the patient's tooth whitening regimen. The VITA Bleachedguide 3D-Master also offers new intermediate shade

tabs that are arranged specifically for measuring bleaching progress.

VITA Bleachedguide 3D-Master Order No. GX69-911 \$103.00



VITA 3D-Master Shade Guide with bleached shade tab

VITA 3D-Master® Shade Guide

Dental professionals around the world have discovered precise and repeatable shade selection and improved communication using the VITA 3D-Master Tooth Guide. Based on true color science, shades can now be selected with incredible accuracy using a fast, efficient, three step process applying Value (lightness), Chroma (saturation), and Hue. VITA offers 29 shade tabs, including 3 bleached tabs in one autoclavable guide. The 3D-Master shade system is supported by VITA toothcolored restorative and prosthetic materials used for the veneering of metal ceramics, densely-sintered zirconia and alumina, as well as indirect laboratory composites, denture teeth and direct micro-hybrid composites.

VITA 3D-Master Tooth Guide With Bleached Cartridge	
Order No. GX69-910	\$103.00

VITA 3D-Master® Reorder Items

3D-Master Bleached Tooth Cartridge OM1, OM2, OM3	
Order No. GX69-950	\$22.60
3D-Master Tooth Cartridge 1M1, 1M2	
Order No. GX69-951	\$17.60
3D-Master Tooth Cartridge 2L1.5, 2L2.5	
Order No. GX69-952	\$17.60
3D-Master Tooth Cartridge 2M1, 2M2, 2M3	
Order No. GX69-953	\$22.60
3D-Master Tooth Cartridge 2R1.5, 2R2.5	
Order No. GX69-954	\$17.60
3D-Master Tooth Cartridge 3L1.5, 3L2.5	
Order No. GX69-955	\$17.60
3D-Master Tooth Cartridge 3M1, 3M2, 3M3	
Order No. GX69-956	\$22.60
3D-Master Tooth Cartridge 3R1.5, 3R2.5	
Order No. GX69-957	\$17.60
3D-Master Tooth Cartridge 4L1.5, 4L2.5	
Order No. GX69-958	\$17.60
3D-Master Tooth Cartridge 4M1, 4M2, 4M3	
Order No. GX69-959	\$22.60
3D-Master Tooth Cartridge 4R1.5, 4R2.5	
Order No. GX69-960	\$17.60
3D-Master Tooth Cartridge 5M1, 5M2, 5M3	
Order No. GX69-961	\$22.60



\$1,675.00

VITABLOCS® Guide 3D-Master®

Improves the shade selection by matching Vitablocs shades faster. The new Vitablocs shade tabs offer a realistic shade and translucency that's relevant to the final restoration's appearance. All the shade tabs are fabricated from Vitablocs material. The new VITABLOCS Guide is based on the original VITA 3D-Master system, reproducing the three natural tooth shade dimensions of value, Chroma and Hue. VITA future-oriented for more than 80 years

Order No. GX69-920	\$103.00
Uruer No. GX69-920	\$103.00

VITA Simulate Color Indicator

Ideal for preparation shade choice and comparisons.

Order No. GX70-430 . \$73.25

VITA Improved Classical Shade Guide

Order No. G68-990 . \$114.00

VITA has improved the VITA Classical Shade Guide to include additional features and benefits:

- The new holder is made with a stronger material which allows the entire guide (holder and tabs) to be autoclaved.
- Wider metal tabs increase the retention of the tab into the slots in the holder.

Improved Classical Reorder Shade Tabs

os		\$10.20
ade	Order No.	Shade
	G68-978	B4
	G68-979	C1
	G68-980	C2
.5	G68-981	C3
	G68-982	C4
	G68-983	D2
	G68-984	D3
	G68-985	D4
	ade 5	G68-978 G68-979 G68-980 5 G68-981 G68-982 G68-983 G68-984

Original Classical Reorder Shade Tabs

•				
Original Class	ical Tabs		\$10.20	
Order No.	Shade	Order No.	Shade	
G68-950	A1	G68-954	A4	
G68-958	B4	G68-962	C4	
G68-951	A2	G68-955	B1	
G68-959	C1	G68-963	D2	
G68-952	A3	G68-956	B2	
G68-960	C2	G68-964	D3	
G68-953	A3.5	G68-957	B3	
G68-961	C3	G68-965	D4	



VITABLOCS Guide 3D-Master



VITA Simulate Color Indicator





VITA VM_®13

VITA's VM13 ceramic veneering porcelain offers outstanding physical properties, handling characteristics and, due to the homogeneous distribution of glasses in the fine microstructure, clinical wear characteristics that mimic those of enamel.

VM13 has a CTE of 13.1-13.6 x 10-6K-1 and is ideal for veneering metal frameworks that have a coefficient of thermal expansion in the range of 13.8–15.2 x 10-6K⁻¹. It is also an excellent material for laminated veneers.

The revolutionary formulation of VITA VM13:

- Low wear rate
- Offers optimum firing stability, which preserves shades even after several firings.
- Has a low firing temperature and reduced shrinkage
- Offers excellent shade cover with the new higher chroma, higher opacity Base Dentines
- Matches dentition with accurate, consistent 3D-Master and Classical shades

Dental ceramists will marvel at its modeling characteristics and its structural and color stability through its firing cycles. Grinding and polishing has been made easier and more effective than most other porcelain systems. The esthetic build-up technique has been improved over historical systems by introducing a simpler, more logical approach to shading/layering.

VM_®13 Basic Kit - 12g Powder Opaque

Order No. VX50-0101

. \$1,405.00 VITA VM 13 Wash Opaque Powder

- VITA VM 13 Opaque Powder
- VITA VM 13 Sun Opaque
- VITA VM 13 Base Dentine Powder VITA VM 13 Chroma Plus Powder
- VITA VM 13 Sun Dentine
- VITA VM 13 Enamel Powder (ENL, END)
- VITA VM 13 Neutral Powder NT
- VITA VM 13 Window Powder WIN
- VITA VM 13 Corrective Powder
- VITA VM 13 Opaque Fluid, 50ml
- Modeling Liquid, 50ml
- Modisol Pen
- Crown Stand (A & B)
- Firing Tray G
- 3D-Master Tooth Guide
- VITA VM 13 Shade Indicator Chart
- Working Instructions


VM_®13 Basic Kit - 12g Paste Opaque

Order No. VX50-0102.

- VITA VM 13 Wash Opaque Paste (7g)
- VITA VM 13 Paste Opaque (5q)
- VITA VM 13 Base Dentine Powder
- VITA VM 13 Chroma Plus Powder VITA VM 13 Enamel Powder ENL, END
- VITA VM 13 Neutral Powder NT
- VITA VM 13 Window Powder WIN VITA VM 13 Corrective Powder
- VITA VM 13 Paste Fluid, 15ml
- Modeling Liquid, 50ml
- Modisol Pen
- Paste Opaque Brush
- Spatula
- Crown Stand (A and B)
- Firing Tray G
- 3D-Master Tooth Guide
- VITA VM13 Shade Indicator Chart
- Working Instructions

VM_®13 Basic Kit (Small) - 12g Powder Opaque

Order No. VX50-0103... VITA VM 13 Wash Opaque Powder



- VITA VM 13 Base Dentine Powder
- VITA VM 13 Chroma Plus Powder
- VITA VM 13 Enamel Powder ENL, END
- VITA VM 13 Neutral Powder NT
- VITA VM 13 Window Powder WIN
- VITA VM 13 Corrective Powder
- VITA VM 13 Opaque Fluid, 50ml
- Modeling Liquid, 50ml
- Modisol Pen
- Crown Stand (A and B)
- Firing Tray G
- 3D-Master Tooth Guide
- VITA VM13 Shade Indicator Chart
- Working Instructions

VM_®13 Basic Kit (Small) - 12g Paste Opaque

Order No. VX50-0104... . \$1,260.00

VITA VM13 Wash Opaque Paste (7g)

VITA VM13 Paste Opaque (5g)

VITA VM13 Base Dentine Powder

VITA VM13 Chroma Plus Powder

VITA VM13 Enamel Powder ENL, END

VITA VM13 Neutral Powder NT

VITA VM13 Window Powder WIN

VITA VM13 Corrective Powder

VITA VM13 Paste Fluid, 15ml

Modeling Liquid, 50ml

Modisol Pen

Paste Opaque Brush

Spatula

Crown Stand (A and B)

Firing Tray G

3D-Master Tooth Guide

VITA VM13 Shade Indicator Chart

Working Instructions





VM_®13 Buildup Kit, 12g

der No.	VX50-0005	\$6
26	VITA VM13 Transpa Dentine Powder	
1	Modeling Liquid, 50ml	Day.
1	Working Instructions	-



VM_®13 Buildup Kit (Small), 12g

Order No. VX50-0006\$398.00		
15	VITA VM13 Transpa Dentine Powder	
1	Modeling Liquid, 50ml	
1	Working Instructions	

VM_®13 Bleached Shade Kit 12g Powder Opaque

Order No.	VX50-0007	\$281.00
1	VITA VM13 Opaque Powder (OP0)	
3	VITA VM13 Base Dentine Powder	No.
3	VITA VM13 Transpa Dentine Powder	All
1	VITA VM13 Enamel Light ENL	1000
1	VITA VM13 Neutral Powder NT	
1	VITA VM13 Window Powder WIN	
1	Modeling Liquid, 50ml	
1	Opaque Fluid, 50ml	
1	3D-Master Bleach Guide	

VM_®13 Bleached Shade Kit, Paste Opaque

Working Instructions

Order No	VX50-0008	\$340.00
1	VITA VM13 Paste Opaque (5q) (OPO)	φυτο.σο
3	VITA VM13 Base Dentine Powder	
3	VITA VM13 Transpa Dentine Powder	
1	VITA VM13 Enamel Light ENL	
1	VITA VM13 Neutral Powder NT	
1	VITA VM13 Window Powder WIN	
1	Modeling Liquid, 50ml	
1	VITA VM13 Paste Fluid, 15ml	
1	Paste Opaque Brush	
1	Spatula	

VM₀12 Professional Vit 12a

6

3D-Master Bleach Guide

Working Instructions

V IVI® I 3	Professional Kit, 12g	
Order No.	VX50-0009	\$1,745.00
11	VITA VM13 Effect Chroma	
11	VITA VM13 Effect Enamel	
3	VITA VM13 Mamelon	The same of the sa
3	VITA VM13 Effect Pearl	A SECTION A
5	VITA VM13 Effect Opal	1 Sept.
6	VITA VM13 Effect Liner	
3	VITA VM13 Color Indicators	
1	Working Instructions	

VM_®13 Professional Kit (Small), 12g

VX50-0010	\$575.00
VITA VM13 Effect Chroma	
VITA VM13 Mamelon	
VITA VM13 Effect Pearl	
VITA VM13 Effect Opal	
VITA VM13 Effect Enamel	
VITA VM13 Color Indicator	
Working Instructions	
	VITA VM13 Effect Chroma VITA VM13 Mamelon VITA VM13 Effect Pearl VITA VM13 Effect Opal VITA VM13 Effect Enamel VITA VM13 Color Indicator

VM_®13 Gingiva Kit **12g Powder Opaque**

rder No.	VX50-0017	\$246.00
5	VITA VM13 Gingiva	
2	VITA VM13 Gingiva Opaque	
1	VITA VM13 Color Indicator	A STATE OF THE PARTY OF THE PAR
1	Working Instructions	Aller to

VM_®13 Gingiva Kit Paste Opaque

Order No.	VX50-0018	\$283.00
5	VITA VM13 Gingiva	
2	VITA VM13 Gingiva Paste Opaque (5g)	
1	VITA VM13 Color Indicator	
1	Working Instructions	

VM_®13 Margin Kit, 12g

Order No.	. VX50-0019	\$254.00
6	VITA VM13 Margin	
1	Color Indicator	
1	Working Instructions	
/B.B. 44) F((4 01 V)4	

VM_®13 Effect Opal Kit

Order No.	VX50-0020	\$ FP (
5	VITA VM13 Effect Opal	
1	Shade Indicator	
1	Working Instructions	

VM_®13 Sun Kit Order No. VX50-0021.

3	Sun Opaque	
3	Sun Dentine	
1	Shade Indicator	
1	Working Instructions	

..\$**FPO**

VM_®13 One Shade Kits, 12g

Order No.	Description	
VX50-0011	Shade 2M2 Powder Opaque\$204	.00
VX50-0012	Shade 3M2 Powder Opaque\$204	.00
VX50-0013	Shade 2M2 Paste Opaque\$324	.00
VX50-0014	Shade 3M2 Paste Opaque\$324	.00
VX50-0020	Effect Opal Kit - NEW	
VX50-0021	Sun Kit - NEW	

VM®13 Reorder Materials

VM 13 Margir Order No. VX50-2261 VX50-2264 VX50-2265 VX50-2267 VX50-2268 VX50-2269	n Porcelain, 12g Shade MAR 1 MAR 4 MAR 5 MAR 7 MAR 8 MAR 9		\$38.20
VM13 Wash C)paques		
Order No.	Description		****
VX50-2359	Paste Wash Opaque		
VX50-2259 VX51-2259	Powder Wash Opaqu Powder Wash Opaqu		
	paques 5g		\$76.25
Order No.	Shade		
VX50-2350	OPO		
VX50-2351	OP1		
VX50-2352	OP2		
VX50-2353 VX50-2354	OP3 OP4		
VX50-2354 VX50-2355	0P5		
VX30-2333	013		
	Opaque Paste, 5g		\$76.25
Order No.	Shade		
VX50-2371	GOL (Gingiva Opaque		
VX50-2373	GOD (Gingiva Opaqu	e Dark)	
VM13 Gingiya	Opaque Powder, 12g		\$34.00
Order No.	Shade		φυτ.υυ
VX50-2271	GOL (Gingiva Opaque	e Light)	
VX50-2273	GOD (Gingiva Opaqu		
			*
	paque Powder, 12g		\$34.00
Order No.	Shade		
VX50-2281 VX50-2282	CO1 CO2		
VX50-2282 VX50-2283	CO3		
V//30 ZZ03	000		
VM13 Color O	paque Paste, 5g		\$76.25
Order No.	Shade		
VX50-2381	CO1		
VX50-2382	CO2		
VX50-2383	CO3		
VM13 Opaque	Powders		
			\$34.00
12g Order No.	•	Shade	
VX50-2250		OP0	
VX50-2251	VX51-2251	OP1	
VX50-2252	VX51-2252	OP2	
VX50-2253 VX50-2254	VX51-2253 VX51-2254	OP3 OP4	
VX50-2254 VX50-2255	v/\J1-ZZJ4 —	0P5	
V/100 ZZ00		0.0	

Technical Help Desk: 888-249-1640 USA | 800-324-6224 Canada

12g Order No.		250g Order No.	
VX50-2031	_	_	0M1
VX50-2032	_	_	0M2
VX50-2033			0M3
VX50-2034	VX51-2034	_	1M1
VX50-2034 VX50-2035	VX51-2034 VX51-2035		1M2
VX50-2035 VX50-2036	VX51-2035 VX51-2036		2L1.5
		_	2L1.5 2L2.5
VX50-2037	VX51-2037	_	
VX50-2038	VX51-2038		2M1
VX50-2039	VX51-2039	VX52-2039	2M2
VX50-2040	VX51-2040	VX52-2040	2M3
VX50-2041	VX51-2041	_	2R1.5
VX50-2042	VX51-2042	_	2R2.5
VX50-2043	VX51-2043	_	3L1.5
VX50-2044	VX51-2044	_	3L2.5
VX50-2045	VX51-2045		3M1
VX50-2046	VX51-2046	VX52-2046	3M2
VX50-2047	VX51-2047	VX52-2047	3M3
VX50-2048	VX51-2048	_	3R1.5
VX50-2049	VX51-2049	_	3R2.5
VX50-2050	—	_	4L1.5
VX50-2051	_	_	4L2.5
VX50-2051 VX50-2052	VX51-2052	_	4M1
VX50-2052 VX50-2053	VX51-2052 VX51-2053		4M2
VX50-2053 VX50-2054	VX51-2053 VX51-2054	_	4M3
	VA31-2034	_	
VX50-2055	_	_	4R1.5
VX50-2056	_	_	4R2.5
VX50-2057	_	_	5M1
VX50-2058			5M2
VX50-2059	_	_	5M3
VM13 Dentine T			
12g			\$34.
			\$430. Shade
	50g Order No.	-	
VX50-2061	_	_	0M1
VX50-2062	_		0M2
VX50-2063	_	_	0M3
WED 0004			1M1
	VX51-2064	_	
VX50-2065	VX51-2065	_	1M2
VX50-2065 VX50-2066	VX51-2065 VX51-2066	_ _ _	1M2 2L1.5
VX50-2065 VX50-2066 VX50-2067	VX51-2065 VX51-2066 VX51-2067	_ _ _	1M2 2L1.5 2L2.5
VX50-2065 VX50-2066 VX50-2067	VX51-2065 VX51-2066	_ _ _ _	1M2 2L1.5
VX50-2065 VX50-2066 VX50-2067 VX50-2068	VX51-2065 VX51-2066 VX51-2067		1M2 2L1.5 2L2.5
VX50-2065 VX50-2066 VX50-2067 VX50-2068 VX50-2069	VX51-2065 VX51-2066 VX51-2067 VX51-2068		1M2 2L1.5 2L2.5 2M1
VX50-2065 VX50-2066 VX50-2067 VX50-2068 VX50-2069 VX50-2070	VX51-2065 VX51-2066 VX51-2067 VX51-2068 VX51-2069		1M2 2L1.5 2L2.5 2M1 2M2
VX50-2065 VX50-2066 VX50-2067 VX50-2068 VX50-2069 VX50-2070 VX50-2071	VX51-2065 VX51-2066 VX51-2067 VX51-2068 VX51-2069 VX51-2070 VX51-2071		1M2 2L1.5 2L2.5 2M1 2M2 2M3 2R1.5
VX50-2064 VX50-2065 VX50-2066 VX50-2067 VX50-2068 VX50-2069 VX50-2070 VX50-2071 VX50-2072	VX51-2065 VX51-2066 VX51-2067 VX51-2068 VX51-2069 VX51-2070 VX51-2071 VX51-2072		1M2 2L1.5 2L2.5 2M1 2M2 2M3 2R1.5 2R2.5
VX50-2065 VX50-2066 VX50-2067 VX50-2068 VX50-2069 VX50-2070 VX50-2071 VX50-2072 VX50-2073	VX51-2065 VX51-2066 VX51-2067 VX51-2068 VX51-2069 VX51-2070 VX51-2071 VX51-2072 VX51-2073		1M2 2L1.5 2L2.5 2M1 2M2 2M3 2R1.5 2R2.5 3L1.5
VX50-2065 VX50-2066 VX50-2067 VX50-2068 VX50-2069 VX50-2070 VX50-2071 VX50-2072 VX50-2073 VX50-2074	VX51-2065 VX51-2066 VX51-2067 VX51-2068 VX51-2069 VX51-2070 VX51-2071 VX51-2072 VX51-2073 VX51-2074		1M2 2L1.5 2L2.5 2M1 2M2 2M3 2R1.5 2R2.5 3L1.5 3L2.5
VX50-2065 VX50-2066 VX50-2067 VX50-2068 VX50-2069 VX50-2070 VX50-2071 VX50-2072 VX50-2073 VX50-2074 VX50-2075	VX51-2065 VX51-2066 VX51-2067 VX51-2068 VX51-2069 VX51-2070 VX51-2071 VX51-2072 VX51-2073 VX51-2074 VX51-2075	VX52-2070 — — — — —	1M2 2L1.5 2L2.5 2M1 2M2 2M3 2R1.5 2R2.5 3L1.5 3L2.5 3M1
VX50-2065 VX50-2066 VX50-2067 VX50-2068 VX50-2069 VX50-2070 VX50-2071 VX50-2072 VX50-2073 VX50-2074 VX50-2075 VX50-2076	VX51-2065 VX51-2066 VX51-2067 VX51-2068 VX51-2069 VX51-2070 VX51-2071 VX51-2072 VX51-2073 VX51-2074 VX51-2075 VX51-2076	VX52-2070 — — — — — VX52-2076	1M2 2L1.5 2L2.5 2M1 2M2 2M3 2R1.5 2R2.5 3L1.5 3L2.5 3M1 3M2
VX50-2065 VX50-2066 VX50-2067 VX50-2068 VX50-2069 VX50-2070 VX50-2071 VX50-2072 VX50-2073 VX50-2074 VX50-2075 VX50-2076 VX50-2077	VX51-2065 VX51-2066 VX51-2067 VX51-2068 VX51-2070 VX51-2071 VX51-2072 VX51-2073 VX51-2074 VX51-2075 VX51-2076 VX51-2076	VX52-2070 — — — — —	1M2 2L1.5 2L2.5 2M1 2M2 2M3 2R1.5 2R2.5 3L1.5 3L2.5 3M1 3M2 3M3
VX50-2065 VX50-2066 VX50-2067 VX50-2068 VX50-2069 VX50-2070 VX50-2071 VX50-2072 VX50-2073 VX50-2074 VX50-2075 VX50-2076 VX50-2077	VX51-2065 VX51-2066 VX51-2067 VX51-2068 VX51-2070 VX51-2071 VX51-2072 VX51-2073 VX51-2074 VX51-2075 VX51-2076 VX51-2077 VX51-2077	VX52-2070 — — — — — VX52-2076	1M2 2L1.5 2L2.5 2M1 2M2 2M3 2R1.5 2R2.5 3L1.5 3L2.5 3M1 3M2 3M3 3R1.5
VX50-2065 VX50-2066 VX50-2067 VX50-2068 VX50-2069 VX50-2070 VX50-2071 VX50-2072	VX51-2065 VX51-2066 VX51-2067 VX51-2068 VX51-2070 VX51-2071 VX51-2072 VX51-2073 VX51-2074 VX51-2075 VX51-2076 VX51-2076	VX52-2070 — — — — — VX52-2076	1M2 2L1.5 2L2.5 2M1 2M2 2M3 2R1.5 2R2.5 3L1.5 3L2.5 3M1 3M2 3M3

VX50-2082 VX50-2083 VX50-2084 VX50-2085 VX50-2086 VX50-2087 VX50-2088 VX50-2089	VX51-2082 VX51-2083 VX51-2084 ————————————————————————————————————		4M1 4M2 4M3 4R1.5 4R2.5 5M1 5M2 5M3
VM13 Window Order No. VX50-2181 VX51-2181 VX52-2181	Description		\$120.00
50g 250g 12g Order No. VX50-2191		250g Order No. VX52-2191 VX52-2192	\$120.00 \$436.00
VM13 Neutral Order No. VX50-2201 VX51-2201 VX52-2201	Description 12g 50g		\$120.00
VM13 Effect E Order No. VX50-2091 VX50-2092 VX50-2093 VX50-2094 VX50-2095 VX50-2096	namel 12g Shade EE1 EE2 EE3 EE4 EE5 EE6	Order No. VX50-2097 VX50-2098 VX50-2099 VX50-2100 VX50-2101	\$38.20 Shade EE7 EE8 EE9 EE10 EE11
VM13 Effect C Order No. VX50-2121 VX50-2122 VX50-2123 VX50-2124 VX50-2125 VX50-2126	Shade	Order No. VX50-2127 VX50-2128 VX50-2129 VX50-2130 VX50-2131	Shade
VM13 Effect Li Order No. VX50-2141 VX50-2142 VX50-2143	Shade EL1	Order No. VX50-2144 VX50-2145 VX50-2146	Shade EL4 EL5
Order No. VX50-2151 VX50-2152	Shade EP1		\$38.20

			\$38.20
Order No. VX50-2171	Shade		
VX50-2172			
VX50-2173			
VX50-2174	E04		
VX50-2175	E05		
VM13 Gingiva	a, 12g		\$38.20
Order No.	Shade		
VX50-2231	G1		
VX50-2232	G2 G3		
VX50-2233	G3		
VX50-2234	G4		
VX50-2235	G5		
VM13 Chroma	a Plus, 12g		\$38.20
Order No.	Shade		
VX50-2241	CP1		
VX50-2242	CP2		
VX50-2243	CP3		
VX50-2244	CP4		
	CP5		
VM13 Mamel	on,12g		\$38.20
Order No.			
VX50-2161			
VX50-2162	MM2		
VX50-2162 VX50-2163	MM2 MM3		
VM13 Correct	ive. 12a		\$38.20
Order No.			
VX50-2221	COR1 Incisal	White	
VX50-2222	COR1 Incisal COR2 Body		
VX50-2223	COR3 Neck		
VM13 Sun De	ntine. 12a		\$38.20
			,
Order No. VX50-2361	Sunlight		
VX50-2362	Sunrise		
VX50-2363	Sunset		
VM13 Sun Or	nague, 12g		\$34.00
Order No.	Shade		
VX50-2291	Yellowish		
VX50-2292	Orang		
VX50-2293	•		
VM_12 C	lassical Sh	ados	
	asic Shade		\$1,475.00
			Opaque A1-D4
ı vvasn U	uduue	ın	UUAUUE AT-174

Opaque A1-D4 Sun Dentine Wash Opaque Base Dentine A1-D4 3 Enamel ENL, END Chroma PlusL Neutral NT Window Modisol Pen Corrective Firing Tray G Shade indicator Modeling Liquid, 50ml 1 Crown Stands (A,B,C) 1 Classical shade guide

VIVI® IJ CIASSICAI JIIAUGS	VM _® 13	Classical	Shades
----------------------------	--------------------	-----------	--------

* ivi	assivai ona	iucs	
V50-0002	VM13 Basic Ki	t Paste Classical	\$2 140 00
V50-0003	VIVITS BASIC KI	t Small Classical	\$900.00
V50-0004	VM13 Basic Ki	t Sml. Paste Classic	\$1,195.00
V50-0005	VM13 Build Up	Kit Classical	\$426.00
V50-0006		Kit Small Classical	
V50-0007		l Color Kit Paste	
V50-0008		I Color Kit Powder	
V50-0011	VM13 One Col	or Kit A2	\$182.00
V50-0012		or Kit A3	
V 30 0012	VIVITO ONG OOI	01 Kit A0	ψ102.00
\	10		404.00
Order No.	Shade	Order No.	Shade
V50-2401	A1	V50-2409	B4
V50-2402	A2	V50-2410	C1
	Λ <u>2</u>		C2
V50-2403	A3	V50-2411	
V50-2404	A3.5	V50-2412	C3
V50-2405	A4	V50-2413	C4
V50-2406	B1	V50-2414	D2
V50-2407	B2	V50-2415	D3
V50-2408	B3	V50-2416	D4
VM13 Opaque	e, 50g		\$120.00
Order No.	Shade	Order No.	Shade
V51-2401	A1		B4
		V51-2409	
V51-2402	A2	V51-2410	C1
V51-2403	A3	V51-2411	C2
V51-2404	A3.5	V51-2412	C3
V51-2405	A4	V51-2413	C4
V51-2406	B1	V51-2414	D2
V51-2407	B2	V51-2415	D3
V51-2408	В3	V51-2416	D4
VM13 Onague	Pasta Sa		\$76.25
	Shade	Order No.	Shade
V50-2301	A1	V50-2309	B4
V50-2302	A2	V50-2310	C1
V50-2303	A3	V50-2311	C2
V50-2304	A3.5	V50-2312	C3
			C4
V50-2305	A4	V50-2313	
V50-2306	B1	V50-2314	D2
V50-2307	B2	V50-2315	D3
V50-2308	В3	V50-2316	D4
\/N/12 Sun On	ianiia Paeto En		¢ 7ፎ ን፫
			φ/0.∠3
Order No.	Shade		
VX50-2391	Yellowish		
VX50-2392	Orang		
VX50-2393	Redish		
V/100 2000	Houisii		
\/\/112 D D	ontino 10-		#04.00
Order No.	Shade	Order No.	Shade
V50-2431	A1	V50-2439	B4
V50-2432	A2	V50-2440	C1
	A3	V50-2441	C2
V50-2433			
V50-2434	A3.5	V50-2442	C3
V50-2435	A4	V50-2443	C4
V50-2436	B1	V50-2444	D2
V50-2437	B2	V50-2445	D3
V50-2438	B3	V50-2446	D4

VM13 Base Den	tine, 50g		\$120.00
Order No.		Order No.	Shade
V51-2431	A1	V51-2439	B4
V51-2432	A2	V51-2440	C1
V51-2433	A3	V51-2441	C2
V51-2434	A3.5	V51-2442	C3
V51-2435	A4	V51-2443	C4
V51-2436	B1	V51-2444	D2
V51-2437	B2	V51-2445	D3
V51-2438	B3	V51-2446	D4
VM13 Transpa [Dentine 12g		\$34.00
Order No.			Shade
V50-2451	A1	V50-2459	B4
V50-2452	A2	V50-2460	C1
V50-2453	A3	V50-2461	C2
V50-2454	A3.5	V50-2462	C3
V50-2455	A4	V50-2463	C4
V50-2456	B1	V50-2464	D2
V50-2457	B2	V50-2465	D3
V50-2458	B3	V50-2466	D4
VM13 Transpa [Dentine 50g		\$120.00
Order No.	Shade	Order No.	Shade
V51-2451	A1	V51-2459	B4
V51-2452	A2	V51-2460	C1
V51-2453	A3	V51-2461	C2
V51-2454	A3.5	V51-2462	C3
V51-2455	A4	V51-2463	C4
V51-2456	B1	V51-2464	D2
V51-2457	B2	V51-2465	D3
V51-2458	B3	V51-2466	D4
uni I £1⊛MV	iids and Acc	essories	
o = iqt			

V68-909 V69-909 VX50-2320 VX50-2321 V68-907	VM Modeling Liquid, 50 ml	\$78.50 \$21.70 \$78.50
VM Color Indica	itors	\$71.25
Order No.	Description	
GX78-414	Effect Chroma, 11 colors	
GX78-415	Effect Liners, EL1 – EL6	
	Effect Pearl	
	Mamelon, 15 colors	
GX78-416	Effect Enamel, 11 colors	
GX78-417	Gingiva, 7 colors	
GX78-418	VM Margin Color Indicator, 6 colors	
GX78-419	Professional Kit Small, 15 colors	
	Effect Chroma, EC1, 4, 6, 8, 9	
	Mamelon, MM2	
	Effect Pearl, EP1	
	Effect Opal, E02	
	Effect Enamel, EE1, 3, 7, 8, 9, 10, 11	



. \$22.70

CERAMICS

TA VM13	EFFEC.	T LINER	
	L1	Snow	White shade for 0-group (bleached
E	L2	Cream	Beige
	L3	Tabac	Brown
E	L4	Golden Fleece	Yellow
E	L5	Papaya	Orange
E	L6	Sesame	Green-yellow
A VM13	EFFEC	T CHROMA	
	C1	Ghost	For white effects
E	C2	Linen	Sand-beige
E	C3	Pale Banana	Light-yellow
E	C4	Lemon Drop	Delicate lemon yellow
E	C5	Golden Rod	Light-orange
E	C6	Sunflower	Orange
E	C7	Light Salmon	Pink
E	C8	Toffee	Beige-brown
	C9	Doe	Brown
	C10	Larch	To reproduce green-brown effects
	C11	Gravel	To reproduce green-grey effects
TA VM13	MAME	LON	
	MM1	Ecru	Beige
	MM2	Mellow Buff	Warm yellow-brown
N	MM3	Peach Puff	Delicate orange

METAL CERAMICS

	IVIIVIZ	IVICIIOVV DUIT	vvaiiii ychovi broviii
	MM3	Peach Puff	Delicate orange
VITA VM	13 GINGI	/A	
	G1	Rose	Old rose
	G2	Nectarine	Orange-red
	G3	Pink Grapefruit	Pink
	G4	Rosewood	Brown-red
	G5	Cherry Brown	Black-red
	GOL	Light Flesh	Light-red

VITA VM13 CORRECTIVE			
	COR1	Neutral	Neutral
	COR2	Sand	Beige
	COR3	Ochre	Brown

Dark-red

GOD

VITA VM13 SUN DENTINE NEW

10

Dark Flesh

Sunlight	For a sunnier and warmer shade result, to be used in the dentine and body area	
Sunrise		
Sunset		

VITA VM13 SUN OPAQUE NEW			
	Yellowish	for mixing the respective opaque shades, in three different nuances	
	Orange		
	Reddish		

VITA VN	113 EFFE	CT ENAMEL	
	EE1	Mint Cream	For whitish incisal areas, margin ridges and triangular protuberances
	EE2	Pastel	For ridges and wedges at the labial surfaces
	EE3	Misty Rose	For pink-translucent discoloratio in the cervical and/or enamel are
	EE4	Vanilla	For ridges and wedges
	EE5	Sun Light	Yellowish-translucent to intensif depth in the approximal area
	EE6	Navajo	For reddish-translucent discoloration, especially in older teeth
	EE7	Golden Glow	Orange-translucent to intensify t effect of depth
	EE8	Coral	Red-translucent for cervical area and/or severe discolorations
	EE9	Water Drop	For subcoating in the enamel are for bluish incisal areas
	EE10	Silver Lake Blue	For bluish wedges and ridges
	EE11	Drizzle	For greyish wedges to increase to effect of depth
VITA VI	113 EFFE	CT PEARL	
	EP1	Pearl	Shade nuance in pastel-yellow
	EP2	Pearl Blush	Shade nuance in pastel-orange
	EP3	Pearl Rose	Shade nuance in pastel-rose
VITA VN	113 FFFF	CT OPAL New	
	EE01	Neutral	Neutral, universally suitable
	E02	Whitish	Distinctive opal effect, for whiti
			wedges and triangular
	E03	Bluish	Bluish
	E04	Blue	Blue
	E05	Dark Violet	Dark Violet
VITA VI	113 MAR	GIN	
	M1	Icy Beige	White
	M4	Wheat	Yellow
	M5	Amber	Amber
	M7	Seashell	Light-Beige
	M8	Tan	Pastel Brown
	M9	Beach	Light-Orange
VITA VIN	112 CUDO	OMA PLUS	0 0
VIIA VIV	CP1	lvory	lvory
	CP2	Almond	Beige
	CP3	Moccasin	Light Orange-Brown
	CP4	Caramel	Orange
	014	Caraniel	orange

VITA VMK®Master Porcelain

This new third generation VITA VMK porcelain offers the lab technician, the dentist and the patient exceptional physical properties, handling characteristics and, due to its fine structure, clinical wear characteristics that mimic natural tooth enamel

From simple minimal powder buildups to the more complex cases with individual characteristics, VITA VMK Master's optimum firing temperature and chemical and physical properties make it suitable for veneering nearly all types of dental alloys, including non-precious. Working with the new VITA VMK Master offers some important advantages, including excellent working properties, improved modeling characteristics, easy and fast handling and minimal shrinkage. Accurate and outstanding esthetic results with traditional layering.

VMK Master offers a selection of additional porcelains, including Luminaries to enhance fluorescence, Translucent porcelains to reproduce the translucency of natural tooth enamel, Neck Modifiers to simulate the cervical third, Gingival porcelains to simulate tissue color, and Opal and Pearl modifiers to achieve true youthful and bleached teeth.

Available in VITA 3D-Master and VITA Classical A1 – D4 shades. A CTE of 13.1 – 13.6 x 10.6 K-1, ideal for veneering metal frameworks that have a coefficient of thermal expansion (CTE) in the range of 13.8 – 15.2 x 10⁻⁶.

. \$1,850.00

\$2,250.00

VMK_®Master 3D-Master **Standard Set, Powder 12g**

Order No.	VX70-001	
1	Wash Onaque	WC

- Opaque, OP1-OP5
- Opaque Dentine, 1M1-5M3
- Dentine, 1M1-5M3
- Enamel EN1, EN3
- Translucent, T4 Window, WIN
- Cervical, CE1-CE3
- Neck, N1-N5
- Corrective, COR1, COR2
- 50ml VITA Modelling Fluid
- 50ml VITA Modelling Fluid RS
- 50ml VITA Opaque Fluid
- Shade indicator 3D-MASTER
- VITA Linearguide 3D-MASTER
- Accessories

VMK_®Master 3D-Master **Standard Set, Paste 12g**

Order No. VX70-002...

- Wash Opaque, WO
- Opaque, OP1-OP5
- Opaque Dentine, 1M1-5M3
- Dentine, 1M1-5M3 Enamel EN1, EN3
- Translucent, T4
- Window, WIN
- Cervical, CE1-CE3
- Neck, N1-N5
- Corrective, COR1, COR2
- 50ml VITA Modelling Fluid
- 50ml VITA Modelling Fluid RS
- 50ml VITA Paste Opaque Liquid Shade indicator 3D-MASTER
- VITA Linearguide 3D-MASTER
- Accessories

VMK_®Master 3D-Master 10-Color Set, Powder 12g

Order No. VX70-003.. . \$1,095.00

- Wash Opaque, WO
- Opaque, OP1-OP4
- Opaque Dentine, 1M1-5M3
- Dentine (1M2, 2L1.5, 2M1, 2M2, 2R2.5, 3L1.5, 3M1,
 - 3M2, 3M3, 4M2)
- Enamel EN1, EN3 Translucent T4
- Window WIN
- Cervical CE1-CE3
- Neck N1-N5
- Corrective COR1, COR2
- 50 ml VITA Modelling Fluid
- 50 ml VITA Modelling Fluid RS
- 50 ml VITA Opaque Fluid
- Shade indicator 3D-MASTER
- VITA Linearguide 3D-MASTER
- Accessories

VMK_®Master 3D-Master 10-Color Set, Paste 12q

Order No. VX70-004. . \$1,495.00

- Wash Opaque, WO
- Opaque, OP1-OP5 5
- Opaque Dentine, 1M1-5M3
- Dentine (1M2, 2L1.5, 2M1, 2M2, 2R2.5, 3L1.5, 3M1, 3M2, 3M3, 4M2)
- Enamel EN1, EN3
- Translucent T4
- Window WIN
- Cervical CE1-CE3
- Neck N1-N5
- Corrective COR1, COR2
- 50 ml VITA Modelling Fluid
- 50 ml VITA Modelling Fluid RS
- 50 ml VITA Paste Opaque Liquid
- Shade indicator 3D-MASTER
- VITA Linearguide 3D-MASTER
- Accessories



Orders: 800-828-3839 USA | 800-263-4778 Canada | www.vident.com

CO1

CO2

CO3

Gold

Brown

Lilac

Orange

Brown

Lilac

. \$69.95

VMK®Master 3D-Master Bleached Color Set, Powder 12g

Order No. VX70-007. \$239.00 Opaque, OPO 1 3 Opaque Dentine, 0M1-0M3 Dentine, 0M1-0M3 Enamel, EN1 Translucent, T4 VITA Modelling Fluid RS, 50ml VITA Opaque Fluid, 50ml Bleached Shade Guide, 0M1-0M3

VMK_®Master 3D-Master **Bleached Color Set, Paste 12g**

Order No. VX70-008. . \$279.00 Opaque, OPO 3 Opaque Dentine, 0M1-0M3 Dentine, 0M1-0M3 3 Enamel, EN1 1 Translucent, T4 VITA Modelling Fluid RS, 50ml 1 VITA PASTE Opaque Liquid, 50ml Bleached Shade Guide, 0M1-0M3

VMK_®Master 3D-Master or Classical **Gingiva Set, Powder 12g**

Order No. V70-010... . \$199.00 2 Opaque, GOL, GOD 5 Gingiva, G1-G5

VMK_®Master 3D-Master or Classical **Gingiva Set, Paste 12g**

Margin/Gingiva shade guide

Order No. V70-011. \$259.00 2 Opaque, GOL, GOD Gingiva, G1-G5 Margin/Gingiva shade guide

VMK®Master 3D-Master Starter Set, Powder (2M2, 3M2) 12g

VITA Modelling Fluid RS, 50ml

VITA Opaque Fluid, 50ml

Shade tabs, 2M2, 3M2

12

Order No. VX70-005. . \$99.00 Opaque, OP2, OP3 2 Opaque Dentine, 2M2, 3M2 Dentine, 2M2, 3M2 Enamel, EN1 Translucent, T4

VMK®Master 3D-Master Starter Set, Paste (2M2, 3M2) 12g

Order No. VX70-006.. . \$139.00 Opaque, OP2, OP3 Opaque Dentine, 2M2, 3M2 Dentine, 2M2, 3M2 Enamel, EN1 Translucent, T4 VITA Modelling Fluid RS, 50ml VITA Paste Opaque Liquid, 50ml Shade tabs, 2M2, 3M2

VMK _® N	Master Classical Standard S	et, Powder 12g
Order No.	V70-001	\$1,550.00
1	Wash Opaque, WO	
16	Opaque, A1-D4	
16	Opaque Dentine, A1-D4	
16	Dentine, A1-D4	
2	Enamel, EN1, EN2	
1	Translucent, T4	
1	Window, WIN	
3	Cervical CE1-CE3	
5	Neck, N1-N5	
2	Corrective, COR1, COR2	
1	VITA Modelling Fluid, 50ml	
1	VITA Modelling Fluid RS, 50ml	
1	VITA Opaque Fluid, 50ml	
1	Shade indicator VITAPAN classical	
1	VITAPAN Classical shade guide	
_	Accessories	

VMK®Master Classical Standard Set, Paste 12g

Order No.	V70-002	\$2,050.00
1	Wash Opaque, WO	
16	Opaque, A1-D4	
16	Opaque Dentine, A1-D4	537
16	Dentine, A1-D4	200
2	Enamel, EN1, EN2	6588V.
1	Translucent, T4	
1	Window, WIN	-4
3	Cervical, CE1-CE3	
5	Nock N1 NE	

Neck, N1-N5 Corrective, COR1, COR2 VITA Modelling Fluid, 50ml VITA Modelling Fluid RS, 50ml

VITA PASTE Opaque Liquid, 15ml

Shade indicator VITAPAN classical VITAPAN Classical shade guide

Accessories

VMK®Master Classical Starter Set, Powder 12g Order No. V70-005..

2 Opaque, A2, A3 Opaque Dentine, A2, A3

Dentine, A2, A3

Enamel, EN1

Translucent, T4

VITA Modelling Fluid RS, 50ml

VITA Opaque Fluid, 50ml

Shade tabs A2, A3



St

Starter S	aster Classical Set, Paste 12g	
2 2 2 1 1 1	70-006	\$139.00
VMK _® Ma	aster Additional Set, 12g	
8 8 6 5 1 1 1 1 1	Po-009	\$999.00
VMK Master	Classical Shade Indicator (Fold-out)	¢70.00
VMK Master	os-940 3D-Master Shade Indicator (Fold-out) (68-941	
Order No. G6	Color Indicator, Translucent 68-942	\$72.95
Order No. G6	Color Indicator, Margin/Gingiva	\$72.95
Order No. G6	Color Indicator, Dentine Modif/Luminary 68-944 Color Indicator, Standard	\$72.95
	88-945	\$72.95

VMK®Master 3D-Master Reorder Materials

Order No. V70 VMK Master V Order No. V71 VMK Master V	Wash Opaque Pow -259 Wash Opaque Past	der, 50g	75.95
VMK Master (Opaque Powder, 12	'g	\$21.95
Order No.	Shade	Order No.	Shade
VX70-250	OP0	VX70-253	OP3
VX70-251	OP1	VX70-254	OP4
VX70-252	OP2	VX70-255	OP5
VMK Master (Opaque Powder, 50	lg	\$75.95
Order No.	Shade	Order No.	Shade
VX71-250	OP0	VX71-253	OP3
VX71-251	OP1	VX71-254	OP4
VX71-252	OP2	VX71-255	OP5

Technical Help Desk: 888-249-1640 USA | 800-324-6224 Canada

	jue Pasie, 5g		
Order No.		Order No.	Shade
VX70-350	OP0	VX70-353	OP3
VX70-351	OP1	VX70-354	OP4
VX70-352	OP2	VX70-355	OP5
VMV Moster Open	uus Dontino 12a		¢ 21 0E
Order No	Shade		
		Order No.	Shade
VX70-061	0M1	VX70-076	3M2
VX70-062	OM2	VX70-077	3M3
VX70-063	0M3	VX70-078	3R1.5
VX70-064	1M1	VX70-079	3R2.5
VX70-065	1M2	VX70-080	4L1.5
VX70-066	2L1.5	VX70-081	4L2.5
VX70-067	2L2.5	VX70-082	4M1
VX70-068	2M1	VX70-083	4M2
VX70-069	2M2	VX70-084	4M3
VX70-070	2M3	VX70-085	4R1.5
VX70-071	2R1.5	VX70-086	4R2.5
VX70-072	2R2.5	VX70-087	5M1
VX70-073	3L1.5	VX70-088	5M2
VX70-074	3L2.5	VX70-089	5M3
VX70-075	3M1		
VMK Master Opac	ue Dentine 50a		\$75 95
Order No.	Shade		Shade
VX71-061	0M1	VX71-076	3M2
VX71-062	OM2	VX71-077	3M3
VX71-063	0M3	VX71-078	3R1.5
VX71-064	1M1	VX71-079	3R2.5
VX71-065	1M2	VX71-080	4L1.5
VX71-066	2L1.5	VX71-081	4L2.5
VX71-067	2L2.5	VX71-081	4M1
VX71-007 VX71-068	2M1	VX71-082 VX71-083	4M2
VX71-068 VX71-069	2M2	VX71-063 VX71-084	4M3
VX71-009 VX71-070	2M3	VX71-004 VX71-085	4R1.5
VX71-070 VX71-071	2R1.5	VX71-065 VX71-086	4R1.5 4R2.5
VX71-071 VX71-072	2R2.5	VX71-080 VX71-087	5M1
VX71-072 VX71-073		VX71-087 VX71-088	5M2
	3L1.5		
VX71-074 VX71-075	3L2.5 3M1	VX71-089	5M3
		O-d N-	
Order No.	Shade	Order No.	Shade
VX70-031	0M1	VX70-046	3M2
VX70-032	OM2	VX70-047	3M3
VX70-033	OM3	VX70-048	3R1.5
VX70-034	1M1	VX70-049	3R2.5
VX70-035	1M2	VX70-050	4L1.5
VX70-036	2L1.5	VX70-051	4L2.5
VX70-037	2L2.5	VX70-052	4M1
VX70-038	2M1	VX70-053	4M2
VX70-039	2M2	VX70-054	4M3
VX70-040	2M3	VX70-055	4R1.5
VX70-041	2R1.5	VX70-056	4R2.5
VX70-042	2R2.5	VX70-057	5M1
VX70-043	3L1.5	VX70-058	5M2
VX70-044	3L2.5	VX70-059	5M3
VX70-045	3M1		

VMK Master Opaque Paste, 5g.



)entine, 50g		\$75.95	VMK Master (Opaque Paste, 5g		\$69.9
Shade	Order No.	Shade	Order No.	Shade	Order No.	Shade
0M1	VX71-046	3M2	V70-301	A1	V70-309	B4
0M2	VX71-047	3M3	V70-302	A2	V70-310	C1
						C2
						C3
						C4
						D2
						D2
			V7U-3U8	ВЗ	V/U-316	D4
			\	D 1 4	•	A 04 6
						Shade
						B4
			V70-432	A2	V70-440	C1
3L2.5	VX71-059	5M3	V70-433	A3	V70-441	C2
3M1			V70-434	A3.5	V70-442	C3
			V70-435	A4	V70-443	C4
namel, 12g		\$21.95	V70-436	B1	V70-444	D2
						D3
						D4
	V/0 133	LIVO	V/0 400	БО	V/O ++O	DŦ
LINZ			VMK Master (Inaque Dentine 50	Πα	\$75.0
namel 50a		\$75.95			-	Shade
•						B4
						C1
	V/1-133	LINO				
EINZ						C2
						C3
ster Classica	ı l					C4
/lotoriolo						D2
nateriais						D3
Opaque Powder, 12	'a	\$21.95	V71-438	B3	V71-446	D4
			VMK Master [•		\$21.9
			Order No.	Shade	Order No.	Shade
			V70-451	A1	V70-459	B4
			V70-452	A2	V70-460	C1
			V70-453	A3	V70-461	C2
				A3.5	V70-462	C3
						C4
						D2
B3	V70-416	D4				D3
						D3 D4
)paque Powder, 50		\$75.95	v / U-430	טט	v / U-400	υ 4
Shade	Order No.	Shade	\/\ /I/ \\ / ==±=== F	Contine FOS		Φ7 Γ (
A1	V71-409					
A2	V71-410	C1				Shade
A3	V71-411	C2	V71-451	A1	V71-459	B4
A3.5	V71-412	C3	V71-452	A2	V71-460	C1
		C4	V71-453	A3	V71-461	C2
	\/ /1_/11/2	UT	1/74 454	A3.5	V71-462	C3
A4	V71-413		V71-454	A3.3	V/1-40Z	U3
A4 B1	V71-414	D2			V71-462 V71-463	
A4 B1 B2	V71-414 V71-415	D2 D3	V71-455	A4	V71-463	C4
A4 B1	V71-414	D2				
	Shade 0M1 0M2 0M3 1M1 1M2 2L1.5 2L2.5 2M1 2M2 2M3 2R1.5 2R2.5 3L1.5 3L2.5 3M1 Inamel, 12g Shade EN1 EN2 Inamel, 50g Shade EN1 EN2 Ster Classica Materials Dpaque Powder, 12 Shade A1 A2 A3 A3.5 A4 B1 B2 B3 Dpaque Powder, 50 Shade A1	Shade Order No. 0M1 VX71-046 0M2 VX71-047 0M3 VX71-048 1M1 VX71-049 1M2 VX71-050 2L1.5 VX71-051 2L2.5 VX71-052 2M1 VX71-053 2M2 VX71-054 2M3 VX71-055 2R1.5 VX71-055 2R2.5 VX71-056 2R2.5 VX71-057 3L1.5 VX71-058 3L2.5 VX71-059 3M1 V70-193 EN2 Shade Order No. EN1 V70-193 EN2 Shade Order No. EN1 V71-193 EN2 Shade Order No. A1 V70-409 A2 V70-410 A3 V70-411 A3.5 V70-412 A4 V70-413 B1 V70-415 B3 V70-416 Opaque Powder, 50g Shade Order No. A1 V70-415 B3 V70-416 Opaque Powder, 50g Shade Order No. A1 V70-415 B3 V70-416	0M1 VX71-046 3M2 0M2 VX71-047 3M3 0M3 VX71-048 3R1.5 1M1 VX71-049 3R2.5 1M2 VX71-050 4L1.5 2L1.5 VX71-051 4L2.5 2L2.5 VX71-052 4M1 2M1 VX71-053 4M2 2M2 VX71-054 4M3 2M3 VX71-055 4R1.5 2R1.5 VX71-056 4R2.5 2R2.5 VX71-057 5M1 3L1.5 VX71-058 5M2 3L2.5 VX71-059 5M3 3M1 ** ** *** ** ** *** ** ** *** ** ** *** ** ** *** ** ** *** ** ** *** ** ** *** ** ** *** ** **	Shade Order No. Shade Order No. 0M1 VX71-046 3M2 V70-301 0M2 VX71-047 3M3 V70-302 0M3 VX71-048 3R1.5 V70-303 1M1 VX71-050 41.5 V70-305 1M2 VX71-051 4L2.5 V70-306 2L1.5 VX71-051 4L2.5 V70-306 2L2.5 VX71-052 4M1 V70-307 2M1 VX71-053 4M2 V70-308 2M2 VX71-054 4M3 WK 2M3 VX71-055 4R1.5 VMK Master Offer No. 2R1.5 VX71-056 4R2.5 Order No. 2R2.5 VX71-057 5M1 V70-431 3L1.5 VX71-058 5M2 V70-432 3L2.5 VX71-059 5M3 V70-433 3M1 V70-432 V70-435 3M1 V70-435 V70-436 Shade Order No. Shade V70-437 EN1 <td> Shade</td> <td> Shade</td>	Shade	Shade

Order No. V70-121 V70-122 V70-123 V70-124 VMK Master 0	Dentine Modifier, 12g Shade DM1 DM2 DM3 DM4 Opal Translucent (OT1),		Shade DM5 DM6 DM7 DM8
VMK Master F	-171 Pearl Translucent (PLT1) -151), 12g	
	ranslucent, 12g		
Order No.	Shade	Order No.	Shade
V70-091	T1	V70-095	T5
V70-092	T2	V70-096	T6
V70-093	T3	V70-097	T7
V70-094	T4	V70-098	T8
VMK Master T	ranslucent, 50g		
	094		\$75.95
	Vindow (WIN), 12g		
	-181		\$21.95
	Vindow (WIN), 50g		
Order No. V71	-181		\$75.95
\/h 41/ h 4			404.05
	Cervical, 12g		
	Shade	Order No.	Shade
V70-201	CE1	V70-203	CE3
V70-202	CE2		
VMK Mastar N	Neck, 12g		¢21 Q5
Order No.	Shade	Order No.	Shade
V70-241	N1	V70-244	N4
V70-242	N2	V70-245	N5
V70-243	N3		
VMK Master (Corrective, 12g		\$24.95
Order No.	Shade	Order No.	Shade
V70-221	COR1	V70-222	COR2
170 221	00111	170 222	00112
VMK Master I	uminary, 12g		\$49 00
Order No.	Shade	Order No.	Shade
V70-141	LM1	V70-144	LM4
V70-141 V70-142	LM2	V70-144 V70-145	LM5
V70-143	LM3	V70-146	LM6
VMK Master (Gingiva, 12g		\$24.95
Order No.	Shade	Order No.	Shade
V70-231	G1	V70-234	G4
V70-232	G2	V70-235	G5
V70-232	G3	¥70 200	GO
V/U 200	uu		

Technical Help Desk: 888-249-1640 USA | 800-324-6224 Canada

VMK Master Mary Order No. V70-260 V70-261 V70-262	gin, 12g Shade MN M1 M2	\$24.95 Order No. V70-263 V70-264 V70-265	Shade M3 M4 M5
VMK Master Ging Order No. V70-271	jiva Opaque, 1 Shade GOL	2g Order No. V70-273	\$21.95 Shade GOD
VMK Master Ging Order No. V70-371 VITA Modeling Flu	Shade GOL	aste, 5g Order No. V70-373	\$ 69.95 Shade GOD
	0		\$ 18.60
0	0		\$ 68.00
Order No. V68-91:	3		\$ 24.70
	3		\$ 79.50
	1		\$ 18.60
	1		\$ 68.00
VITA Paste Opaqu Order No. V68-92			\$ 29.90

Titanium Porcelain

VITA Titanium Porcelain is matched to the specific material properties of titanium. The inherent versatile properties of titanium provide an incredibly biocompatible, corrosionresistant coping that is lightweight, hypo-allergenic, mechanically strong (Vicker's hardness of 180-250, degree 1-4) and has low thermal conductivity. VITA Titanium Porcelain provides the proper coefficient of thermal expansion for titanium copings (9.6 x 10-6 K-1) and is



available in both VITA 3D-Master and Classical shades.

Titanium Kits

METAL CERAMICS

Titanium 3D-Master Starter Kit (3M2)	
Order No. VX16-001	\$199.00
Titanium 3D-Master Starter Kit (2M2)	
Order No. VX16-003	\$199.00
Titanium 3D-Master Standard Set, 12g	
Order No. VX17-002	\$1,899.00
Titanium Classical Standard Set, 12g	
Order No. V17-002	\$1,495.00
Master Kit*	
Order No. V17-003	\$TBA
A2 Starter Kit*	
Order No. V17-004	\$TBA
A3 Starter Kit*	
Order No. V17-005	\$TBA

Titanium Materials

16

	um Opaque Powde	rs	\$32.00
50g*			\$100.00
	50g Order No.		
VX17-2501	VX18-2501	OP1	
VX17-2502	VX18-2502	OP2	
VX17-2503	VX18-2503	0P3	
VX17-2504	VX18-2504	OP4	
VX17-2505	VX18-2505	OP5*	
3D-Master Titani	um Opaque Dentin	e Powder, 12g	\$32.00
Order No.	Shade	Order No.	Shade
VX17-2331	1M1	VX17-2344	3M3
VX17-2332	1M2	VX17-2345	3R1.5
VX17-2333	2L1.5	VX17-2346	3R2.5
VX17-2334	2L2.5	VX17-2347	4L1.5
VX17-2335	2M1	VX17-2348	4L2.5
VX17-2336	2M2	VX17-2349	4M1
VX17-2337	2M3	VX17-2350	4M2
VX17-2338	2R1.5	VX17-2351	4M3
VX17-2339	2R2.5	VX17-2352	4R1.5
VX17-2340	3L1.5	VX17-2353	4R2.5
VX17-2341	3L2.5	VX17-2354	5M1

VX17-2342	3M1	VX17-2355	5M2
VX17-2343	3M2	VX17-2356	5M3
3D-Master Titaniu	m Dentine Powders		
12g			\$32.00
12g Order No.	50g Order No.	Shade	
VX17-2361	VX18-2361	1M1	
VX17-2362	VX18-2362	1M2	
VX17-2363	_	2L1.5	
VX17-2364	VX18-2364*	2L2.5	
VX17-2365	VX18-2365	2M1	
VX17-2366	VX18-2366	2M2	
VX17-2367	VX18-2367	2M3	
VX17-2368	_	2R1.5	
VX17-2369	VX18-2369*	2R2.5	
VX17-2370	VX18-2370	3L1.5	
VX17-2371	VX18-2371	3L2.5	
VX17-2372	VX18-2372	3M1	
VX17-2373	VX18-2373	3M2	
VX17-2374	VX18-2374	3M3	
VX17-2375	VX18-2375	3R1.5	
VX17-2376	VX18-2376	3R2.5	
VX17-2377	_	4L1.5	
VX17-2378	_	4L2.5	
VX17-2379	VX18-2379	4M1	
VX17-2380	VX18-2380	4M2	
VX17-2381	VX18-2381	4M3	
VX17-2382	_	4R1.5	
VX17-2383	_	4R2.5	
VX17-2384	_	5M1	
VX17-2385	_	5M2	
VX17-2386	_	5M3	
Titanium Classical	/3D Enamel, 12g		\$32.00
Order No.	Shade		
VX17-2071	EN1*		
VX17-2072	EN2		
VX17-2073	EN3*		
VX17-2074	EN4		

* While	supp	lies	last
---------	------	------	------

			ტა2.00
•	50g Order No.		\$100.00
2g Order No.			
/17-2001	V18-2001	A1	
/17-2002	V18-2002	A2	
/17-2003	V18-2003	A3	
/17-2004	V18-2004	A3.5	
/17-2005	V18-2005*	A4	
/17-2006	V18-2006*	B1	
′17-2007	V18-2007*	B2	
/17-2008	V18-2008	B3	
/17-2009	V18-2009*	B4	
/17-2010	V18-2010*	C1	
/17-2011	V18-2011*	C2	
/17-2012	V18-2012*	C3	
/17-2013	V18-2013*	C4	
/17-2014	V18-2014*	D2	
/17-2015	V18-2015*	D3	
/17-2016	V18-2016*	D4	
	2010	2.	
	ım Dentine Powders		***
	FO OI N		\$100.00
2g Order No.	-		
/17-2051	V18-2051	A1	
/17-2052	V18-2052	A2	
17-2053	V18-2053	A3	
′17-2054	V18-2054	A3.5	
/17-2055	V18-2055*	A4	
/17-2056	V18-2056*	B1	
/17-2057	_	B2	
	V18-2058	B3	
/17-2058		B4	
/17-2058 /17-2059	V18-2059*		
17-2059		C1	
/17-2059 /17-2060	V18-2060*	C1 C2	
/17-2059 /17-2060 /17-2061	V18-2060* V18-2061*	C2	
/17-2059 /17-2060 /17-2061 /17-2062	V18-2060* V18-2061* V18-2062*	C2 C3	
/17-2059 /17-2060 /17-2061 /17-2062 /17-2063	V18-2060* V18-2061* V18-2062* V18-2063*	C2 C3 C4	
/17-2059 /17-2060 /17-2061 /17-2062 /17-2063 /17-2064	V18-2060* V18-2061* V18-2062* V18-2063* V18-2064*	C2 C3 C4 D2	
/17-2059 /17-2060 /17-2061 /17-2062	V18-2060* V18-2061* V18-2062* V18-2063*	C2 C3 C4	
'17-2059 '17-2060 '17-2061 '17-2062 '17-2063 '17-2064 '17-2065 '17-2066	V18-2060* V18-2061* V18-2062* V18-2063* V18-2064* V18-2065 V18-2066*	C2 C3 C4 D2 D3 D4	\$32.00
/17-2059 /17-2060 /17-2061 /17-2062 /17-2063 /17-2064 /17-2065 /17-2066	V18-2060* V18-2061* V18-2062* V18-2063* V18-2064* V18-2065 V18-2066*	C2 C3 C4 D2 D3 D4	\$32.00 Shade
/17-2059 /17-2060 /17-2061 /17-2062 /17-2063 /17-2064 /17-2065 /17-2066 Classical Titaniu Order No.	V18-2060* V18-2061* V18-2062* V18-2063* V18-2064* V18-2065 V18-2066* Im Opaque Dentine, Shade	C2 C3 C4 D2 D3 D4	Shade
/17-2059 /17-2060 /17-2061 /17-2062 /17-2063 /17-2064 /17-2065 /17-2066 Classical Titaniu Order No. /17-2031	V18-2060* V18-2061* V18-2062* V18-2063* V18-2064* V18-2065 V18-2066* Im Opaque Dentine, Shade A1	C2 C3 C4 D2 D3 D4 12g Order No. V17-2039	Shade B4
/17-2059 /17-2060 /17-2061 /17-2062 /17-2063 /17-2064 /17-2065 /17-2066 Classical Titaniu Order No. /17-2031 /17-2032	V18-2060* V18-2061* V18-2062* V18-2063* V18-2064* V18-2065 V18-2066* Im Opaque Dentine, Shade A1 A2	C2 C3 C4 D2 D3 D4 12g Order No. V17-2039 V17-2040	Shade B4 C1
717-2059 717-2060 717-2061 717-2062 717-2063 717-2065 717-2066 8 Classical Titaniu 9 Order No. 717-2031 717-2032 717-2033	V18-2060* V18-2061* V18-2062* V18-2063* V18-2064* V18-2065 V18-2066* Im Opaque Dentine, Shade A1 A2 A3	C2 C3 C4 D2 D3 D4 12g Order No. V17-2039 V17-2040 V17-2041	Shade B4 C1 C2
717-2059 717-2060 717-2061 717-2062 717-2063 717-2065 717-2066 8 Classical Titaniu 9 Order No. 717-2031 717-2032 717-2034	V18-2060* V18-2061* V18-2062* V18-2063* V18-2064* V18-2065 V18-2066* Im Opaque Dentine, Shade A1 A2 A3 A3.5	C2 C3 C4 D2 D3 D4 12g Order No. V17-2039 V17-2040 V17-2041 V17-2042	Shade B4 C1 C2 C3
17-2059 17-2060 17-2061 17-2062 17-2063 17-2065 17-2066 18assical Titaniu 17-2031 17-2032 17-2033 17-2034 17-2035	V18-2060* V18-2061* V18-2062* V18-2063* V18-2064* V18-2065 V18-2066* Im Opaque Dentine, Shade A1 A2 A3 A3.5 A4	C2 C3 C4 D2 D3 D4 12g Order No. V17-2039 V17-2040 V17-2041 V17-2042 V17-2043	Shade B4 C1 C2 C3 C4
717-2059 717-2060 717-2061 717-2062 717-2063 717-2065 717-2066 8 Classical Titaniu 9 Order No. 717-2031 717-2032 717-2034 717-2034 717-2036	V18-2060* V18-2061* V18-2062* V18-2063* V18-2065 V18-2066* Im Opaque Dentine, Shade A1 A2 A3 A3.5 A4 B1	C2 C3 C4 D2 D3 D4 12g Order No. V17-2039 V17-2040 V17-2041 V17-2042 V17-2043 V17-2044	Shade B4 C1 C2 C3 C4 D2
/17-2059 /17-2060 /17-2061 /17-2062 /17-2063 /17-2065 /17-2066 Classical Titaniu Order No. /17-2031 /17-2032 /17-2034 /17-2035 /17-2036 /17-2036 /17-2037	V18-2060* V18-2061* V18-2062* V18-2063* V18-2065 V18-2066* Im Opaque Dentine, Shade A1 A2 A3 A3.5 A4 B1 B2	C2 C3 C4 D2 D3 D4 12g Order No. V17-2039 V17-2040 V17-2041 V17-2042 V17-2043 V17-2044 V17-2045	Shade B4 C1 C2 C3 C4 D2 D3
/17-2059 /17-2060 /17-2061 /17-2062 /17-2063 /17-2065 /17-2066 Classical Titaniu Order No. /17-2031 /17-2032 /17-2034 /17-2035 /17-2036 /17-2036 /17-2037	V18-2060* V18-2061* V18-2062* V18-2063* V18-2065 V18-2066* Im Opaque Dentine, Shade A1 A2 A3 A3.5 A4 B1	C2 C3 C4 D2 D3 D4 12g Order No. V17-2039 V17-2040 V17-2041 V17-2042 V17-2043 V17-2044	Shade B4 C1 C2 C3 C4 D2
717-2059 717-2060 717-2061 717-2062 717-2063 717-2064 717-2065 717-2066 61assical Titaniu 717-2031 717-2032 717-2034 717-2035 717-2036 717-2037 717-2038	V18-2060* V18-2061* V18-2062* V18-2063* V18-2065 V18-2066* Im Opaque Dentine, Shade A1 A2 A3 A3.5 A4 B1 B2 B3	C2 C3 C4 D2 D3 D4 12g Order No. V17-2039 V17-2040 V17-2041 V17-2042 V17-2043 V17-2044 V17-2045 V17-2046	Shade
717-2059 717-2060 717-2061 717-2062 717-2063 717-2064 717-2065 717-2066 717-2031 717-2032 717-2034 717-2035 717-2037 717-2038 717-2038 717-2038 717-2038 717-2038	V18-2060* V18-2061* V18-2062* V18-2063* V18-2065 V18-2066* Im Opaque Dentine, Shade A1 A2 A3 A3.5 A4 B1 B2 B3 Opaque, 12g	C2 C3 C4 D2 D3 D4 12g Order No. V17-2039 V17-2040 V17-2041 V17-2042 V17-2043 V17-2044 V17-2045 V17-2046 Order No.	Shade
717-2059 717-2060 717-2061 717-2062 717-2063 717-2064 717-2065 717-2066 Classical Titaniu 717-2031 717-2032 717-2034 717-2035 717-2036 717-2038 Titanium Color Colorder No. 717-2021	V18-2060* V18-2061* V18-2062* V18-2063* V18-2064* V18-2066* Im Opaque Dentine, Shade A1 A2 A3 A3.5 A4 B1 B2 B3 Opaque, 12g	C2 C3 C4 D2 D3 D4 12g Order No. V17-2039 V17-2040 V17-2041 V17-2042 V17-2043 V17-2045 V17-2046 Order No. V17-2025	Shade
/17-2059 /17-2060 /17-2061 /17-2062 /17-2063 /17-2064 /17-2065 /17-2066 Classical Titaniu Order No. /17-2031 /17-2032 /17-2034 /17-2035 /17-2036 /17-2037 /17-2038	V18-2060* V18-2061* V18-2062* V18-2063* V18-2065 V18-2066* Im Opaque Dentine, Shade A1 A2 A3 A3.5 A4 B1 B2 B3 Opaque, 12g	C2 C3 C4 D2 D3 D4 12g Order No. V17-2039 V17-2040 V17-2041 V17-2042 V17-2043 V17-2044 V17-2045 V17-2046 Order No.	Shade
717-2059 717-2060 717-2061 717-2062 717-2063 717-2064 717-2065 717-2066 Classical Titaniu 717-2031 717-2032 717-2034 717-2035 717-2036 717-2038 Titanium Color Colorder No. 717-2021	V18-2060* V18-2061* V18-2062* V18-2063* V18-2064* V18-2066* Im Opaque Dentine, Shade A1 A2 A3 A3.5 A4 B1 B2 B3 Opaque, 12g	C2 C3 C4 D2 D3 D4 12g Order No. V17-2039 V17-2040 V17-2041 V17-2042 V17-2043 V17-2045 V17-2046 Order No. V17-2025	Shade

Order No. V18-2071	Shade EN1		
V18-2073	EN3		
	nslucent, 12g		\$
Order No.	Shade		
V17-2081	T1 Milky		
V17-2082	T2 Opal		
V17-2083	T3 Pearl		
V17-2084	T4 Neutral		
V17-2085	T5 Polar		
V17-2086	T6 Blue		
	nslucent (T4 Shade), 50g		0.4
Order No. V1	8-2084 Clear		\$1
Titanium Wir	ndow, 12g (17-2122		\$
	vical, 12g		\$
Order No.	Shade		
V17-2091	CE1 Light Orange		
V17-2092	CE2 Light Yellow		
V17-2093	CE3 Gold		
V17-2094	CE4 Clear		
	ensive, 12g		
Order No.	Shade	Order No.	Shac
V17-2111	IN1	V17-2116	IN6
V17-2112	IN2	V17-2117	IN7
V17-2113	IN3	V17-2118	IN8
V17-2114	IN4	V17-2119	IN9
V17-2115	IN5	V17-2120	IN10
Titanium Gin	giva, 12g		\$
Order No.	Shade		
V17-2141	G1 Light		
V17-2142	G2 Dark		
Titanium Ma	rgin Porcelain, 12g		\$
Order No.	Shade	Order No.	Shac
V17-2600	MAR N	V17-2604	MAF
V17-2601	MAR 1	V17-2605	MAF
V17-2602	MAR 2	V17-2606	MAF
V17-2603	MAR 3		
Titanium Bor			
Order No.	Description		_
V15-100	Paste Bonder, 7g		
V15-101	Powder Bonder, 7.5g		
V16-101	Powder Bonder, 36g .		\$
Titanium Sha			6 4
Lirder No. (3)	(17-906		\$1

* While supplies last



VITA VM®9

ALL-CERAMICS AND PRESSABLES

VITA VM9 features desirable characteristics while satisfying the unique requirements of yttrium-stabilized zirconia. This 3D-Master and Classical shaded porcelain is compatible with VITA PM9 pressable ceramics, VITA's In-Ceram YZ Cubes for CEREC and many of today's zirconia substructure materials.

VM_®9 Basic Kit, 12g

Order No.	NX98-0101	\$1,450.00
3	VITA VM9 Chroma Plus	
26	VITA VM9 Base Dentine	1
2	VITA VM9 Enamel	
1	VITA VM9 Neutral	The state of the s
1	VITA VM9 Window	1 0000
3	VITA VM9 Corrective	L. Car

1 VITA VM Modeling Liquid, 50ml

Modisol Pen

1 Firing Tray G

1 Fibrous pad firing support

1 VITA VM9 shade indicator

1 3D-Master shade guide

1 Working instructions

3 Sun Dentine SD1-SD3

VM_®9 Basic Kit, 12g (Small)

1 Working instructions

VITA VM9 shade indicator

3D-Master shade guide

VM®9 Buildup Kit, 12g



VM®9 Buildup Kit, 12g (Small)

VITA VM Modeling Liquid, 50ml

Order No. I	NX98-0006	\$436.00
15	VITA VM9 Transpa Dentine	

VM_®9 Professional Kit, 12g

VM®9 Professional Kit, 12g (Small)

Order No.	NX98-0014.		\$550.00
5	VITA VM9	Effect Chroma	
7	VITA VM9	Effect Enamel	
1	VITA VM9	Mamelon	
1	VITA VM9	Effect Pearl	
1	VITA VM9	Effect Opal	
1	VITA VM9	Shade Indicators	

VM_®9 Gingiva Kit, 12g

Order No.	NX98-0015\$226.00	
5	VITA VM9 Gingiva	
1	VITA VM9 Shade Indicator	

VM®9 Bleached Shade Kit, 12g

Order No.	NX98-0007	\$338.00
3	VITA VM9 Base Dentine	
3	VITA VM9 Transpa Dentine	
1	VITA VM9 Enamel ENL (light)	
1	VITA VM9 Neutral	100
1	VITA VM9 Window	
1	VITA VM Modeling Liquid, 50ml	
1	3D-Master bleached shade guide	
1	Working instructions	

VM_®9 2M2 One Shade Kit, 12g

Order No.	NX98-0011	\$133.00
1	VITA VM9 Base Dentine	
1	VITA VM9 Transpa Dentine	
1	VITA VM9 Enamel ENL (light)	
1	VITA VM9 Neutral	4
1	VITA VM9 Window	Tallar)
1	VITA VM Modeling Liquid, 50ml	
1	Working instructions	

VM®9 Margin Kit, 12g

Order No. NX	98-0020	\$254.00
1 N	/11 White	
1 N	14 Yellow	
1 N	15 Amber	
1 N	И7 Light Beige	
1 N	18 Pastel Brown	
1 N	19 Light Orange	

VM®9 Effect Opal Kit

Working instructions

der No.	NX98-0021 FPO	
5	Effect Opal	
1	Shade indicator	

VITA VM®9 Low Fusing Add-On Powders

The low-fusing (780° C) VITA VM9 Add-On materials are particularly well suited for individualizing (layering technique) substructure-free VITA PM9 restorations. This material is based on the proven VITA VM9 veneering material, which enables the technician to reproduce shade nuances in a fast manner to achieve top-quality, esthetic restorations.

VITA VM 9 Add-On:

- Economic efficiency: The excellent stability during layering, minimal shrinkage and outstanding shade reproduction ensure economic processing.
- Natural shade effect: Thanks to its beautiful translucency, enamel-like refraction behavior, efficient polishability and excellent bond to the VITA PM9 base structure, VITA VM 9 Add-On materials create a natural interplay of colors.
- Reliable bonding (adhesion): VITA VM 9 Add-On materials have been perfectly matched with the VITA PM9 ceramic and easily achieve excellent, distortion-free sintered bond to the VITA PM9 base structure. As a result, permanent loadability and strength is ensured.

VM_®9 Add-On Kit (for VITA PM_®9)

Order No.	N98-0035	\$329.00
1	VM9 Add-On ADD1 (Transparent), 12g	
1	VM9 Add-On ADD2 (Light Enamel), 12g	
1	VM9 Add-On ADD3 (Dark Enamel), 12g	
1	VM9 Add-On ADD4 (White Transparent), 12g	
1	VM9 Add-On ADD5 (Yellow Transparent), 12g	
1	VM9 Add-On ADD6 (Orange Translucent), 12g	
1	VM9 Add-On ADD7 (Red Translucent), 12g	
1	VM9 Add-On ADD8 (Blue Translucent), 12g	
1	VITA VM9 Glaze LT, 12g	
1	VITA VM Modeling Liquid, 50ml	
1	VITA AKZENT Fluid, 20ml	
1	Fibrous pad firing support	
1	A&E Brush, No. 3/0	
1	VITA VM9 Add-On Shade Indicator	
1	Working Instructions VITA PM9, 1450	

VM_®9 Add-On Powders

VM9 Add-On Po	wders	. \$34.00
Order No.	Description	
NX98-2291	VM9 Add-On ADD1 (Transparent)	
NX98-2292	VM9 Add-On ADD2 (Light Enamel)	
NX98-2293	VM9 Add-On ADD3 (Dark Enamel)	
NX98-2294	VM9 Add-On ADD4 (White Transparent)	
NX98-2295	VM9 Add-On ADD5 (Yellow Transparent)	
NX98-2296	VM9 Add-On ADD6 (Orange Transluscent)	
NX98-2297	VM9 Add-On ADD7 (Red Transluscent)	
NX98-2298	VM9 Add-On ADD8 (Blue Transluscent)	

Technical Help Desk: 888-249-1640 USA | 800-324-6224 Canada

VM9 3D-Master Reorder Materials

\/\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\	T D	40	004.00
	er Transpa Denti	ne, 12g	
Order No.	Shade	Order No.	Shade
NX98-2061	0M1	NX98-2076	3M2
NX98-2062	0M2	NX98-2077	3M3
NX98-2063	0M3	NX98-2078	3R1.5
NX98-2064	1M1	NX98-2079	3R2.5
NX98-2065	1M2	NX98-2080	4L1.5
NX98-2066	2L1.5	NX98-2081	4L2.5
NX98-2067	2L2.5	NX98-2082	4M1
NX98-2068	2M1	NX98-2083	4M2
NX98-2069	2M2	NX98-2084	4M3
NX98-2070	2M3	NX98-2085	4R1.5
NX98-2071	2R1.5	NX98-2086	4R2.5
NX98-2072	2R2.5	NX98-2087	5M1
NX98-2073	3L1.5	NX98-2088	5M2
NX98-2074	3L2.5	NX98-2089	5M3
NX98-2075	3M1	14/100 2000	OIVIO
14/130 2073	JIVII		
VM9 3D-Maste	er Base Dentine	, 12g	\$34 00
Order No.	Shade	Order No.	Shade
NX98-2031	0M1	NX98-2046	3M2
NX98-2032	0M2	NX98-2047	3M3
NX98-2033	0M3	NX98-2048	3R1.5
NX98-2034	1M1	NX98-2049	3R2.5
NX98-2035	1M2	NX98-2050	4L1.5
NX98-2036	2L1.5	NX98-2051	4L1.5
NX98-2037	2L1.5 2L2.5	NX98-2052	4L2.5 4M1
NX98-2038	2M1	NX98-2053	4M2
NX98-2039	2M2	NX98-2054	4M3
NX98-2040	2M3	NX98-2055	4R1.5
NX98-2041	2R1.5	NX98-2056	4R2.5
NX98-2042	2R2.5	NX98-2057	5M1
NX98-2043	3L1.5	NX98-2058	5M2
NX98-2044	3L2.5	NX98-2059	5M3
NX98-2045	3M1		
\/\ 40 T			# 400.00
	lentine, 50g	0 N	
Order No.	Shade	Order No.	Shade
NX99-2061	0M1	NX99-2074	3L2.5
NX99-2062	0M2	NX99-2075	3M1
NX99-2063	0M3	NX99-2076	3M2
NX99-2064	1M1	NX99-2077	3M3
NX99-2065	1M2	NX99-2078	3R1.5
NX99-2066	2L1.5	NX99-2079	3R2.5
NX99-2067	2L2.5	NX99-2082	4M1
NX99-2068	2M1	NX99-2083	4M2
NX99-2069	2M2	NX99-2084	4M3
NX99-2070	2M3	NX99-2087	5M1
NX99-2071	2R1.5	NX99-2088	5M2
NX99-2072	2R2.5	NX99-2089	5M3
NX99-2073	3L1.5		



	ne, 50g				r Effect Pearl, 12g		\$38.20
Order No.	Shade	Order No.	Shade	Order No.	Description		
NX99-2031	0M1	NX99-2044	3L2.5	NX98-2151	EP1		
NX99-2032	0M2	NX99-2045	3M1	NX98-2152	EP2		
NX99-2033	0M3	NX99-2046	3M2	NX98-2153	EP3		
NX99-2034	1M1	NX99-2047	3M3	10/100 2100	LIO		
				\/\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\\	40		400.00
NX99-2035	1M2	NX99-2048	3R1.5		ne, 12g		\$38.20
NX99-2036	2L1.5	NX99-2049	3R2.5	Order No.	Shade		
NX99-2037	2L2.5	NX99-2052	4M1	NX98-2361	Sunlight		
NX99-2038	2M1	NX99-2053	4M2	NX98-2362	Sunrise		
NX99-2039	2M2	NX99-2054	4M3	NX98-2363	Sunset		
NX99-2040	2M3	NX99-2057	5M1	14/100 2000	Odrisot		
				\/MO OD M+-	M		ΦΩΩ ΩΩ
NX99-2041	2R1.5	NX99-2058	5M2		r Mamelon, 12g		\$38.20
NX99-2042	2R2.5	NX99-2059	5M3	Order No.	Description		
NX99-2043	3L1.5			NX98-2161	MM1	NX98-2174	E05
				NX98-2162	MM2	NX98-2175	E05
VM9 Fnamels 50)g		\$120.00	NX98-2163	MM3		
Order No.	Shade		ψ120.00	10/100 2100	IVIIVIO		
				\/MO OD M+-	F((+ 01 10-		ΦΩΩ ΩΩ
NX99-2181	Window, WIN				r Effect Opal, 12g		\$38.20
NX99-2191	Enamel ENL (for I			Order No.	Description		
NX99-2192	Enamel END (for a	dark shades)		NX98-2171	E01		
NX99-2201	Neutral NT			NX98-2172	E02		
				NX98-2173	E03		
\/M0 Enamala 25	50g		¢426 00	NX98-2174	E04		
	•		\$430.00				
Order No.	Shade			NX98-2175	E05		
NX100-2181	Window, WIN						
NX100-2191	Enamel ENL (for I	ight shades)		VM9 3D-Maste	r Effect Chroma, 12g		\$38.20
NX100-2192	Enamel END (for a	dark shades)		Order No.	Shade	Order No.	Shade
NX100-2201	Neutral NT	,		NX98-2121	EC1	NX98-2127	EC7
10/100 2201	riodiai iii			NX98-2122	EC2	NX98-2128	EC8
\/N40.0D.N4+	F		Φ04.00				
	Enamels, 12g		\$34.00	NX98-2123	EC3	NX98-2129	EC9
Order No.	Shade			NX98-2124	EC4	NX98-2130	EC10
NX98-2181	Window, WIN			NX98-2125	EC5	NX98-2131	EC11
NX98-2191	Enamel ENL (for I	ight shades)		NX98-2126	EC6		
NX98-2192	Enamel END (for o	-					
NX98-2201	Neutral NT	aam onaaoo,		VM9 3D-Master	r Chroma Plus, 12g		\$38.20
NX98-2221	Corrective, COR1			Order No.	Desc.	Order No.	Desc.
INAGO-ZZZ I	Corrective, Cont						
				NX98-2241	CP1	NX98-2244	CP4
VM9 3D-Master I	Effect Enamel, 12g .		\$38.20	NX98-2242	CP2	NX98-2245	CP5
Order No.	Shade	Order No.	Shade	NX98-2243	CP3		
NX98-2091	EE1	NX98-2097	EE7				
NX98-2092	EE2	NX98-2098	EE8	VM9 Gingiya 1	2g		\$38.20
NX98-2093	EE3	NX98-2099	EE9	Order No.	Shade	Order No.	Shade
NX98-2094	EE4	NX98-2100	EE10	NX98-2231	G1	NX98-2234	G4
NX98-2095	EE5	NX98-2101	EE11	NX98-2232	G2	NX98-2235	G5
NX98-2096	EE6			NX98-2233	G3		
VM9 3D-Master	Effect Liner, 12g		\$114.00	VM9 Margin, 12	2g		\$35.10
Order No.	Shade	Order No.	Shade	Order No.	Shade		
NX98-2141	EL1	NX98-2144	EL4	NX98-2261	M1 White		
NX98-2142	EL2	NX98-2145	EL5	NX98-2264	M4 Yellow		
NX98-2143	EL3	NX98-2146	EL6f	NX98-2265	M5 Amber		
				NX98-2267	M7 Light Beige		
VM9 Corrective,	12g		\$38.20	NX98-2268	M8 Pastel Brown		
Order No.	Shade			NX98-2269	M9 Light Orange		
NX98-2221	COR1				- 5		
NX98-2222	COR2						
NX98-2223	COR3						

VM9 Bonder, 12	<u>2g*</u>		\$34.00
Order No.	Desc.	Order No.	Desc.
NX98-2210	EB0	NX98-2213	EB3
NX98-2211	EB1	NX98-2214	EB4
NX98-2212	EB2	NX98-2215	EB5
VM9 Bonder Pa	ste, 5g*		\$80.50
VM9 Bonder Pa Order No.	ste, 5g* Desc.	Order No.	\$80.50 Desc.
	-		_
Order No.	Desc.	Order No.	Desc.
Order No. NX98-2310	Desc. EB0	Order No. NX98-2313	Desc. EB3

* While supplies last

VM®9 Classical Materials

VM9 Basic Kit Powder Classical

Order No. I	N98-0121	\$ FPO
16	Opaque, A1-D4	
16	Base Dentine, A1-D4	
3	Sun Dentine	
5	Chroma Plus	
2	Enamel ENL, END	
1	Neutral NT	
1	Window	
3	Corrective	
1	Modisol Pen	
1	Modeling Liquid, 50ml	
1	Firing Tray GI	

VM9 Basic Kit Small Classical

Shade indicator Classical shade guide

Fibrous firing pad supports

Order No.	N98-012	<u> </u>	 	 \$ F I	PC
		. .			

VM9 Base Der	ntine Classical 12g		\$34.00
Order No.	Shade	Order No.	Shade
N98-2431	A1	N98-2439	B4
N98-2432	A2	N98-2440	C1
N98-2433	A3	N98-2441	C2
N98-2434	A3.5	N98-2442	C3
N98-2435	A4	N98-2443	C4
N98-2436	B1	N98-2444	D2
N98-2437	B2	N98-2445	D3
N98-2438	B3	N98-2446	D4
VM9 Base Der	ntine Classical 50g		\$120.00
Order No.	Shade	Order No.	Shade
N99-2431	A1	N99-2439	B4
N99-2432	A2	N99-2440	C1
N99-2433	A3	N99-2441	C2
N99-2434	A35	N99-2442	C3
N99-2435	A4	N99-2443	C4
N99-2436	B1	N99-2444	D2
N99-2437	B2	N99-2445	D3
N99-2438	B3	N99-2446	D4

VM9 Transpa	Dentine Classical	12g	\$34.00
	Shade	Order No.	
N98-2451	A1	N98-2459	B4
N98-2452	A2	N98-2460	C1
N98-2453	A3	N98-2461	C2
N98-2454	A3.5	N98-2462	C3
N98-2455	A4	N98-2463	C4
N98-2456	B1	N98-2464	D2
N98-2457	B2	N98-2465	D3
N98-2458	B3	N98-2466	D4
VM9 Transpa	Dentine Classical	50g	\$120.00
VM9 Transpa Order No.		50g Order No.	\$120.00 Shade
	Shade	-	
Order No.	Shade	Order No.	Shade
Order No. N99-2451	Shade A1	Order No. N99-2459	Shade B4
Order No. N99-2451 N99-2452	Shade A1 A2	Order No. N99-2459 N99-2460	Shade B4 C1
Order No. N99-2451 N99-2452 N99-2453	Shade A1 A2 A3	Order No. N99-2459 N99-2460 N99-2461	Shade B4 C1 C2
Order No. N99-2451 N99-2452 N99-2453 N99-2454	Shade A1 A2 A3 A3.5	Order No. N99-2459 N99-2460 N99-2461 N99-2462	Shade B4 C1 C2 C3
Order No. N99-2451 N99-2452 N99-2453 N99-2454 N99-2455	Shade A1 A2 A3 A3.5 A4	Order No. N99-2459 N99-2460 N99-2461 N99-2462 N99-2463	Shade B4 C1 C2 C3 C4
Order No. N99-2451 N99-2452 N99-2453 N99-2454 N99-2455 N99-2456	Shade A1 A2 A3 A3.5 A4 B1	Order No. N99-2459 N99-2460 N99-2461 N99-2462 N99-2463 N99-2464	Shade B4 C1 C2 C3 C4 D2

VM Modeling Liquids

M Modeling Liqui Order No. V69-909 Effect Bonder Flui Order No. V68-910 Paste Fluid for VM	gd, 250ml gd, 50ml* g	\$78.50 \$27.90
Order No. GX78-414 GX78-415 GX78-416 GX78-417	rs	\$71.25
GX78-419	Professional Kit Small, 15 colors	



VITAVM_®9 EFFECT Porcelains

ALL-CERAMICS AND PRESSABLES

VITA VIV	VITA VM9 EFFECT LINER				
	EL1	Snow	White shade for 0-group (bleached)		
	EL2	Cream	Beige		
	EL3	Tabac	Brown		
	EL4	Golden Fleece	Yellow		
	EL5	Papaya	Orange		
	EL6	Sesame	Green-yellow		

VITA VM	VITA VM9 EFFECT CHROMA				
	EC1	Ghost	For white effects		
	EC2	Linen	Sand-beige		
	EC3	Pale banana	Light-yellow		
	EC4	Lemon Drop	Tender lemon yellow		
	EC5	Golden Rod	Light-orange		
	EC6	Sunflower	Orange		
	EC7	Light salmon	Pink		
	EC8	Toffee	Beige-brown		
	EC9	Doe	Brown		
	EC10	Larch	To reproduce green-brown effects		
	EC11	Gravel	To reproduce green-grey effects		

SUN DENTINE NEW			
	Light Yellow	for a sunnier and warmer shade	
	Light Orange		
	Orange-Red		

VITA VM	VITA VM9 MAMELON			
	MM1	Ecru	Beige	
	MM2	Mellow Buff	Warm yellow-brown	
	MM3	Peach Puff	Tender orange	

VITA VM	VITA VM9 EFFECT PEARL				
	EP1	Pearl	Shade in pastel-yellow		
	EP2	Pearl Blush	Shade in pastel-orange		
	EP3	Pearl Rose	Shade in pastel-rose		

EE1	Mint Cream	For whitish incisal areas, marginal
LLI	Willit Glean	ridges and triangular tubers
EE2	Pastel	For ridges and wedges at the labia surfaces
EE3	Misty Rose	For pink-translucent discoloration in the cervical and/or enamel area
EE4	Vanilla	For ridges and wedges
EE5	Sun Light	Yellowish-translucent to intensify depth in the approximal area
EE6	Navajo	For reddish-translucent discoloration, especially in older teeth
EE7	Golden Glow	Orange-translucent to intensify the effect of depth
EE8	Coral	Red-translucent to intensify the effect of depth
EE9	Water Drop	For subcoating in the enamel area for bluish incisal areas
EE10	Silver Lake Blue	For bluish wedges and ridges
EE11	Drizzle	For greyish wedges to increase the effect of depth

VITA VM9 EFFECT OPAL NEW				
	E01	Neutral	Neutral, universally suitable	
	E02	Whitish	Distinctive opal effect, for whitish wedges/triangular protuberances in occlusal area	
	E03	Bluish	Bluish	
	E04	Blue	Blue	
	E05	Dark Violet	Dark Violet	

VITA VM9 CORRECTIVE			
	E03	Opal Bluish	Bluish
	COR1	Incisal	Neutral
	COR2	Body	Beige
	COR3	Neck	Brown

VITA PM®9

VITA PM9 is VITA's latest restorative material in the VM line of small particle, 3D-Master shaded porcelain systems with low wear rates. PM9 is an extremely versatile press-to-zirconia system that enables technicians to consistently produce beautiful, high-strength restorations.

Pressable VITA PM9 is based on the proven fine-structure ceramic VITAVM9 and is used for pressing to yttria-stabilized ZrO2 frameworks in the CTE range of $10.4 - 10.6 \times 10^{-6}$, such as VITA In-Ceram YZ and other zirconium oxide substructure materials for natural looking highstrength anterior and posterior crowns and bridges.

VITA PM9 may also be pressed without the support of a zirconia substructure for the fabrication of single or multi-surface inlays, onlays, partial crowns, veneers and anterior crowns using the staining/glazing or layering technique. Restorations made from VITA PM9 can be easily characterized using the layering technique with the entire range of VITA VM9 ceramic materials or individualized (staining technique) with VITA Akzent or VITA Shading Paste stains. Each of the ten shaded pellets in the system is available in two levels of translucency.

PM®9 Accessory Kit

<i>O</i> rder No. VX36-9004		\$399.00
56	PM9 Investment, 100g	
1	Measuring Cup, 60ml	

PM9 Investment Liquid, 900ml

PM9 Investment System 200g, 3 parts

PM9 Disposable Plungers

Working Instructions, VITA PM9

Working Instructions, Investment

PM®9 Translucent Pellet Kit

order No.	VX36-9005\$649.00
10	PM9 Press Pellets "T" (0 M1P, 0 M2P, 1 M1P,

1 M2P, 2 M1P, 2 M2P, 2 M3P, 3 M1P, 3 M2P, 3 M3P), 5pcs each

VITA PM9 Shade Indicator, 10 Shade Tabs "T"

Working Instructions VITA PM9

PM®9 High-Translucent Pellet Kit

Order No. '	VX36-9006\$649.00	
	PM9 Press Pellets "HT" (0 M2P, 1 M1P, 1 M2P, 2 M2P, 3	
	M2P, ENOP, EN1P, EN2P, ENLP, ENDP) 5pcs each	

VITA PM9 Shade Indicator, 9 Shade Tabs "HT"

Working Instructions VITA PM9

2 VITA PM Working Instructions

PM®9 Sample Set

Order No.	VX36-9003	\$199.00
1	VITA PM9, 1M1P-T, 2pcs	
1	VITA PM9, ENLP-HT, 2pcs	
1	VITA PM9, EN1P-HT, 2pcs	
2	VITA VM9, ADD-ON, ADD 2+8, 12ml	
1	VITA VM Modeling Liquid, 50ml	
1	VITA Glaze LT Paste, 4g	
4	PM Investment Material, 100ml	
1	PM Investment Material Liquid, 100ml	

VM_®9 ADD-ON Kit

Order No. N98-0035	\$329.00
1 VM9 Add-On ADD1 (Transparent), 12g	
1 VM9 Add-On ADD2 (Light Enamel), 12g	
1 VM9 Add-On ADD3 (Dark Enamel), 12g	
1 VM9 Add-On ADD4 (White Transparent), 12g	
1 VM9 Add-On ADD5 (Yellow Transparent), 12g	
1 VM9 Add-On ADD6 (Orange Translucent), 12g	
1 VM9 Add-On ADD7 (Red Translucent), 12g	
1 VM9 Add-On ADD8 (Blue Translucent), 12g	
1 VITA VM9 Glaze LT, 12g	
 VITA VM Modeling Liquid, 50ml 	
1 VITA AKZENT Fluid, 20ml	
1 Fibrous Pad Firing Support	
1 A&E Brush, No. 3/0	

VITA VM9 Add-On Shade Indicator


Working Instructions VITA PM9

PM®9 Reorder Materials

VM9 Add-On Powe	ders		\$34.00
Order No.	Description		
NX98-2291	VM9 Add-On AE	D1 (Transparent)	
NX98-2292	VM9 Add-On AE	DD2 (Light Enamel)	
NX98-2293	VM9 Add-On AE	DD3 (Dark Enamel)	
NX98-2294	VM9 Add-On AE	DD4 (White Transparent)	
NX98-2295	VM9 Add-On AE	DD5 (Yellow Transparent)	
NX98-2296	VM9 Add-On AE	DD6 (Orange Transluscent)	
NX98-2297	VM9 Add-On AE	DD7 (Red Transluscent)	
NX98-2298	VM9 Add-On AD	DD8 (Blue Transluscent)	
PM9 Press Pellets	. 2g		. \$84.50
Order No.	Shade	Description	
VX36-9402	0M2P-HT	P-HT, High Translucent, 5pk	
VX36-9411	1M1P-HT	P-HT, High Translucent, 5pk	
VX36-9412	1M2P-HT	P-HT, High Translucent, 5pk	
VX36-9422	2M2P_HT	P-HT High Translucent 5nk	

VX36-9411	1M1P-HT	P-HT, High Translucent, 5pk
VX36-9412	1M2P-HT	P-HT, High Translucent, 5pk
VX36-9422	2M2P-HT	P-HT, High Translucent, 5pk
VX36-9432	3M2P-HT	P-HT, High Translucent, 5pk
VX36-9450	ENOP-HT	P-HT, High Translucent, 5pk
VX36-9451	EN1P-HT	P-HT, High Translucent, 5pk
VX36-9452	EN2P-HT	P-HT, High Translucent, 5pk
VX36-9453	ENLP-HT	P-HT, High Translucent, 5pk
VX36-9454	ENDP-HT	P-HT, High Translucent, 5pk
VX36-9201	0M1	P-T, Translucent, 5pk
VX36-9202	0M2	P-T, Translucent, 5pk
VX36-9211	1M1	P-T, Translucent, 5pk
VX36-9212	1M2	P-T, Translucent, 5pk
VX36-9221	2M1	P-T, Translucent, 5pk
VX36-9222	2M2	P-T, Translucent, 5pk
VX36-9223	2M3	P-T, Translucent, 5pk
VX36-9231	3M1	P-T, Translucent, 5pk
VX36-9232	3M2	P-T, Translucent, 5pk
VX36-9233	3M3	P-T, Translucent, 5pk
VX36-9301	0M1	P-O, Opaque, 5pk
VX36-9302	0M2	P-O, Opaque, 5pk
VX36-9311	1M1	P-O, Opaque, 5pk
VX36-9312	1M2	P-O, Opaque, 5pk
VX36-9321	2M1	P-O, Opaque, 5pk
VX36-9322	2M2	P-O, Opaque, 5pk
VX36-9323	2M3	P-O, Opaque, 5pk

ALL-CERAMICS AND PRESSABLES



PM9 Investment System, 3 parts	
Order No. J36-001	\$115.00
PM9 Investment, 50g x100g	
Order No. J36-002	\$159.00
PM9 Investment Liquid, 900ml	
Order No. J36-003	\$40.00
PM9 Disposable Plungers, 50 count	
Order No. J36-004	\$135.00
PM9 Firing Paste,1 Syringe, 12g	
Order No. A36-001	\$40.00
PM9 Firing Paste, 3 Syringes, 12g each	
Order No. A36-002	\$92.00
PM9 Color Indicator O (Opaque)	
Order No. GX36-901	\$71.2
PM9 Color Indicator T (Translucent)	
Order No. GX36-902	\$71.2
PM9 Color Indicator HT (High Translucent)	
Order No. GX36-903	\$71.2
VITA Glaze LT, Powder, 7.5g	
Order No. S28-201	\$49.00
VITA Glaze LT, Paste, 4g	
Order No. S28-202	\$39.00

Additional Accessories

VITA In-Ceram Interspace Varnish, 30ml	
Order No. A87-254	\$63.00
VITA In-Ceram Interspace Thinner, 30ml	
Ordor No. A87-255	¢ 50 በበ

VITA VM_®7

VITA VM7, the first of the VITA VM fine-grain family, is made to veneer alumina substructures. With an average expansion (CTE) of 7, and formulated to accurately reproduce VITA 3D-Master shades directly out of the bottle, it is ideal to use with VITA In-Ceram Alumina, VITA In-Ceram 2000 AL, Spinell and Zirconia.



After VITA VM7



Dr. Ed McLaren

VM_®7 Basic Kit, 12g

	. 0	
Order No. I	NX77-0001	\$1,230.00
3	VITA VM7 Effect Liner (EL2, EL3, EL4)	
3	VITA VM7 Chroma Plus (CP2, CP3, CP4)	
26	VITA VM7 Base Dentine	
2	VITA VM7 Enamel	600
1	VITA VM7 Neutral	
1	VITA VM7 Window	ALDER .
3	VITA VM7 Corrective (COR1, COR2, COR3)	40
1	VITA VM Veneering Liquid, 50ml	
1	Firing Tray G	
1	Fibrous pad firing support	
1	VITA VM7 shade indicator	
1	3D-Master shade guide	
1	Working Instructions	
1	Modisol Pen	

VM_®7 Basic Kit, 12g (Small)

Order No.	NX77-0002\$895.00
3	VITA VM7 Chroma Plus
15	VITA VM7 Base Dentine
2	VITA VM7 Enamel
1	VITA VM7 Neutral
1	VITA VM7 Window
1	VITA VM7 Corrective
1	VITA VM Veneering Liquid, 50ml
1	Firing Tray G
1	Fibrous pad firing support
1	VITA VM7 shade indicator
1	3D-Master shade guide
1	Working Instructions

VM_o7 Duildup Vit 12a

VIVI®/ E	Buildup Kit, 12g	
Order No. 26 1 1	NX77-0003 VITA VM7 Transpa Dentine VITA VM Modeling Liquid, 50ml Working Instructions	\$515.00
VM®7 E	Buildup Kit, 12g (Small)	
Order No. 15 1	NX77-0004 VITA VM7 Transpa Dentine VITA VM Modeling Liquid, 50ml Working Instructions	\$319.00

VM®7 Bleached Shade Kit, 12g

Order No.	NX77-0005	\$306.00
1	VITA VM7 Effect Liner	
3	VITA VM7 Base Dentine	
3	VITA VM7 Transpa Dentine	distr.
1	VITA VM7 Enamel (ENL), 12g	115
1	VITA VM7 Neutral, 12g	4
1	VITA VM7 Window, 12g	
1	VITA VM Modeling Liquid, 50ml	
1	Bleached shade guide group 0M	
1	Working Instructions	

VM®7 One Color Kit, 12g (3M2)

Order	No.	. NX77-0006	\$195.00
	1	VITA VM7 Base Dentine	
	1	VITA VM7 Transpa Dentine	
	1	VITA VM7 Enamel (ENL)	
	1	VITA VM7 Neutral	
	1	VITA VM7 Window	
	1	VITA VM7 Effect Liner	
	1	VITA VM Modeling Liquid, 50ml	
		Working Instructions	

VM_®7 One Color Kit, 12g (2M2)

Order No.	NX77-0010	\$195.00
1	VITA VM7 Base Dentine	
1	VITA VM7 Transpa Dentine	
1	VITA VM7 Enamel ENL Dentine	Die A
1	VITA VM7 Neutral	637
1	VITA VM7 Window	
1	VITA VM7 Effect Liner	
1	VITA VM Veneering Liquid, 50ml	
1	Working Instructions	

VM_®7 Professional Kit 12n

VITA VM7 Effect Opal

Working Instructions

VITA VM7 shade indicators

V IVI®/	riviessiviiai kit, izy	
Order No.	NX77-0007	\$1,330.00
11	VITA VM7 Effect Chroma	
11	VITA VM7 Effect Enamel	
6	VITA VM7 Effect Liner EL5, EL6	
3	VITA VM7 Mamelon	· A
3	VITA VM7 Effect Pearl	and the same

VM®7 Professional Kit, 12g (Small)

Order No.	NX77-0008\$570.00
5	VITA VM7 Effect Chroma
7	VITA VM7 Effect Enamel
1	VITA VM7 Mamelon
1	VITA VM7 Effect Pearl
1	VITA VM7 Effect Opal
3	VITA VM7 shade indicator
1	Working Instructions

VM®7 Gingiva Kit, 12g

\$242.00
600

VM_®7 Margin Kit 12g

Order No.	NX77-0011	\$214.00
6	VITA VM7 Margin	
1	VITA VM7 shade indicator	
1	Liquid	
	Working Instructions	1997

VM®7 Effect Opal

Order No. I	NX77-0012\$ FPO
5	Effect Opal
1	Shade Indicator
1	Working Instructions

VM®7 Reorder Materials

NX77-1038

NX77-1039

NX77-1040

NX77-1041

VM Veneering Material Liquid				
	Order No.V68-909, 50ml			
	Order No.V69-909	, 250ml		\$78.50
	VM Color Indicato	rs		\$71.25
	Order No.	Description		
	GX78-414	Effect Chroma, 11	colors	
	GX78-415		I/Mamelon, 15 colors	
	GX78-416	Effect Enamel, 11		
	GX78-417	Gingiva, 7 colors	001010	
	GX78-419	Professional Kit Sr	nall, 15 colors	
	VM7 Base Dentine	e Powders		
				¢3/I UC
	12g Order No.		Shade	ψ120.00
	NX77-1031	og order No.	0M1	
	NX77-1031 NX77-1032		0M2	
	NX77-1032 NX77-1033	_	0M3	
		— NIV70 1004		
	NX77-1034	NX78-1034	1M1	
	NX77-1035	NX78-1035	1M2	
	NX77-1036		2L1.5	
	NX77-1037	_	2L2.5	

NX78-1038

NX78-1039

NX78-1040



2M1

2M2

2M3

2R1.5

NX77-1042

NX//-1U42		2R2.5		VIVI / Effect Ena	amel, 12g
NX77-1043	NX78-1043	3L1.5		Order No.	Shade
NX77-1044	NX78-1044	3L2.5		NX77-1091	EE1
NX77-1045	NX78-1045	3M1		NX77-1092	EE2
NX77-1045		3M2			EE3
	NX78-1046			NX77-1093	
NX77-1047	NX78-1047	3M3		NX77-1094	EE4
NX77-1048	NX78-1048	3R1.5		NX77-1095	EE5
NX77-1049	NX78-1049	3R2.5		NX77-1096	EE6
NX77-1050	_	4L1.5			
NX77-1051	_	4L2.5		VM7 Effect Lin	er. 12a
NX77-1052	NX78-1052	4M1		Order No.	Shade
NX77-1053	NX78-1053	4M2		NX77-1141	EL1
NX77-1054	NX78-1054	4M3		NX77-1142	EL2
NX77-1055		4R1.5		NX77-1143	EL3
NX77-1056	_	4R2.5			
NX77-1057	_ _ _	5M1		VM7 Effect Pea	arl, 12g
NX77-1058		5M2		Order No.	Shade
NX77-1059	_	5M3		NIX77-1151	FP1
147.77 1000		OIVIO		NX77-1151	EP2
VM7 Transpa D	antina Davidara			NX77-1153	EP3
			Φ04.00	INV/1-1139	EF3
			\$120.00	VM7 Mamelon	, 12g
12g Order No.	50g Order No.	Shade		Order No.	Shade
NX77-1061		0M1		NX77-1161	MM1
NX77-1062		0M2		NX77-1162	MM2
NX77-1063	_	0M3		NX77-1163	MM3
NX77-1064	NX78-1064	1M1		14/1// 1100	IVIIVIO
				\/\ \/7 \ \\ Fff==+ \ \(\) =	al 10a
NX77-1065	NX78-1065	1M2		VM7 Effect Opa	
NX77-1066	_	2L1.5		Order No.	
NX77-1067	_	2L2.5		NX77-1171	E01
NX77-1068	NX78-1068	2M1		NX77-1172	E02
NX77-1069	NX78-1069	2M2		NX77-1173	E03
NX77-1070	NX78-1070	2M3			
NX77-1071	_	2R1.5		VM7 Window	
NX77-1072	_	2R2.5		Order No. NX7	7 1101 12a
NX77-1073	NX78-1073	3L1.5		Order No. NX7	3-1181, 50g
NX77-1074	NX78-1074	3L2.5			
NX77-1075	NX78-1075	3M1		VM7 Enamel El	
NX77-1076	NX78-1076	3M2		Order No. NX7	7-1191, 12g
NX77-1077	NX78-1077	3M3		Order No. NX7	3-1191, 50g
NX77-1078	NX78-1078	3R1.5			
NX77-1079	NX78-1079	3R2.5		VM7 Enamel El	ND
NX77-1080		4L1.5		Order No. NX7	
NX77-1081	— NIV/70 4000	4L2.5		Order No. NX7	5-1192, 50g
NX77-1082	NX78-1082	4M1			_
NX77-1083	NX78-1083	4M2		VM7 Neutral N	
NX77-1084	NX78-1084	4M3		Order No. NX7	7-1201, 12g
NX77-1085		4R1.5		Order No. NX78	3-1201, 50g
NX77-1086		4R2.5			
NX77-1087		5M1		VM7 Chroma P	lus 12a
NX77-1088		5M2		Order No.	Shade
NX77-1089	_	5M3		NX77-1242	CP2
				NX77-1243	CP3
VM7 Effect Chr	oma, 12g		\$40.20	NX77-1244	CP4
Order No.	Shade	Order No.	Shade		
NX77-1121	EC1	NX77-1127	EC7	VM7 Gingiva, 1	2g
NX77-1122	EC2	NX77-1128	EC8	Order No.	Shade
NX77-1123	EC3	NX77-1129	EC9	NX77-1231	G1
NX77-1124	EC4	NX77-1123	EC10	NX77-1231	G2
NX77-1125	EC5	NX77-1131	EC11	NX77-1233	G3
NX77-1126	EC6				

2R2.5

VM7 Effect Enamel, 12g.

.. \$40.20

.. \$98.00 Shade

.. \$40.20

.. \$40.20

... \$40.20

.. \$34.00 .. \$120.00

.. \$34.00 .. \$120.00

.. \$34.00 .. \$120.00

.. \$34.00 .. \$120.00

.. \$40.20

... \$40.20

Shade G4 G5

E04

E05

Shade

EE7

EE8

EE9 EE10

EE11

EL4

EL5 EL6

Order No.

NX77-1097

NX77-1098

NX77-1099 NX77-1100 NX77-1101

Order No. NX77-1144

NX77-1145

NX77-1146

NX77-1174

NX77-1175

Order No.

NX77-1234 NX77-1235

VITAVM _® 7	EFFECT	Porcelains
-----------------------	---------------	-------------------

VITA VM7 EFFECT LINER			
	EL1	Snow	White shade for 0-group (bleached)
	EL2	Cream	Beige
	EL3	Tabac	Brown
	EL4	Golden Fleece	Yellow
	EL5	Papaya	Orange
	EL6	Sesame	Green-yellow

VITA VM7 EFFECT CHROMA			
	EC1	Ghost	For white effects
	EC2	Linen	Sand-beige
	EC3	Pale banana	Light-yellow
	EC4	Lemon Drop	Tender lemon yellow
	EC5	Golden Rod	Light-orange
	EC6	Sunflower	Orange
	EC7	Light salmon	Pink
	EC8	Toffee	Beige-brown
	EC9	Doe	Brown
	EC10	Larch	To reproduce green-brown effects
	EC11	Gravel	To reproduce green-grey effects

ITA VI	17 EFFEC	T CHROMA	
	EC1	Ghost	For white effects
	EC2	Linen	Sand-beige
	EC3	Pale banana	Light-yellow
	EC4	Lemon Drop	Tender lemon yellow
	EC5	Golden Rod	Light-orange
	EC6	Sunflower	Orange
	EC7	Light salmon	Pink
	EC8	Toffee	Beige-brown
	EC9	Doe	Brown
	EC10	Larch	To reproduce green-brown effects
	EC11	Gravel	To reproduce green-grey effects

VITA VM	7 EFFECT	ENAMEL	
	EE1	Mint Cream	For whitish incisal areas, marginal ridges and triangular tubers
	EE2	Pastel	For ridges and wedges at the labial surfaces
	EE3	Misty Rose	For pink-translucent discoloration in the cervical and/or enamel area
	EE4	Vanilla	For ridges and wedges
	EE5	Sun Light	Yellowish-translucent to intensify depth in the approximal area
	EE6	Navajo	For reddish-translucent discoloration, especially in older teeth
	EE7	Golden Glow	Orange-translucent to intensify the effect of depth
	EE8	Coral	Red-translucent to intensify the effect of depth
	EE9	Water Drop	For subcoating in the enamel area for bluish incisal areas
	EE10	Silver Lake Blue	For bluish wedges and ridges
	EE11	Drizzle	For greyish wedges to increase the effect of depth

VITA VM7 MAMELON			
	MM1	Ecru	Beige
	MM2	Mellow Buff	Warm yellow-brown
	MM3	Peach Puff	Tender orange

VITA VM7 EFFECT PEARL			
	EP1	Pearl	Shade in pastel-yellow
	EP2	Pearl Blush	Shade in pastel-orange
	EP3	Pearl Rose	Shade in pastel-rose

VITA VM	7 EFFECT	OPAL NEW	
	E01	Neutral	Neutral, universally suitable
	E02	Whitish	Distinctive opal effect, for whitish wedges/triangular protuberances in occlusal area
	E03	Bluish	Bluish
	E04	Blue	Blue
	E05	Dark Violet	Dark Violet

VITA VM	7 CORRECTIVE		
	E03	Opal Bluish	Bluish

VITA VI	A7 GINGIVA		
	COR1	Incisal	Neutral
	COR2	Body	Beige
	COR3	Neck	Brown
	G1	Rose	Old rose
	G2	Nectarine	Orange-red
	G3	Pink Grapefruit	Pink
	G4	Rosewood	Brown-red
	G5	Cherry Brown	Black-red

VITA In-Ceram®

The proven In-Ceram family of all-ceramic substructure porcelains includes three materials possessing differing properties. The phenomenal success rate (99%+) of the In-Ceram family of materials is unsurpassed by any other all-ceramic system. They do not require enamel/dentine bonding and their varying levels of translucency enable ceramists to deliver restorations that satisfy every esthetic restorative need.

In-Ceram Alumina is indicated for use with single unit anterior or posterior crowns and three-unit anterior bridges. This popular substructure ceramic offers a high level of translucency and a flexural strength of 500 MPa.

In-Ceram Spinell is indicated for single unit anterior crowns, with a very high level of translucency and a flexural strength of 350 MPa.

In-Ceram Zirconia, with its very high strength, is indicated for single unit posterior and three-unit posterior bridges. In-Ceram Zirconia has a moderate level of translucency and a flexural strength of over 700 MPa.

In-Ceram Alumina, In-Ceram 2000 AL, In-Ceram Spinell and In-Ceram Zirconia are also available in machinable block form for CEREC inLab and other milling systems.

In-Ceram Sprint offers the lab the ability to control the production of Alumina crowns (no waiting for copings in the mail). Though limited to single Alumina crowns, the end result is the same outstanding strength and esthetics of In-Ceram Alumina. In-Ceram Sprint is ready to sinter in about one hour and easily comes off the die in just twenty minutes.





Dr. Ed McLaren

In-Ceram_® Alumina

28

Alumina Slip Powo Order No. V98-110 Alumina Special P		\$645.00
		\$225.00
Alumina Additive		
Order No. V98-117		\$39.20
Alumina/Zirconia I	Mixing Liquid, 20btls, 5ml	
Order No. V98-113		\$185.00
*Alumina Infiltrati	on Glass, 25g	\$340.00
Order No.	Shade	
V98-251	AL1 (A1, B1, C2, D2)	
V98-252	AL2 (A2, A3, B2, B3, D3)	
V98-253	AL3 (A3.5, A4, B4, C1, D4)	
V98-254	AL4 (C3, C4)	

In-Ceram® Spinell Assortment Kit

V98-004\$3	2,095.00
Spinell Glass Powder, 25g	
Spinell Slip Powder, 200g	
Special Plaster, 20g	
Mixing Liquid, 2.85ml	
Duplicating Rings, 28mm	
Vacuum Mixing Beakers (one rubber stopper)	
Brushes, Nos. 1 and 3	
Ice Cube Bags	
Each of the following: Acrylic Beaker, Measuring C	Cylinder,
Marking Pencil, Pencil Refill Pack (Red 12 ea), wor	king
instructions.	
	Spinell Glass Powder, 25g Spinell Slip Powder, 200g Special Plaster, 20g Mixing Liquid, 2.85ml Duplicating Rings, 28mm Vacuum Mixing Beakers (one rubber stopper) Brushes, Nos. 1 and 3 Ice Cube Bags Each of the following: Acrylic Beaker, Measuring 0 10ml, Glass Spatula, IC Brush, IC3 (Synthetic), Car Marking Pencil, Pencil Refill Pack (Red 12 ea), wor

In-Ceram_® Spinell

Order No. V98			\$645.00
	ell Slip Liquid, 16 vial -121 -2 85ml	is, 2.85mi	\$137 00
	r (all 3 materials)		φ107.00
Order No. V98-111			
*In-Ceram Spinell Infiltration Glass Powders, 25g\$340.00			
Order No.	Powder No.	Color	
V98-311	S11	White Gray	
V98-312	S12	Brownish Gray	
V98-313	S13	Yellow	
V98-314	S14	Brown	

In-Ceram® Zirconia Assortment Kit

Order No. \	/98-005	\$1,675.00
1	Zirconia Glass Powder, 25g	
1	Zirconia Slip Powder, 300g	
20	Alumina/Zirconia Mixing Liquid, 5ml	
1	Zirconia Additive, 5ml	
1	Alumina/Spinell/Zirconia shade indicator	

*In-Ceram® Zirconia

Infiltration Glass Powder, 25g	\$340.00			
Powder No.				
Z21				
Z22				
Z23				
Z24				
ler, 300g				
S	\$645.00			
5ml				
¹	\$38.20			
Alumina/Zirconia Slip Mixing Liquid, 20btls, 5ml				
3	\$185.00			
l 3 materials)				
	\$225.00			
	Powder No. Z21 Z22 Z23 Z24 ler, 300g 55ml Slip Mixing Liquid, 20btls, 5ml			

* While supplies last

In-Ceram® Accessories

*In-Ceram Duplicating Rings 2pk 25mm

Order No. A87-141	\$10.00
In-Ceram Spacer Varnish, 30ml	
Order No. A87-254	\$63.00
In-Ceram Thinner, 30ml	
Order No. A87-255	\$59.00
In-Ceram Marking Pen (lead not incl.)	
Order No. A41-452	\$21.00
In-Ceram Marking Pencil Leads, 12pcs	
Order No. A41-453	\$9.00
In-Ceram Vacuum mixing beaker w/ stopper	
Order No. A87-256	\$61.00
In-Ceram Mixing beaker (no stopper)	
Order No. A87-260	\$21.00
*In-Ceram Ice Cube Freezing Bags	
Order No. A87-262	\$19.00
In-Ceram Glass Mixing Spatula	
Order No. A51-263	\$15.00
IC3 Slipping Brush	
Order No. A51-213	\$14.00
IC4 Slipping Brush	
Order No. A51-211	\$20.00
In-Ceram Testing Liquid, 6ml	400.00
Order No. V98-116	\$29.90
In-Ceram Glass Color Indicator Alumina/Spinell/Zirconia	ቀ ሰብ ፓር
Order No. G98-909	\$83.50
In-Ceram Plaster Mixing Cup	Φ0.00
Order No. A87-257	\$9.00
In-Ceram Storage Cup w/lid Order No. A87-258	ቀ ፫ ሰር
	\$5.00
In-Ceram Pipette Order No. A87-259	¢1./ ∩(
Honeycomb Tray W, 2pk	\$14.00
Order No A61-216	¢/11.20
Platinum Pins, 6pk	φ41.Δι
Order No. A61-215	\$288 00
Platinum Foil, 60mm X 50mm	φ200.00
Order No. M98-125	\$/197 NC
Platinum Foil, 96mm diameter (round)	ψ437.00
Order No. M98-126	\$815 NC
In-Ceram Porcelain Tray, 2pk	ψυ ι J.UC
Order No. A61-212	\$169 00
Celay Alumina Optimizer, 10g	ψ100.00
Order No. V54-020	\$240 00
Celay Zirconia Optimizer, 8g	\$2 10.00
Order No. V54-025	\$297 nr
	4207.00

In-Ceramat III_® T Touchscreen Furnace

The Inceramat "high temp" sintering furnace is used for the In-Ceram technique as well as all prolonged high temperature sintering and glass infiltration requirments.

Features/Benefits

- High-heat quartz spiral muffle (rated 1200°c)
- 100 user-modifiable programs (10x In-Ceram preinstalled)
- Self-calibrates automatically after every 100 hours of use
- Firing table moves forward for easy access
- Fold-away twin cooling table
- Stainless steel construction
- Supplied complete with second muffle as back up
- 2 year warranty on parts and labor
- 1 year warranty on muffle

Order No. E98-051\$7,425.00

Vitasonic Mixer

The Vitasonic is a special frequency controlled ultrasonic with highly precise timing to ensure that correct flow and maximium thixotropic properties are achieved with the VITA In-Ceram slip mixing process. 2 year warranty.

Order No. E98-004\$690.00

* While supplies last



VITA Akzent Stains

Developed by Claude Sieber, Akzent can be used in both metal-ceramic and all-ceramic restorations. These stains enable the ceramist or dentist to achieve a more life-like three-dimensional effect.

Features/Benefits

- A selection of 20 colors based upon the spectrum of colors found in natural teeth.
- Akzent Stains demonstrate fluorescence as observed in natural
- Only completely safe metal oxides are used in Akzent Stains
- The stains are extremely fine-grained, resulting in excellent wear rate.
- The fine grain structure helps eliminate porosity on the surface of the restoration, which helps prevent plaque accumulation and

biocompatibility with the gingival area.

 When used with the Akzent finishing agent, a semi-translucent glazed effect is achieved, resulting in excellent light transmission.

VITA Akzent Stain Kit

Order No.	S28-004	\$1,170.00
20	Akzent Stains, 3g	
1	Akzent Finishing Agent 5g	

- Akzent Glaze, 5g
- Akzent Stain Liquid, 20ml
- Akzent Mixing Palette
- Glass Spatula—E30
- Sable Stain Brush #0
- Akzent Color Indicators
- Conversion Chart
- Working Instructions

VITA Akzent Stain Reorder Materials

Akzent Stains, 3g refill\$46.40		
Order No.	Description	
S28-760	Akz 01 (Birch)	
S28-761	Akz 02 (Mellow Yellow)	
S28-762	Akz 03 (Sun Kiss)	
S28-763	Akz 04 (Inca Gold)	
S28-764	Akz 05 (Gold Earth)	
S28-765	Akz 06 (Curry)	
S28-766	Akz 07 (Adobe)	
S28-767	Akz 08 (Almond)	
S28-768	Akz 09 (Burnt Clay)	
S28-769	Akz 10 (Muted Clay)	
S28-770	Akz 11 (Tabac)	
S28-771	Akz 12 (Redwood)	
S28-772	Akz 13 (Shak)	
S28-773	Akz 14 (Jet Black)	
S28-774	Akz 15 (Brown Stone)	
S28-775	Akz 16 (Seal Brown)	
S28-776	Akz 17 (Niagara)	
S28-777	Akz 18 (Fumo 01)	
S28-778	Akz 19 (Fumo 02)	
S28-779	Akz 20 (Fumo 03)	

Order No.	Description	
S28-780	Akzent Glaze, 5g	\$52.75
S28-784	Akzent Glaze, 30g	. \$239.00
S28-781	Akzent Finishing Agent, 5g	\$67.00
S28-128	Akzent Stain Liquid, 20ml	\$28.90
A51-240	Akzent Mixing Plate	. \$101.00
A51-157	Ceramic Mixing Palette	\$56.00
A51-261	Glass Spatula – E30	\$12.00
G68-420	Akzent Color Indicator #1 (1-10)	\$76.25
G68-421	Akzent Color Indicator #2 (11-20)	\$76.25

VITA Akzent Glaze Spray

Easy to apply, spray ceramic powder for simple, time-saving glazing of all ceramic restorations with a sintering temperature of 850°C.

Akzent Glaze Spray 125ml, 15g	
Order No. S28-790	\$122.00
Patterson Order No. 706-9628	
VITA LT Glaze Powder, 7.5g	
Order No. S28-201	\$FPO
Patterson Order No. 706-9776	
VITA LT Glaze Paste, 4g	
Order No. S28-201	\$122.00
Patterson Order No. 706-9784	

VITA Interno Stains

The Interno effect powders, developed by Claude Sieber, enable the ceramist to reproduce the

internal structures found in the interior of a natural tooth. Interno powders are for use in both VITA metal-ceramic and all-ceramic restorations.

Features/Benefits

- Increased luminescence due to intensified light transmission
- High fluorescence to increase the effects of light and depth
- Interno powders are particularly suitable for washing-in to the Dentine and enamel area
- They can be mixed with translucent powders for layering to increase the chromatic effect
- Used to create the illusion of depth when space is lacking

VITA Interno Effect Powder Kit

Working Instructions

Orders: 800-828-3839 USA | 800-263-4778 Canada | www.vident.com

Order No	. V41-0001	\$880.00
12	Interno Effect Powders, 3g	
1	Interno Liquid, 20ml	
1	Mixing Palette	
1	Glass Spatula – E30	
1	Stain Brush – 310	
1	Interno Shade Sample Disc	

VITA Interno Effect Powder - Reorder Materials

VITA Interno Stair	ns, 3g refill	\$59.00
Order No.	Description	
V41-1001	Int 01 (White Smoke)	
V41-1002	Int 02 (Sand)	
V41-1003	Int 03 (Lemon Chrome)	
V41-1004	Int 04 (Sunflower)	
V41-1005	Int 05 (Coral Blush)	
V41-1006	Int 06 (Gold Earth)	
V41-1007	Int 07 (Gold Misty Lilac)	
V41-1008	Int 08 (Marina)	
V41-1009	Int 09 (Cream Pearl)	
V41-1010	Int 10 (Burnt Orange)	
V41-1011	Int 11 (Clove)	
V41-1012	Int 12 (Splinter)	
Internal liquid 20	ml	
Interno Liquid, 20		Φ4E 40
	20	
Mixing Palette	0	¢101 00
		\$101.00
Glass Spatula, E3	u 1	¢12.00
		Φ12.00
Red Sable Kolinsk	0	ቀ20.00
		\$29.00
	•	ዕ ርብ 75
OTUEL INO. 008-92	<i>/</i>	გეგ./ე
Interno Shade Sai Order No. G68-92	mple Disc 7	\$58.75

VITA Shading Paste

These fluorescent fine-grained staining pastes are self-glazing and designed specifically for CEREC porcelain. Ready-to-use paste system does not require mixing a powder and a liquid.

VITA Shading Paste Classic Kit

Order No. \	V55-500\$835.00
8	Shading Pastes, 4g
1	VITA Glaze Paste, 7g
1	VITA Stain Liquid, 15ml
1	Shade Sample Disc
1	Mixing Palette
1	Agate Mixing Spatula
1	Stain Brush with Tips
1	Flat Squirrel Hairbrush

VITA Shading Paste - Classic Reorder Materials

Shading Paste Sta	ins, 4g refills		. \$66.00
Order No.	Description		
V55-501	White		
V55-502	Orange		
V55-503	Ochre		
V55-504	Khaki		
V55-505	Olive		
V55-506	Medium Brow	vn	
V55-507	Reddish Brow	'n	
V55-508	Mauve		
VITA Glaze Paste,	7g		
Order No. V55-515)		\$111.00
VITA Stain Liquid,	15ml		
Order No. V55-520)		. \$45.40

VITA Shading Paste 3D-Master Kit

Order No.	VX55-006	. \$925.00
9	Shading Pastes, 4g	
1	VITA Glaze Paste, 7g	
1	VITA Akzent Finishing Agent, 5g	
1	VITA Shading Paste Liquid, 15ml	
1	Shade Sample Blade	
1	Mixing Palette	
1	Glass Spatula	
1	Stain Brush with Tips	
1	Flat Squirrel Hairbrush	
1	Working Instructions	

VITA Shading Paste, 3D-Master Reorder Materials

Shading Paste Sta	ins, 4g refills	. \$66.00
Order No.	Description	
VX55-521	Birch	
VX55-522	Mellow Yellow	
VX55-523	Curry	
VX55-524	Fumo 1	
VX55-525	Fumo 2	
VX55-526	Fumo 3	
VX55-527	Tabac	
VX55-528	Shak	
VX55-529	Niagara	


VITA Extreme Esthetic Characterization Kit

Order No. NX98-0017 ..

8	VITA Akzent (AKZ01, AKZ02, AKZ03, AKZ04, AKZ09, AKZ12, AKZ16, AKZ17), 3g each
1	Akzent Finishing Agent (AKZ26)
1	Akzent Glaze (AKZ25), 5g
1	VITA Akzent Fluid, 20ml
1	VM9 Window (WIN), 12g
1	VM9 Enamlel Light (ENL), 12g
1	VM9 Corrective (COR1), 12g
1	VM9 Neutral (NT), 12g
1	VM9 Effect Enamel (EE9), 12g

VM9 Effect Enamel (EE1), 12g

VITA Glass Mixing Rod





. \$1.045.00

VITA VM®LC

INDIRECT COMPOSITE

VITA VMLC is a light curing microparticle paste resin for the indirect manufacture of fixed and removable restorations. The fine-particle inorganic filler, the size of a few nanometers, and its unique distribution, allow high translucency and natural scattering of light. VMLC exhibits excellent moldability, grindability and polishability which provides high esthetics. VMLC can be built up over both metal and resin frameworks and is indicated for veneering of full and partial crowns, bridges, inlays and onlays. It can also be used to individualize VITA Physiodens and Vitapan denture teeth.

VITA VM_®LC Basic Kit - 3D-Master

Order No. KX44-0003. \$2,695.00

- VITA VMLC Opaque, 10g 26
- VITA VMLC Chroma Plus CP1-CP5, 4g
- VITA VMLC Base Dentine 1M1-5M3, 4g
- 2 VITA VMLC Enamel ENL, END, 4g
- VITA VMLC Neutral NT, 4g
- VITA VMLC Window WIN, 2g
- VITA VMLC Opaque Liquid, 5ml
- VITA VMLC VM bond/LC Opaque Solvent, 10ml
- VITA VMLC Modeling Liquid, 30ml
- VITA VMLC Separator, 30ml
- flat brush, mixing spatula, measuring scoop
- VITA VM shade sample blades "individual"
- VITA Toothquide 3D-Master
- VITA VMLC working Instructions

VITA VM_®LC Basic Kit Small - 3D-Master

Order No. KX44-0004 . \$1,750.00

THE THE THE

- 15 VITA VMLC Opaque, 10q
- VITA VMLC Chroma Plus CP1-CP5, 4g
- 15 VITA VMLC Base Dentine, 4g ea. (1M1, 2M1, 2M3, 3L1.5, 3L2.5, 3M1, M2, 3M3, 3R1.5, 3R2.5, 4M1, 4M3, 5M2)
- VITA VMLC Enamel ENL, END, 4g
- VITA VMLC Neutral NT, 4g
- VITA VMLC Window WIN, 2g
- VITA VMLC Opaque Liguid, 10ml
- VITA VMLC VM bond/LC Opaque Solvent, 10ml
- VITA VMLC Modeling Liquid, 30ml
- VITA VMLC Separator, 30ml
- flat brush, mixing spatula, measuring scoop
- VITA VM shade sample blades"individual'
- VITA Toothguide 3D-Master
- VITA VMLC working Instructions
- VITA Toothquide 3D-Master

VITA VM_®LC Two Color Kit 3D-Master

Order No. KX44-0011

.\$FP0

- OPAQUE 2M2, 3M2, 10g
- BASE DENTINE 2M2, 3M2, 4g
- ENAMEL ENL, 4g
- CHROMA PLUS CP2, CP3, 4q
- Opaque Liquid, 5ml
- Modelling Liquid, 10ml
- Brush
- Black Mixing Palette Porcelain
- Measuring Spoon
- Mixing Spatula
- Shade Indicator Stick 2M2
- Shade Indicator Stick 3M2
- Working Instructions

VITA VM_®LC Bleached Color Kit

Order No. KX44-0008 \$450.00

- Opaque 0M1, 0M3, 10g
- Base Dentine 0M1, 0M3, 4g
- Enamel ENL, 4g
- Neutral NT, 4g
- Window WIN, 2g
- VITA VMLC Window WIN, 2g
- Opaque Liquid, 5ml
- Solvent for VITAVM Bond/LC Opaque, 10ml
- VITA VMLC Modeling Liquid, 30ml
- VITA VMLC Separator, 30ml
- flat brush, mixing spatula, measuring scoop
- Bleached Shade Guide Group 0M
- VITA VMLC Working Instructions

VITA VM_®LC Bleached Color Add-On Kit

Order No. KX44-0009 . \$160.00

- Opaque 0M1, 0M3, 10g
- Base Dentine 0M1, 0M3, 4g
- Enamel ENL, 4g
- Bleached Shade Guide Group 0M VITA VMLC Working Instructions

VITA VM_®LC Professional Kit

Order No. KX44-0005 \$397.00

- Color Opaque CO1,CO2, CO3, 10g
- VITA VMLC Effect Liner EL2-EL6,4q
- Effect Enamel EE3, EE6, EE9, 4g
- VITA VM shade sample blades "individual"
- VITA VMLC working Instructions

VITA VM_®LC Paint Kit

Order No. KX44-0007 \$357.00

- 10 Paint PT 1-19 2g
- WINDOW, 10g
- VITA VM shade sample blades "individual"
- VITA VMLC working Instructions

VITA VM®LC Gingiva Kit

Order No. KX44-0006 \$276.00

- Gingiva Opaque GOL, 10g
- Gingiva G1,G2,G4,G5
- VITA VM shade sample blades"individual"
- VITA VMLC working Instructions

VITA VM®LC Reorder Materials

VMLC Opaque	Refills 10g		\$42.30
Order No.	Shade	Order No.	Shade
KX44-4001	0M1	KX44-4016	3M2
KX44-4003	0M3	KX44-4017	3M3
KX44-4004	1M1	KX44-4018	3R1,5
KX44-4005	1M2	KX44-4019	3R2,5
KX44-4006	2L1,5	KX44-4020	4L1,5
KX44-4007	2L2,5	KX44-4021	4L2,5
KX44-4008	2M1	KX44-4022	4M1
<x44-4009< td=""><td>2M2</td><td>KX44-4023</td><td>4M2</td></x44-4009<>	2M2	KX44-4023	4M2
X44-4010	2M3	KX44-4024	4M3
X44-4011	2R1,5	KX44-4025	4R1,5
<x44-4012< td=""><td>2R2,5</td><td>KX44-4026</td><td>4R2,5</td></x44-4012<>	2R2,5	KX44-4026	4R2,5
(X44-4013	3L1,5	KX44-4027	5M1
(X44-4014	3L2,5	KX44-4028	5M2
<x44-4015< td=""><td>3M1</td><td>KX44-4029</td><td>5M3</td></x44-4015<>	3M1	KX44-4029	5M3
VMI C Raco Do	ntino Pofillo Aa		¢52.75
Viville base be Order No.	Shade	Order No.	\$55.75 Shade
X44-4031	OM1	KX44-4046	3M2
X44-4031 X44-4033	0M3	KX44-4047	3M3
(X44-4034	1M1	KX44-4047 KX44-4048	3R1,5
<x44-4034 <x44-4035< td=""><td>1M2</td><td>KX44-4049</td><td>3R2,5</td></x44-4035<></x44-4034 	1M2	KX44-4049	3R2,5
(X44-4036	2L1,5	KX44-4050	3112,5 4L1,5
X44-4030 X44-4037	2L1,5 2L2,5	KX44-4051	4L1,5 4L2,5
(X44-4037 (X44-4038	2L2,5 2M1	KX44-4052	4L2,5 4M1
(X44-4039	2M2	KX44-4053	4M2
(X44-4035 (X44-4040	2M3	KX44-4054	4M3
<x44-4040 <x44-4041< td=""><td>2R1,5</td><td>KX44-4055</td><td>4R1,5</td></x44-4041<></x44-4040 	2R1,5	KX44-4055	4R1,5
X44-4041 X44-4042		KX44-4056	4R2,5
X44-4042 X44-4043	2R2,5	KX44-4057	4nz,5 5M1
	3L1,5		
(X44-4044	3L2,5	KX44-4058	5M2
(X44-4045	3M1	KX44-4059	5M3
VMLC Effect Li	ner Refills 4g		\$53.75
Order No.	Shade	Order No.	Shade
KX44-4142	EL2	KX44-4145	EL5
KX44-4143	EL3	KX44-4146	EL6
KX44-4144	EL4		

Technical Help Desk: 888-249-1640 USA | 800-324-6224 Canada

Order No. KX44-4191 KX44-4192	Refills 4g Shade ENL END		\$53.75			
KX44-4201	NT					
VMLC Paint Ref Order No. KX44-4401 KX44-4403 KX44-4405	ills 2g Shade PT1 PT3 PT5	Order No. KX44-4413 KX44-4415 KX44-4417	\$33.00 Shade PT13 PT15 PT17			
KX44-4408 KX44-4409 KX44-4412	PT8 PT9 PT12	KX44-4419 KX44-4181	PT19 WIN			
Order No. KX44-4271 KX44-4281	Shade Gingiva GOL CO1	Order No. KX44-4282 KX44-4283	Shade CO2 CO3			
Order No. KX44-4231 KX44-4232	Shade G1 G2	Order No. KX44-4234 KX44-4235	Shade G4 G5			
VMLC Effect En Order No. KX44-4093 KX44-4096 KX44-4099	amels Refills 4g Shade EE3 EE6 EE9		\$53.75			
VMLC Modeling	Effect Liquid Refills VMLC Modeling Liquid 30ml					
			\$79.50			
	902	 10ml	\$12.40			
VMLC Solvent for Bond/LC Opaque 10ml Order No. V69-903\$5.00						
VMLC Separato						

Order No. V69-904



\$28.00

VITA VM_®LC Basic Kit Classical

- OPAQUE A1, A2, A3, A3.5, B3, D3, 10g
- BASE DENTINE A1, A2, A3, A3.5, B3, D3, 4g
- ENAMEL ENL, END, 4g
- NEUTRAL NT, 4g
- WINDOW WIN, 2g
- CHROMA PLUS CP1, CP4, 4g
- EFFECT LINER EL2, EL3, 4g 2
- Opaque Liquid, 5ml
- Solvent for VITA VM Bond/LC Opaque, 10ml
- Modelling Liquid, 30ml
- Separator for VITA VM LC, 30ml
- Brush

Order No. K44-0011

INDIRECT COMPOSITE

Order No. K44-0004

- Black Mixing Palette Porcelain
- Measuring Spoon
- Mixing Spatula
- VITA VM Shade Indicator "Individual"
- VITA Classical A1-D4 Shade Indicator
- Working Instructions

VITA VM_®LC Two Color Kit Classical

OPAQUE A2, A3, 10 g BASE DENTINE A2, A3, 4 q

- ENAMEL ENL, 4g
- CHROMA PLUS CP4, 4g
- EFFECT LINER EL3, 4g
- Opaque Liquid, 5ml
- Modelling Liquid, 10 ml
- Brush

K44-4445

34

- Black Mixing Palette Porcelain
- Measuring Spoon
- Mixing Spatula
- Shade indicator stick A2
- Shade indicator stick A3
- Working Instructions VMLC VITA Classical A-D shades, 10g

D3

K44-4401	A1		
K44-4402	A2		
K44-4403	A3		
K44-4404	A3.5		
K44-4408	B3		
K44-4415	D3		
VITA VM LC BAS	SE DENTINE, 4g		.\$ FPO
K44-4431	A1	ENL	
K44-4432	A2	ENL	
1144-4402	AZ	LINL	
K44-4433	A3	ENL	
K44-4433	A3	ENL	

END


VM Bond Kit

.\$**FPO**

.\$**FPO**

.\$**FPO**

Order No. KX44-0010 **SFPO**



- Bond LIQUID for VITA VM LC 5 ml
- Solvent for VM Bond/LC Opaque 10 i
- Mixing pad, 100 sheets pk
- Brush-holder, white
- Measuring spoon, mixing spatula
- Disposable brush tips, 50 pcs pk
- Instructions for use

VITA VM® Bond Reorder Materials

V69-910	VM Bond Powder 12g
V69-911	VM Bond Mixing liquid 5 ml
V69-903	VM Solvent for VM Bond/LC Opag

Cleaning Solution 10 ml

Measuring spoon/mixing spatula (2-sided)

VITA VM®CC Cold Cure Resin

VITA VMCC cold cure polymer resin is ideal for producing beautiful long term temporary restorations in VITA 3D Master shades. VITA VMCC is an unfilled, acrylate based polymer which polishes to a high luster. VMCC offers four minutes of working time and self cures completely in ten minutes. VMCC is available in 18 3D-Master shades. 15 Classical shades and one bleach shade (0M1). All shades are available in 30 and 100g powder.

VMCC can be used with any traditional technique for creating provisional restorations and may be placed with any temporary cement. When shade customized temporaries are required, VMCC is ideal because in addition to the dentine shades, there are also enamels necks, and translucents available.

Indications:

- For all applications involving temporary restorations and extra-oral repairs
- Full crowns and bridges
- Affixing shell crowns
- Affixing teeth to model casts
- Covering attachments/joints on fixed removable restorations

VITA VM®CC Basic Kit (16) Order No. KX45-0001

. \$770.00

- Dentine powders, 30g ea. (0M1,1M1,2L1.5, 2M1, 2M2, 2M3, 3L1.5, 3L2.5,3M1, 3M2, 3M3, 3R1.5, 3R2.5, 4L1.5 4M2, 5M1
- Enamel powders ENL, END, 30g
- Transparent powders WIN, 30g
- Liquid, 100ml
- Measuring Spoon
- VITA Toothquide 3D-Master
- VITA VMCC working Instructions

VITA VM®CC Basic Kit (5)

Order No. KX45-0002 \$317.00

- Dentine powders, 30g ea. (1M1, 2M2, 3M2 3R2.5, 4M2)
- VITA VMCC Enamel ENL, END, 4g
- Transparent powders WIN, 30g
- Liquid, 30ml
- Measuring Spoon
- VITA Toothguide 3D-Master
- VITA VMCC working Instructions

VITA VM®CC Professional Kit

\$167.00 Order No. KX45-0003.

5 Effect powders EL2-6, 30g ea.

VITA VM®CC Basic Kit Classical

Order No. K45-0001. \$735.00

- Base Dentine powders A1-D4 (No B1), 30g ea.
- Enamel ENL, END, 30g
- Window WIN, 30a
- VMCC Liquid, 100ml
- Measuring spatula Shade Indicator
- Working instructions

VITA VM®CC Basic Kit Small Classical

Order No. K45-0002. \$346.00

- 6 Base Dentine powders A1,A2,A3,A3.5,B3,D3
- Working instructions

VITA VM®CC Assortment Kit Classical (extension kit for VMCC 3D users)

Order No. K45-0003 \$630.00

- Base Dentine powders A1-D4 (No B1), 30g ea.
- VMCC Liquid, 100ml

Shade Indicator

Measuring Spoon

VITA VM®CC Reorder Materials

VITA VMCC 3D-Master Base Dentine

100g			\$133.00
30g Order No.	100g Order No.	Shade	
KX45-7031	KX46-7031	0M1	
KX45-7034	KX46-7034	1M1	
KX45-7036	KX46-7036	2L1.5	
KX45-7038	KX46-7038	2M1	
KX45-7039	KX46-7039	2M2	
KX45-7040	KX46-7040	2M3	
KX45-7043	KX46-7043	3L1.5	
KX45-7044	KX46-7044	3L2.5	
KX45-7045	KX46-7045	3M1	
KX45-7046	KX46-7046	3M2	
KX45-7047	KX46-7047	3M3	
KX45-7048	KX46-7048	3R1.5	
KX45-7049	KX46-7049	3R2.5	
KX45-7050	KX46-7050	4L1.5	
KX45-7053	KX46-7053	4M2	
KX45-7057	KX46-7057	5M1	
KX45-7191	KX46-7191	ENL	
KX45-7192	KX46-7192	END	
	sical Base Dentine		
30g			\$42.30
100g			\$133.00

K45-7433	K46-7433
K45-7434	K46-7434
K45-7435	K46-7435
K45-7437	K46-7437

30g Order No.

K45-7431

K45-7432

K45-7445

K45-7446

B2 В3 K45-7438 K46-7438 K45-7439 B4 K45-7440 C1 K46-7441 C2 K45-7441 K45-7442 K46-7442 C3K45-7443 C4K45-7444 D2

K46-7445

100g Order No.

K46-7431

K46-7432

Shade

Α1

A2

Α3

A3.5

A4

D3

D4

VITA VMCC Transparent (WIN), 30g Order No. KX45-7181

Order No. KX46-710.

VITA VMCC	Effect Liner, 30g.	\$42.30
Order No.	Shade	
KX45-7142	EL2	
KX45-7143	EL3	
KX45-7144	EL4	
KX45-7145	EL5	

EL6 KX45-7146 VITA VMCC Liquid, 30ml Order No. KX45-730. . \$42.30 VITA VMCC Liquid, 100ml



\$126.00

\$42.30

36

VITABLOCS® RealLife

VITABLOCS RealLife are the first CAD/CAM machinable ceramic blocks to offer three-dimensional esthetics for the fabrication of highly esthetic anterior restorations*. Only suitable for MC XL milling units.



Vident Order No.	Patterson Order No.	Shade	Qty	Size
VX55-4476	707-1749	0M1	5	14/14mm
VX55-4479	707-1756	1M1	5	14/14mm
VX55-4480	707-1715	1M2	5	14/14mm
VX55-4483	707-1764	2M1	5	14/14mm
VX55-4484	707-1723	2M2	5	14/14mm
VX55-4491	707-1731	3M2	5	14/14mm

VITABLOCS® for Rapid Layer Technology

VITABLOCS Mark II and TriLuxe forte in a 40/19 size for highly efficient CAD/CAM veneering of bridge substructures with CEREC/inLab MC XL. *Only suitable for MC XL milling units*.

Call Patterson Dental for Block Orders: 800-873-7683 US • 800-570-6356 Canada

Vident Order No.	Patterson Order No.	Shade	Qty	Size
TriLuxe forte				
VX55-4725	707-3398	1M2	2	40/19mm
VX55-4729	707-3406	2M2	2	40/19mm
VX55-4736	707-3414	3M2	2	40/19mm
Mark II				
VX55-4501	707-2366	1M1	2	40/19mm
VX55-4502	707-2374	1M2	2	40/19mm
VX55-4506	707-2382	2M2	2	40/19mm
VX55-4513	707-2390	3M2	2	40/19mm

VITABLOCS® Mark II, 3D-Master® Shades

Fine particle feldspar ceramic indicated for inlays, onlays, full and partial crowns in the anterior and posterior area. Highly translucent monochromatic blocs guarantee excellent shade matching with the patient's remaining tooth structure. Clinically proven with over 20 years of documented success.*




				_	
Vident Order No.	Patterson Order No.	Shade	Block Type	Qty	Size
VX55-4300	707-3729	0M1	MarkII	10	112
VX55-4378	707-4214	0M1	MarkII	5	114
VX55-4301	707-3737	1M1	MarkII	10	18
VX55-4327	707-3828	1M1	MarkII	10	l10
VX55-4353	707-3919	1M1	MarkII	10	112
VX55-4379	707-4222	1M1	MarkII	5	114
VX55-4302	707-3745	1M2	MarkII	10	18
VX55-4328	707-3836	1M2	MarkII	10	I10
VX55-4354	707-3927	1M2	MarkII	10	112
VX55-4380	707-4230	1M2	MarkII	5	114
VX55-4305	707-3752	2M1	MarkII	10	18
VX55-4331	707-3844	2M1	MarkII	10	l10
VX55-4357	707-3935	2M1	MarkII	10	112
VX55-4383	707-4248	2M1	MarkII	5	114
VX55-4306	707-3760	2M2	MarkII	10	18
VX55-4332	707-3851	2M2	MarkII	10	I10
VX55-4358	707-3943	2M2	MarkII	10	112
VX55-4384	707-4255	2M2	MarkII	5	114
VX55-4307	707-3778	2M3	MarkII	10	18
VX55-4333	707-3869	2M3	MarkII	10	110
VX55-4359	707-3950	2M3	MarkII	10	112
VX55-4385	707-4263	2M3	MarkII	5	114
VX55-4312	707-3786	3M1	MarkII	10	18
VX55-4338	707-3877	3M1	MarkII	10	110
VX55-4364	707-3968	3M1	MarkII	10	112
VX55-4390	707-4271	3M1	MarkII	5	114
VX55-4313	707-3794	3M2	MarkII	10	18
VX55-4339	707-3885	3M2	MarkII	10	110
VX55-4365	707-3976	3M2	MarkII	10	112
VX55-4391	707-4289	3M2	MarkII	5	114
VX55-4314	707-3802	3M3	MarkII	10	18
VX55-4340	707-3893	3M3	MarkII	10	l10
VX55-4366	707-3984	3M3	MarkII	10	112
VX55-4392	707-4297	3M3	MarkII	5	l14
VX55-4320	707-3810	4M2	MarkII	10	18
VX55-4346	707-3901	4M2	MarkII	10	110
VX55-4372	707-3992	4M2	MarkII	10	112
VX55-4398	707-4305	4M2	MarkII	5	114

VITABLOCS® Mark II, Classical Shades

Fine particle nethaline feldspar ceramic. Indicated for crowns, veneers, inlays, and onlays.*

Vident Order No.	Patterson Order No.	Shade	Block Type	<i>Qty</i>	Size
V55-4301	707-3141	A1	Mark	10	18
V55-4311	707-3174	A1	MarkII	10	l10
V55-4321	707-3208	A1	MarkII	10	l12
V55-4341	707-3315	A1	MarkII	5	114
V55-4302	707-3158	A2	MarkII	10	18
V55-4312	707-3182	A2	MarkII	10	l10
V55-4322	707-3216	A2	Markll	10	l12
V55-4342	707-3323	A2	MarkII	5	114
V55-4303	707-3166	A3	MarkII	10	18
V55-4313	707-3190	A3	MarkII	10	l10
V55-4323	707-3224	A3	MarkII	10	l12
V55-4343	707-3331	A3	MarkII	5	114
V					

^{*} Mandrel fits both compact and MC XL milling units



VITABLOCS® TriLuxe forte

Polychromatic TriLuxe forte feldspar ceramic blocks feature fine nuances of shade transition from the incisal edge to the cervical area. With more intense color saturation and increased fluorescence, TriLuxe forte blocs provide unsurpassed shade effects and a chameleon blend with natural dentition.*

Vident Order No.	Patterson Order No.	Shade	Block Type	Оty	Size	
VX55-4780	707-5781	1M2	TriLuxe Forte	10	12mm	
VX55-4784	707-5799	2M2	TriLuxe Forte	10	12mm	
VX55-4791	707-5807	3M2	TriLuxe Forte	10	12mm	
VX55-4754	707-5716	1M2	TriLuxe Forte	5	14/12mm	
VX55-4758	707-5724	2M2	TriLuxe Forte	5	14/12mm	
VX55-4765	707-5732	3M2	TriLuxe Forte	5	14/12mm	
VX55-4750	707-5633	1M2	TriLuxe Forte	5	14/14mm	
VX55-4751	707-5641	2M2	TriLuxe Forte	5	14/14mm	
VX55-4752	707-5658	3M2	TriLuxe Forte	5	14/14mm	
New Classical Shad	des					
V55-4751	707-5815	A1	TriLuxe Forte	10	12mm	
V55-4752	707-5823	A2	TriLuxe Forte	10	12mm	
V55-4753	707-5831	A3	TriLuxe Forte	10	12mm	NEW
V55-4851	707-5740	A1	TriLuxe Forte	5	14/12mm	NE
V55-4852	707-5757	A2	TriLuxe Forte	5	14/12mm	
V55-4853	707-5765	A3	TriLuxe Forte	5	14/12mm	•

5 14/12mm

5 14/12mm

VITABLOCS® TriLuxe

Polychromatic TriLuxe feldspar ceramic blocks feature three levels of tooth shades representing the enamel, dentine and cervical tooth areas within the same bloc.*

Vident Order No.	Patterson Order No.	Shade	Block Type	<i>Qty</i>	Size
VX55-4680	707-4628	1M2	TriLuxe	10	12mm
VX55-4684	707-4636	2M2	TriLuxe	10	12mm
VX55-4691	707-4644	3M2	TriLuxe	10	12mm
VX55-4654	707-4545	1M2	TriLuxe	5	14/12mm
VX55-4658	707-4552	2M2	TriLuxe	5	14/12mm
VX55-4665	707-4560	3M2	TriLuxe	5	14/12mm
VX55-4650	707-4461	1M2	TriLuxe	5	14/14mm
VX55-4651	707-4479	2M2	TriLuxe	5	14/14mm
VX55-4652	707-4487	3M2	TriLuxe	5	14/14mm
New Classical Share	des				
V55-4651	707-4415	A1	TriLuxe	10	12mm
V55-4652	707-4123	A2	TriLuxe	10	12mm
V55-4653	707-4131	A3	TriLuxe	10	12mm
V55-4681	707-4578	A1	TriLuxe	5	14/12mm



VITA CAD-Temp™

707-4586

707-4594

V55-4682

V55-4683

Composite blocks for the milling of high quality, fully and partially anatomical long-term provisional crowns and multi-unit bridges. Available in monocolor and polychromatic multi-color for lasting color, stability and esthetics.*

TriLuxe

TriLuxe

A2

А3

Vident Order No.	Patterson Order No.	Shade	Block Type	Qty	Size
VX59-4002	707-4677	0M1T	CAD-Temp	2	15.5 x 19 x 39mm
VX59-4010	707-4347	0M1T	CAD-Temp	10	15.5 x 19 x 39mm
VX59-4102	707-4685	1M2T	CAD-Temp	2	15.5 x 19 x 39mm
VX59-4110	707-4354	1M2T	CAD-Temp	10	15.5 x 19 x 39mm
VX59-4202	707-4693	2M2T	CAD-Temp	2	15.5 x 19 x 39mm
VX59-4210	707-4362	2M2T	CAD-Temp	10	15.5 x 19 x 39mm
VX59-4302	707-4701	3M2T	CAD-Temp	2	15.5 x 19 x 39mm
VX59-4310	707-4370	3M2T	CAD-Temp	10	15.5 x 19 x 39mm
VX59-4155	707-6748	1M2T	CAD-Temp CT-55 [†]	1	15 x 19 x 55mm
VX59-4255	707-6755	2M2T	CAD-Temp CT-55 [†]	1	15 x 19 x 55mm
VX59-4355	707-6763	3M2T	CAD-Temp CT-55 [†]	1	15 x 19 x 55mm
VX59-4121	707-6581	1M2T	CAD-Temp MultiColor	2	15.5 x 19 x 39mm
VX59-4122	707-6599	1M2T	CAD-Temp MultiColor	10	15.5 x 19 x 39mm
VX59-4221	707-6607	2M2T	CAD-Temp MultiColor	2	15.5 x 19 x 39mm
VX59-4222	707-6664	2M2T	CAD-Temp MultiColor	10	15.5 x 19 x 39mm
VX59-4321	707-6672	3M2T	CAD-Temp MultiColor	2	15.5 x 19 x 39mm
VX59-4322	707-6680	3M2T	CAD-Temp MultiColor	10	15.5 x 19 x 39mm

Call Patterson Dental for Block Orders: 800-873-7683 US • 800-570-6356 Canada



For the processing of VITA CAD-Waxx and CAD-Temp in the compact inLab and CEREC 3 systems (not in the MC XL), the Sirona CAD-Waxx Starter Set with a modified tank and reinforced filter system is required.

- * Mandrel fits both compact and MC XL milling units
- † Only suitable for MC XL milling unit





38



In-Ceram® YZ for inLab

MACHINABLE CERAMICS

Sintering ceramic indicated for crown copings and multi-unit bridge substances in the anterior and posterior area. This ytrium-stabilized zirconium dioxide is the "diamond" among dental ceramics for the construction of high strength, attractive crown and bridge substructures (inVizion restorations). Unsurpassed strength for highly stable, thin substructures. Flexural strength ≥900MPa. Veneer with VITA VM9.*

g, ctable, cabet.	actareor riexarar en en engar =ce	on a ronce mar min time		100
Vident Order No.	Patterson Order No.	Block Type	Qty	Size
V59-4141	707-1327	YZ 14/13	20	13 x 13 x 14mm
V59-4205	707-3000	YZ 20/15	5	14 x 15 x 20mm
V59-4220	707-0972	YZ 20/15	20	14 x 15 x 20mm
V59-4291	707-3026	YZ 20/19	4	14 x 19 x 20mm
V59-4293	707-0980	YZ 20/19	16	14 x 19 x 20mm
V59-4402	707-3042	YZ 40/15	2	14 x 15 x 40mm
V59-4410	707-3059	YZ 40/15	10	14 x 15 x 40mm
V59-4491	707-3067	YZ 40/19	2	14 x 19 x 39mm
V59-4492	707-3075	YZ 40/19	10	14 x 19 x 39mm
V59-551	706-9636	YZ 55/20 flip†	1	15.5 x 19 x 55mm
V59-554	707-3109	YZ 55/20 flip†	4	15.5 x 19 x 55mm
V59-4551	707-3083	YZ 55/20‡	1	15.5 x 19 x 55mm
V59-4554	707-3117	YZ 55/20‡	4	15.5 x 19 x 55mm
V59-4651	707-2077	YZ 65/25‡	1	22 x 25 x 65mm
V59-4653	707-2085	YZ 65/25‡	3	22 x 25 x 65mm
V59-4662	707-1350	YZ 65/40‡	1	17 x 40 x 65mm
V59-4671	707-1343	YZ 85/40‡	1	22 x 40 x 85mm

[†] Only suitable for inLab Compact milling unit. NOT suitable for MC XL

In-Ceram® YZ Pre-Shaded for inLab

Replaces the need for VITA YZ Coloring Liquid to obtain an optimally shaded substructure base. These industrially pre-shaded zirconium oxide blocks reduce valuable working time and permit homogeneous shading of YZ substructures.*

Vident Order No.	Patterson Order No.	Block Type	Qty	Size	Con Con
V59-4142	707-1285	YZ, Color LL1p 14/13	20	13 x 13 x 14mm	4
V59-4240	707-1145	YZ, Color LL1p 20/15	5	14 x 15 x 20mm	
V59-4242	707-0988	YZ, Color LL1p 20/15	20	14 x 15 x 20mm	
V59-4243	707-1160	YZ, Color LL1p 20/19	4	14 x 19 x 20mm	
V59-4245	707-1079	YZ, Color LL1p 20/19	16	14 x 19 x 20mm	
V59-4440	707-1186	YZ, Color LL1p 40/15	2	14 x 15 x 40mm	
V59-4441	707-1194	YZ, Color LL1p 40/15	10	14 x 15 x 40mm	
V59-4443	707-1202	YZ, Color LL1p 40/19	2	14 x 19 x 39mm	
V59-4444	707-1210	YZ, Color LL1p 40/19	10	14 x 19 x 39mm	
V59-540	707-1228	YZ, Color LL1p 55/20 flip†	1	15.5 x 19 x 55mm	
V59-541	707-1236	YZ, Color LL1p 55/20 flip†	4	15.5 x 19 x 55mm	
V59-4540	707-1244	YZ, Color LL1p 55/20‡	1	15.5 x 19 x 55mm	
V59-4541	707-1251	YZ, Color LL1p 55/20‡	4	15.5 x 19 x 55mm	
V59-4640	707-1269	YZ, Color LL1p 65/25‡	1	22 x 25 x 65mm	
V59-4641	707-1277	YZ, Color LL1p 65/25‡	3	22 x 25 x 65mm	
V59-4663	707-1368	YZ, Color LL1p 65/40‡	1	17 x 40 x 65mm	
V59-4675	707-1061	YZ, Color LL1p 85/40‡	1	22 x 40 x 85mm	

Call Patterson Dental for Block Orders: 800-873-7683 US • 800-570-6356 Canada



High molecular translucent acrylpolymer for the construction of PFM metal frameworks for crowns & bridges. Also indicated as a drilling template in implantology. Burns away cleanly during the lost-wax technique. Not for clinical use.*



Vident Order No.	Patterson Order No.	Block Type	Qty	Size
V52-4402	707-2994	CAD-Waxx	2	14 x 15 x 40mm
V52-4410	707-2986	CAD-Waxx	10	14 x 15 x 40mm
_	636-6678	CAD-Waxx starter set		

Note: For the processing of VITA CAD-Waxx and CAD-Temp in the compact inLab and CEREC 3 systems (not in the MC XL), the Sirona CAD-Waxx Starter Set with a modified tank and reinforced filter system is required.

In-Ceram® AL

Partially sintered alumina blocks suitable for tooth colored primary components, crowns, and anterior bridge substructures with up to 3 units. Veneer with VITA VM7.*



Vident Onden Ne	D-++ O N-	Dis als Torra	0+	0:	
Vident Order No.	Patterson Order No.	Block Type	Qty	Size	
V58-4204	707-2879	AL - 20/19	4	14 x 19 x 20mm	
V58-4224	707-2861	AL - 20/19	24	14 x 19 x 20mm	
V58-4402	707-2887	AL - 40/19	2	14 x 19 x 39mm	
V58-4410	707-2895	AL - 40/19	10	14 x 19 x 39mm	

In-Ceram® Classic Spinell for inLab

Features the highest translucency of all oxide ceramics for excellent esthetics. Indicated for anterior, single unit crown copings. 350 MPa.*



Vident Order No.	Patterson Order No.	Block Type	Qty	Size
V55-4355	707-2978	Spinell CS-11	10	10 x 11 x 15mm

In-Ceram® Classic Alumina for inLab


Provides the perfect combination of esthetics and strength. Indicated for crowns and anterior bridge substructures with up to 3 pontics. 450-600 MPa.*



Vident Order No.	Patterson Order No.	Block Type	Qty	Size
V55-4350	707-2960	Alumina CA-12	10	10 x 12 x 15mm
V54-4104	707-2903	Alumina BA-28	2	14 x 15 x 28mm
V55-4349	707-2952	Alumina CA-40	2	14 x 15 x 40mm

In-Ceram® Classic Zirconia for inLab

Features the highest masking power for discolored preparations. Indicated for crown and bridge substructures with up to 3 units in the anterior and posterior area. 650-700 MPa.*



^{*} Mandrel fits both compact and MC XL milling units







[‡] Only suitable for inLab MC XL milling unit

[†] Only suitable for inLab Compact milling unit. NOT suitable for MC XL

[‡] Only suitable for inLab MC XL milling unit

^{*} Mandrel fits both compact and MC XL milling units

VITABLOCS® for KaVo Everest

MACHINABLE CERAMICS

VITABLOCS Mark II for KaVo Everest are made from VITA's fine-structure feldspar ceramic offering outstanding clinical reliability. Indicated for inlays, onlays, veneers and pardial or full coverage crowns. These blocks consist of a dentin and enamel blend of ceramics for highly esthetic results.

Vident Order No.	Description	Qty	Price
VX55-2354	VITABLOCS for KaVo 1M2C – I12	5 pieces	\$118.00
VX55-2358	VITABLOCS for KaVo 2M2C – I12	5 pieces	\$118.00
VX55-2365	VITABLOCS for KaVo 3M2C – I12	5 pieces	\$118.00
VX55-2380	VITABLOCS for KaVo 1M2C – I14	5 pieces	\$169.00
VX55-2384	VITABLOCS for KaVo 2M2C – I14	5 pieces	\$169.00
VX55-2391	VITABLOCS for KaVo 3M2C – I14	5 pieces	\$169.00

VITA CAD-Temp™ DISCS for KaVo Everest

CAD-Temp for KaVo Everest systems is suitable for the fabrication of multi-unit, fully or partially anatomical long-term provisional restorations with up to 2 pontics. Up to 35 restorations can be milled from each disc. An integrated "chip" ensures that the system detects which restorations have already been milled and how many can still be produced. Disc Size 100 x 20mm.

Vident Order No.	Description	Qty	Price
VX58-2101	CAD-Temp for KaVo 1M2T-100 x 20mm	1 piece	\$199.00
VX58-2201	CAD-Temp for KaVo 2M2T-100 x 20mm	1 piece	\$199.00
VX58-2301	CAD-Temp for KaVo 3M2T-100 x 20mm	1 piece	\$199.00

VITA In-Ceram® YZ DISCS for Open CAD/CAM Systems

VITA In-Ceram YZ made of zirconium dioxide partially stabilized with yttrium in DISC format with a circumferential groove. This extra-large geometry is designed for the manufacture of multi-unit bridge substructures and the economical manufacture of single crowns using batch milling in open CAD/CAM Systems. The enlargement factor in all 3 spatial dimensions is clearly imprinted on each DISC.

Vident Order No.	Description	Oty	Price
V59-1401	VITA In-Ceram YZ DISC-98mm x 14mm	1 piece	\$375.00*
V59-1801	VITA In-Ceram YZ DISC-98mm x 18mm	1 piece	\$450.00*
V59-2001	VITA In-Ceram YZ DISC-98mm x 20mm	1 piece	\$500.00*

VITA CAD-Temp™ DISCS for Open CAD/CAM Systems

The provisional material VITA CAD-Temp monoColor available in DISC format with an external diameter of 98mm (including circumferential groove) for open CAD/CAM Systems. VITA CAD-Temp is a CAD/CAM composite material for full and partially anatomical long-term provisional restorations, with the perfect combination of strength (approx. 80 MPa flexural strength) and elasticity (2800 MPa). Available in three shades for exceptional esthetics.

Vident Order No.	Description	Qty	Price	
VX57-2101	CAD-Temp DISC 1M2T-98mm x 20mm	1 piece	\$225.00*	
VX57-2201	CAD-Temp DISC 2M2T-98mm x 20mm	1 piece	\$225.00*	
VX57-2301	CAD-Temp DISC 3M2T-98mm x 20mm	1 piece	\$225.00*	
VX57-2401	CAD-Temp Multi-Colored DISC 1M2T-98mm x 20mm	1 piece.		
VX57-2501	CAD-Temp Multi-Colored DISC 2M2T-98mm x 20mm	1 piece.		
VX57-2601	CAD-Temp Multi-Colored DISC 3M2T-98mm x 20mm	1 piece.		MIEVY
*Auantity Discount P	ricing available Appual Purchase Contracts available Co	ntact Vident for additio	nal information	

Uuantity Discount Pricing available. Annual Purchase Contracts available. Contact Vident for additional information

VITA CAD-Waxx™ DISCS for Open CAD/CAM Systems

CAD-Waxx DISCS are made of a millable, filler-free acrylate polymer which can be used in place of modeling wax. CAD-Waxx in DISC format has an external diameter of 98mm (including circumferential groove) and a height of 20mm for open CAD/CAM Systems. Clear, transparent acrylic gives excellent visibility for easy implant placement. Residue free.

Vident Order No.	Description	Qty	Price
V53-2001	CAD-Waxx DISC-98mm x 20mm	1 piece	\$100.00*

Orders: 800-828-3839 USA | 800-263-4778 Canada | www.vident.com



VITABLOCS® 3D-Master Guide

VITABLOCS 3D-Master Guide . \$103.00 Order No. GX69-920...

VITABLOCS Guide 3D-Master

VITA YZ Finishing Kit

This unique kit features rotating instruments that were created specially to address the needs of zirconium oxide materials. The instruments cut easily without pressure, and create a smooth surface that minimizes the stresses often found in these materials when ground upon.



Order No. B31-001	\$189.00
Round end taper abrasive	
Order No. B31-101	\$40.20
Disc abrasive	
Order No. B31-102	\$46.40
Barrel abrasive	
Order No. B31-103	\$46.40
Inverted cone abrasive	
Order No. B31-104	\$40.20
Round end Cylinder diamond	
Order No. B31-105	\$45.40
Round End Chamfer diamond	
Order No. B31-106	\$45.40

Vident EASY Polishing Sets

The unique EASY polishing system offers efficient and effective finishing of all CAD/CAM ceramics in just 2 steps. These high performance diamond grinders and polishers eliminate the need for a pre-polishing step.



Easy Polishing Sets

A POPULO

For MOD filings and single, freestanding teeth

Order No. B24-100	. \$219.00
EASY Polishing Set, Right Angle (4 grinders & 6 polishers) Order No. B24-200	. \$182.00

EASY Polishing Set. Hand Piece (6 grinders & 6 polishers)

EASY Grinders ar	nd Polishers, Hand Piece\$18.00
Order No.	Description
B24-110	Blue Coarse Cone Grinder, HP
B24-111	Blue Coarse Cylinder Grinder, HP
B24-112	Blue Coarse Wheel Grinder, HP
B24-120	Pink Medium Cone Grinder, HP
B24-121	Pink Medium Cylinder Grinder, HP
B24-122	Pink Medium Wheel Grinder, HP
B24-130	Pink Medium Knife Edge Polisher, HP
B24-131	Pink Medium Flame Polisher, HP
B24-132	Pink Medium Wheel Polisher, HP
B24-140	Grey Fine Knife Edge, HP
B24-141	Grey Fine Flame Polisher, HP
B24-142	Grey Fine Wheel Polisher, HP
EASY Grinders a	nd Polishers, Right Angle\$18.00
Order No.	Description
B24-210	Blue Coarse Knife Edge Grinder, RA
B24-211	Blue Coarse Cone Grinder, RA
B24-220	Pink Medium Knife Edge Grinder, RA
B24-221	Pink Medium Cone Grinder, RA
B24-230	Pink Medium Flame Polisher, RA
B24-231	Pink Medium Cup Polisher, RA
B24-232	Pink Medium Knife Edge Polisher, RA
B24-240	Grey Fine Flame Polisher, RA
B24-241	Grey Fine Cup Polisher, RA

VITA Akzent Glaze Spray

Easy to apply, spray ceramic powder for simple, time-saving glazing of all ceramic restorations with a sintering temperature of 850°C.

Grey Fine Knife Edge Polisher, RA

Akzent Glaze S	pray 125ml	
Order No. S28-	790	\$122.00
Patterson Orde	r No. 706-9628	

Accu-Spray™

B24-242

Reflective spray scan powder specifically designed for lab use with a non-clinical propellant. Accu-Spray enhances the reflective properties of a prepared dental model fabricated with dental stone or scan stone. 200 ml can. Order No. V55-120...

Accu-Liquid™

A liquid version of scan powder that dries to a reflective powder prior to scanning. 10ml bottle specifically for lab use. Accu-Liquid is perfect for scanning implant abutments. Order No. V55-123...

Accu-Wax™

This highly reflective wax provides increased laser-scanning accuracy. Perfect for blocking out undercuts, voids, defects and recontouring die stone prior to scanning. 50g Tin. Order No. V55-122



^{*} Quantity Discount Pricing available. Annual Purchase Contracts available. Contact Vident for additional information.

Accu-Scan Plus™ Universal Scan Stone

Accu-Scan Plus is an exceptional resin-reinforced die stone for conventional and electronic laser scanning applications. This Class IV low expansion, high density die stone produces highly-detailed models for accurate fitting restorations.

Order No. B43-225. . \$101.00

VITA CEREC Powder

CHAIRSIDE PRODUCTS

White contrast powder to eliminate reflections for optical impressions.

Order No. V55-101 Patterson Order No.707-9925

VITA CEREC Spray Head

Application rod can be rotated 360°

Order No. V55-103 . \$56.75 Patterson Order No.707-9940

VITA CEREC Liquid

Temporary adhesive for the CEREC Powder. Remove with water 6ml

Order No. V55-104 . \$42.30 Patterson Order No. 707-9957

VITA CEREC Propellant

CFC-free propellant for applying CEREC Powder to the tooth with the spray head, 250ml

Order No. V55-106. . \$49.50 Patterson Order No. 707-9957



VITA Powder Scan Spray

VITA Powder Scan Spray is a blue colored pigment suspension that perfectly adjusts the optical reflective properties of tooth surfaces for optical scanning. The new formula features a high concentration of powder with high yield for greater efficiency, 75ml

Order Number: V55-117 Patterson Order No. 708-0138

Dentaco Scan Film

An alternative to powder and propellants, Scan Film is applied as a liquid and allowed to dry into a homogeneous powder surface. 10 micron film thickness. 5 ml can.

Order No. V55-113. \$220.00

Patterson Order No. 708-0070

Unit Dose 0.01ml packets. 60 ea.

Order No. V55-114.

Patterson Order No. 708-0088



\$220.00

VITA YZ Coloring Liquid stains unsintered YZ substructures from bright white to a range of 3D-Master and Classic shades. The Light and Medium stains will obtain the base shade while the Intense and Neutral stains will increase or decrease the intensity of the base shade.

VITA YZ Coloring Liquid Set Order No. V59-003

.\$**FPO**

YZ Coloring Liquid Light, 100ml

YZ Colorina Liquid Medium, 100ml

YZ Coloring Liquid Intense, 100ml

YZ Coloring Liquid Neutral, 100ml

Work containers

Acrylic tweezers for Coloring Liquid

YZ Coloring Liquid

VITA YZ Coloring Liquid Reorder Materials

YZ Coloring Liquid, Light, 250ml
YZ Coloring Liquid, Medium, 250ml
YZ Coloring Liquid, Intense, 100ml
YZ Coloring Liquid, Neutral, 100ml
Shade Indicator

VITA Shading Paste

These fluorescent fine-grained staining pastes are self-glazing and designed specifically for CEREC porcelain. Ready-to-use paste system does not require mixing a powder and a liquid

VITA Shading Paste Classic Kit

Order No. V55-500. \$835.00

Shading Pastes, 4g

VITA Glaze Paste, 7g

VITA Stain Liquid, 15ml

Shade Sample Disc

Mixing Palette

Glass Spatula

Stain Brush with Tips

Flat Squirrel Hairbrush

VITA Shading Paste

Classic Reorder Materials

Shading Paste	Stain Refills, 4g	\$66.00
Order No.	Shade	
V55-501	White	
V55-502	Orange	
V55-503	Ochre	
V55-504	Khaki	
V55-505	Olive	
V55-506	Medium Brown	
V55-507	Reddish Brown	
V55-508	Mauve	

VITA Glaze Paste, 7g

Order No V55-515	. \$111.00
VITA Stain Liquid, 15ml	
Order No. V55-520	\$45.40

VITA Shading Paste 3D-Master Kit

Order No.	vX55-006\$925.00	J
9	Shading Pastes, 4g	
1	VITA Glaze Paste, 7g	
1	VITA Akzent Finishing Agent, 5g	
1	VITA Shading Paste Liquid, 15ml	
1	Shade Sample Blade	
1	Mixing Palette	

Glass Spatula Stain Brush with Tips

Flat Squirrel Hairbrush Working Instructions

VITA Shading Paste 3D-Master

Reorder Materials

Shading Paste S	Stain Refills, 4g\$66.00
Order No.	Shade
VX55-521	Birch
VX55-522	Mellow Yellow
VX55-523	Curry
VX55-524	Fumo 1 (Sand)
VX55-525	Fumo 2 (Light Brown)
VX55-526	Fumo 3 (Medium Brown)
VX55-527	Tabac (Tobacco)
VX55-528	Shak (Gray)
VX55-529	Niagara (Blue)

VITA VM9 Esthetic Kit for VITABLOCS

Order No. NX98-0016 . \$ 1.435.00 Patterson Order No. 707-1970

Sample set of VITABLOCS 3D-Master

WINDOW WIN, 12g

NEUTRAL NT, 12q

ENAMEL ENL. END. 12a

EFFECT PEARL EP1, 12q

CORRECTIVE COR1, 12 g

AKZENT Finishing Agent Akz 25, 5g

EFFECT ENAMEL EE1, EE10, 12g

Shading Paste Glaze SP15, 7g

EFFECT OPAL EO2, 12g

EFFECT CHROMA EC1, EC4, 12g

MAMELON MM2, 12q

VITA VM Modelling Liquid, 50ml VITA Shading Paste Liquid, 15ml

Firing Travs

Platinum Pins, 2 pcs

Working Instructions

VITA Extreme Esthetic Characterization Kit

Order No. NX98-0017

\$1,045.00

VITA Akzent (AKZ01, AKZ02, AKZ03, AKZ04, AKZ09,

AKZ12, AKZ16, AKZ17), 3g each

Akzent Finishing Agent (AKZ26)

Akzent Glaze (AKZ25), 5a

VITA Akzent Fluid, 20ml

VM9 Window (WIN), 12g VM9 Enamlel Light (ENL), 12g

VM9 Corrective (COR1), 12g

VM9 Neutral (NT), 12g

VM9 Effect Enamel (EE9), 12a

VM9 Effect Enamel (EE1), 12a

VITA Glass Mixing Rod

Akzent Trav

Stain Brush (#0)

Quick-Peg[™] Custom Firing Pegs

Quick Peg is a paste that is squeezed into/onto the restoration, and placed onto the peg or firing tray. It dries guickly, holds the restoration in place during firing and is removed easily. It can be used with allceramic or porcelain-fused-to-metal restorations. It can also be used for degassing and staining and glazing. It is especially helpful in eliminating the frustrating balancing act experienced placing the crown onto the peg. With over 60 applications per 10cc syringe, Quick-Peg is economical enough to use with every type of restoration.

Quick—Peg 3-pack of 10cc syringes Order No. A47-320.. . \$35.50

VITA CEREC Etch

5% hydrofluoric acid gel for etching ceramic restorations. Not for use etching tooth structure.

Order No. V55-105. \$**FPO** Patterson Order No. 707-9965



. \$ TBA

VITA Simulate Preparation Material — a light curing composite — is used for simulating the shade of the prepared tooth that is being restored. This material is used as an aid in all-ceramic restorative cases in which the final shade of the restoration may be strongly influenced by the shade of the prepared tooth stump. The simulate material makes it easier for the dental technician or CEREC doctor to control the final shade of the restoration and to correct it if necessary before the restoration is tried in with the patient.



VITA Simulate Preparation Color Guide

VITA Simulate® Preparation Material Kit

Order No.	KX43-0001	\$345.00
6	VITA SIMULATE Preparation Material, 4 g	
	0M1S, 1M1S, 2M3S, 3M2S, 4M3S, 5M3S	
1	VITA SIMULATE Insulation Liquid, 16 ml	
150	Application Sticks	
1	VITA SIMULATE Preparation Material Guide	
	(Shade indicator)	
1	Instructions for use	

VITA Simulate® Preparation Material

VITA Simulate	Preparation Mate	rial Powder, 4g	\$56.75
Order No.	Shade	Order No.	Shade
KX43-4031	0M1	KX43-4046	3M2
KX43-4034	1M1	KX43-4054	4M3
KX43-4040	2M3	KX43-4059	5M3
	Preparation Mate	rial Liquid, 16ml	\$45 40
	Applications stick		ψ 10. 10
	• •		\$42.30
	Preparation Color		
Order No. GX7	0-430		\$73.25

Easymatch™ Software Shade Accuracy from Prep to Placement

Easymatch Software saves time and money by predicting the optimal VITA 3D-Master shaded VITABLOCS before placing a milling block into your CEREC or Sirona inLab device. With the new Easymatch software it is easy to predict difficult shade variables that could lead to remakes. Easymatch software displays the best shade solution based on critical shade parameters such as ceramic thickness, prep shade, and prescribed shade. Easymatch works with VITABLOCS Mark II, TriLuxe and TriLuxe Forte.

Use Easymatch software with the VITA Easyshade for highly precise, predictable shade outcomes every time.

Order No. GX10-001	\$795.00
Patterson Order No. 706-8992	



Orders: 800-828-3839 USA | 800-263-4778 Canada | www.vident.com

VITA Physiodens® Anterior and Posterior Teeth

VITA Physiodens teeth are designed to provide the highest esthetics and natural tooth morphology. Each anterior tooth is created using 3 layers of material which provides exciting nuances that provide light reflectiveness matching natural teeth. This multi layering provides opalescence, luminescence and balanced translucency. All teeth are hand made ensuring that each tooth is unique and provides the subtle characterization needed to enhance optimum esthetics. The teeth are made from a ceramic filled (16%) hardened acrylic material which provides strength, chip resistance and wear resistance for long term function. Since Physiodens teeth are the same size as natural teeth this makes the teeth ideal for partials. There is exceptional shade accuracy even after grinding the Vitapan tooth to a shell like thickness.

Physiodens teeth are available in 21 anterior moulds, 8 lower anterior moulds, 6 upper and lower posterior moulds. All Physiodens teeth come in 26 3D-Master, 15 VITA Classical and two bleached shades.



Dentures by Dennis Purinton, CDT

VITA Vitapan® Anterior & Posterior Teeth

Ergonomic, simple, easy to set up, and natural esthetic separate Vitapan from other premium teeth. The hand made 3 layer esthetic teeth provide for individual characterization and a high degree of translucency, florescence and opalescence. The Vitapan material is the same high quality ceramic filled hardened acrylic as the VITA Physiodens. This unique formula provides a high strength, shock and impact resistant tooth which resists chipping even when grinding.

Vitapan teeth offer a large selection of anteriors and posteriors teeth to meet everyday needs. There are 31 anterior uppers and 13 lower anteriors. The lower anteriors are manufactured either in full or abraided anatomy to mimic natural wear. The posteriors are available in either Cuspiform, with a 23-28% cusp angulation or Synoform, with a 5% angulation to meet individual patient needs. The Synoforms are specially made where there is a low degree of intercuspation and their narrow design provides for more tongue room for aging patients.

All Vitapan teeth are available in 26 3D-Master, 15 VITA Classical and two bleached shades.

Synoform Posteriors exhibit the same material physical properties as Vitapan Anteriors and have shallow inter-locking occlusal surfaces.

Synoform teeth are particularly indicative where there is a low degree of intercuspation. Their narrow design provides for more torque room for the aging patient.

- Vitapan Anterior 1 x 6
- Vitapan Posterior 1 x 8
- Vitapan Plus Anterior 1 x 6
- Physiodens Anterior 1 x 6
- Physiodens Posterior 1 x 8
- Lumin Vacuum Porcelain Anterior 1 x 6
 Lumin Vacuum Porcelain Posterior 1 x 8
- Lingoform 1 x 8
- VITA MFT 1 x 6
- VITA MFT 1 x 8

VITAPAN® PLUS Anterior Teeth

VITA has modernized the classical VITAPAN and given it a new, future-oriented design. VITAPAN PLUS is a superior anterior tooth that combines a natural appearance and ideal function. The distinct horizontal and vertical curvature characteristic of the tooth provide optimum support of the lip area and offers a balanced anterior set-up. The anatomy of the tooth matches the design of the palatal surfaces for conserving the patient's phonetics and accustomed tongue sensation. With life-like layering, VITAPAN PLUS features noticeably enhanced esthetics with natural shade, light and reflective properties. To facilitate better dental hygiene and for more realistic gums, the broadened tooth neck area features interdental spacing and a natural papilla form, which offers several possibilities for age-appropriate prosthesis. Produced from excellent material quality, VITAPAN PLUS consists of homogeneously uniform VITA MRP (Microfiller Reinforced Polyacrylic) material throughout all layers and features an outstanding degree of wear resistance.

VITAPAN PLUS is available in 10 upper anterior and 5 lower anterior moulds and 10 Classical and 3D-Master shades.





DENTURETEETH

VITA TOOTH DEALERS

lowa Dental Supply	(800) 362-2072
Nowak Dental Supplies, Inc	(800) 654-7623
The Teeth Connection, LLC	(808) 722-8511
Arnold Dental Supply	(800) 562-6645
Dorado Dental Supply, Inc	(619) 281-4337
Hi-Tec Dental Products	(800) 859-2006
Zahn Dental	(800) 496-9500
Henry Schein	(800) 668-5558
Sinclair	(800) 663-7393
DenPlus	(888) 344-4424
Central Dental	(800) 268-4442
Patterson	(800) 570-6356
Westan	(204) 946-5452

VITA LINGOFORM Posterior Denture Teeth

The aim of a lingualized occlusion denture set up is the stabilization of the prosthetic restoration. The lingualized occlusion concept reduces lateral excursion loading which could minimize off axis stress to the underlying bone structure. In implant supported overdentures, this may minimize destructive lateral forces on the implant.

The LINGOFORM posterior teeth makes it possible to combine excellent function and superior esthetics. The occlusal surfaces in the mandibular teeth have a buccal cuspal inclination of 15° and 20° in the lingual area. In the upper jaw the somewhat deeper occlusal surface relief comes very close to nature without sacrificing function. The ridge lap has been reduced to facilitate faster set up and reduce grinding. The VITA LINGOFORM posteriors are indicated for:

- complete denture prosthetics
- implant supported overdentures
- attachments

48

partial dentures

VITA LINGOFORM posteriors can be combined with VITA PHYSIODENS or VITAPAN anteriors as required.

VITA LINGOFORM posterior teeth sets are exclusively available in 26 3D-MASTER shades, 15 VITA CLASSICAL shades and 1 Bleached shade. LINGOFORM is available in three upper and lower moulds 21L, 22L, 23L and 24L.

VITA® MFT (Multi Functional teeth)

The VITA MFT teeth are a new tooth concept designed to provide excellent function as well as esthetic in a medium priced tooth. The design principle was to create harmony with nature. VITA MFT teeth are indicated for full and partial dentures.

The true-to-nature appearance of the anterior teeth, built up in three layers, meet all the criteria of a restoration designed to satisfy the patients needs:

- Opalescence, luminescence and translucence combined with perfectly matching tooth morphology and an individual mamelon structure, create a play of rich colors in anterior nuances as well as a vivid light reflection in every single tooth.
- The slightly convex design of the labial facettes purposefully enhances the lip contour and thus promotes dentofacial harmony.
- The natural reproduction of the palatal surfaces assist every patient in maintaining their familiar phonetic skills and tongue contact.
- As a further highlight, the white incisal margins support the outstanding esthetic of the VITA MFT anteriors.

It is the unique design of the posteriors that gives VITA MFT its name. This multifunction teeth posterior design lends itself to various occlusal schemes with ease, including lingualized occlusion. The occlusal design is based on the innovative easy-centric principle, i.e. the occlusal surfaces are designed in such a way that, almost automatically, the perfect centric is reached.

- The unique pre-abraded occlusal surfaces of the lower posteriors is designed to minimize horizontal shear forces. This is particularly important with implant supported overdentures.
- The ridge lap has been reduced to minimize the need for grinding and enlarges the surface for safe bonding to the denture base material.

VITA MFT is available in 13 Shades, 0M1, 0M3, 1M1(B1), 2M1 (A1), 2L1.5 (B2), 2M2 (A2), 3L1.5, 3L2.5 (B3), 3M2 (A3), 3R2.5 (A3.5), 4L1.5, 4M2, 5M1, (A4) and 17 anteriors and 3 posterior moulds.

The whole range of indications in a pleasantly small VITA MFT product line. A reduced and harmonic range of shades/moulds makes for increased economic efficiency.

Multifunctional VITA MFT teeth are suited for:

- Full Dentures
- Partial Dentures
- Implant Prosthetics
- Combination Cases



Denture Teeth Mould Guides

Order No.	Description	
D85-020	VITAPAN Anterior Classical, 44 x 6	. \$575.00
D85-040	VITAPAN Anterior A1, 44 x 6	. \$575.00
D85-045	VITAPAN Anterior A2, 44 x 6	. \$575.00
DX85-023	VITAPAN Lingoform Posterior 3D, 12 x 8,	
	10 x 8, 6 x 8	. \$575.00
DX85-040	VITAPAN Anterior 1M1, 44 x 6	. \$575.00
D85-042	Physiodens A1, 29 x 6, 12 x 8	. \$575.00
DX85-022	Physiodens 1M1, 29 x 6, 12 x 8	. \$575.00
DX85-021	Physiodens 2M2, 29 x 6, 12 x 8	. \$575.00
DX85-027	Physiodens/Lingoform 3D-Assorted,	
	29 x 6, 20 x 8	. \$575.00
DX85-030	MFT Anterior & Posterior 2M2	
	17 x 6, 6 x 8	. \$200.00
DX85-031	MFT Anterior & Posterior 2M2	
	17 x 6, 6 x 8	. \$155.00

Denture Cabinets

Order No.	Description
D85-030	Denture cabinet, 2 drawer inserts\$100.00
D85-031	Denture cabinet, 4 drawer inserts\$150.00
D85-032	Denture cabinet, 12 drawer inserts\$380.00

Drop Dead Gorgeous Dentures (DVD or VHS)

This innovative video, featuring Dr. Charles Barotz, demonstrates how to think about, market, sell and fabricate amazingly lifelike dentures.

Order No.	Description	
A47-627	(DVD) or A47-628 (VHS)	\$265.00

BIO-Logical Prosthetics DVD (By Dr. Eugen End)

BIO-Logical Prosthetics is an occlusion concept according to nature's example. This interactive DVD-ROM gives a demonstration of BIO-Logical Prosthetics and its application in the field of complete denture prosthetics.

BIO-Logical Prosthetics DVD	
Order No. A47-630\$	175.00

VITACOLL™ Bond

VITACOLLTM ensures reliable chemical bonding of acrylic teeth to hot-curing and cold-curing resins as well as polymethacrylate injection-molding materials, thus avoiding breakage of the teeth out of the base. The basal surface of the teeth should be roughened prior to placement of the VITACOLL material.

VITACOLL Bond, 100ml	
Order No. D85-251	\$40.00

VITAFOL® H Paste

and saves time during the processing of prefabricated teeth. When VITAFOL H is used, the teeth are protected during devesting, interdental spaces remain clean and gingival characterization is retained.

VITAFOL H with adhesive crystals ensures reliability

VITATUL Lab Kit, 70IIII	
Order No. D85-252	\$42.00





VITACOLL

VITAFOL H Paste

VACUMAT 6000M Ceramic Furnaces

VITA presents a first in the dental industry; one-of-a-kind VITA VACUMAT 6000M modular ceramic furnaces. The unique and innovative concept satisfies any laboratory's firing needs in an efficient and economical way. The heart of the system is the VITA VACUMAT 6000M furnace with its four unique control units, the vPad easy, the vPad comfort, vPad Clinical and the vPad excellence. These feature-filled controllers enable up to four furnaces to be controlled with one control unit. Buy one controller with a firing unit now, then add additional firing units inexpensively as needed.

The new VITA VACUMAT 6000M is a premium microprocessor-controlled, fully automatic furnace used to fire all types of veneering and crystallization ceramics. The unit features a modern, compact and ergonomic design and just like its predecessors, requires considerably less bench space in comparison to conventional ceramic furnaces.

In addition, each of the VITA VACUMAT 6000M firing units can be individualized with colored side panels. The white panels can also serve as a notepad.

A choice of three different icon-driven control units with intuitive symbols are available to program and operate the furnace.

- **The vPad** easy has 200 user-adjustable programs, all the required functions and produces excellent firing results.
- The vPad Comfort has 500 user-adjustable programs, a color touch-screen, a Photo Viewer and can control one or two furnaces simultaneously.
- The vPad Clinical is matched with the requirements of CEREC dentists, featuring a 7-inch touch display that offers clearly arranged menus and enables fast and intuitive operation. All necessary block programs are pre programmed. Up to 2 furnaces can be conveniently controlled.
- The vPad Excellence control unit offers 1,000 customizable programs, is equipped with an 8.5 inch color touch-screen, a Photo Viewer, can control up to four furnaces simultaneously, has individual user accounts and many other useful extras.

VACUMAT 6000M Firing Units

VACUMAT 6000M Firing Unit	
Order No. E60-001	\$2650.00
VACUMAT 6000M Firing Unit - Stainless Steel	
Order No. E60-002	\$3250.00



VACUMAT 6000M Controllers

vPad Easy	
Order No.E60-003	\$850.00
vPad Comfort	
Order No.E60-004	\$1695.00
vPad Clinical	
Order No.E60-040	\$3195.00
vPad Excellence	
Order No.E60-005	\$3195.00

VACUMAT 6000M Vacuum Pumps

Vac 6000M Vacuum Pump 110V	
Order No. E60-006\$	3795.00

VACUMAT 6000M Firing Accessories

Vac6000M Metal Cooling Tray Set	
Order No. E60-007	\$89.00
Side Panel Set, Blue	
Order No. E60-008	\$89.00
Side Panel Set, Pink	
Order No. E60-009	\$89.00
Side Panel Set, Mint Green	
Order No. E60-010	\$89.00
Side Panel Set, Apricot	
Order No. E60-011	\$89.00
Side Panel Set, White	
Order No. E60-012	\$89.00
Side Panel Set, Dark Gray	
Order No. E60-013	\$89.00
Number Set	
Order No. E60-014	\$39.00
Switchbox	
Order No. E60-015	\$379.00
Switchbox Cable	
Order No. E60-016	\$29.00

In-Ceramat III® T Touchscreen Furnace

The Inceramat "high temp" sintering furnace is used for the In-Ceram technique as well as all prolonged high temperature sintering and glass infiltration requirments.

Features/Benefits

- High-heat quartz spiral muffle (rated 1200°c)
- 100 user-modifiable programs (10x In-Ceram pre-installed)
- Automatic calibration at power up
- Firing table moves forward for easy access
- Fold-away twin cooling table
- Stainless steel construction
- Supplied complete with 2nd muffle as back up
- 2 year warranty on parts and labor
- 1 year warranty on muffle

Easyglaze Furnace

The Easyglaze Furnace is ideal for any dentist or lab that is staining, glazing, or performing any firing function that does not require vacuum. Fully automatic, precise and affordable, the Easyglaze Furnace offers a roomy muffle, a 3 inch diameter firing tray and a microprocessor capable of holding temperatures to within ± 2 degrees C. The accurate muffle is wound alloy with high temperature fibrous alumina insulation.

Features/Benefits

- 2 Preset programs for typical VITA and CEREC Glazes
- 3 Customizable user programs
- Automatic Re-glaze Feature
- Operating temperature up to 1000°C
- Controlled temperature ramp-up rate (60-120°C/min.)
- Clamshell design allows the firing tray to remain stationary
- Small size, does not take up much bench-top space
 Simple and user-friendly interface and controls
- Large and clear LCD display
- ETL/CE safety markings
- Retractable cooling tray
- 1 year warranty

Order No. E45-001 \$1,795.00



In-Ceramat III T Touchscreen Furnace



Program example



Easyglaze Furnace



\$11.60

The VITA Zyrcomat T is a completely new programmable atmospheric high-temperature sintering furnace for VITA In-Ceram YZ and In-Ceram AL cubes, and is recommended for use with the Sirona inLab system. Time-saving sintering features include fifty adjustable programs that are easily accessed using the new large color touchscreen.

Features/Benefits

FURNACES

- High-heat firing chamber (rated 1600°C)
- Furnace and control (power) unit is separate for convenient shipping
- Attractive and robust casing consisting of painted sheet steel and stainless
- Processing capabilities of 75 + plus/single/mixed unit restorations (Stacked)
- 5 molybdenum elements for uniform heat distribution (4 with 230V)
- Typical processing time including cooling and removal 5.45 hours
- Voltage loss protection with restart recovery
- 2 year warranty parts and labor (6 months on heating elements)
- 110 voltage (with optional 230v delivery)
- Total weight 130.9 lbs (57 kgs)

Accessories included:

- Furnace tweezers
- Sintering crucible
- Sintering dish
- 1 pack ZrO2 sphere pearls
- Allen wrenches
- Operating instructions
- Automatic calibration at power up
- Upgradeable software via USB
- Multiple user languages
- 40 programmable programs

Order No. E59-014. .. \$11.990.00

VITA_® **Zyrcomat T Furnace Accessories**

Order No.	Description	
38002	Sintering Pearls	\$128.00
38010	Sintering Crucible Short	\$585.00
38006	Sintering Dish	\$535.00

Super Steam Cleaner

The Vident Super Steam Cleaner provides continuous demineralized steam cleaning. Its large tank capacity allows for extended use between fillups.

- Simple to use and install.
- Fitted with three pressure safety devices
- Large half gallon tank capacity
- For cleaning models through to completed work
- Degrease of frameworks prior to ceramics
- Cleaning attachments and implant components
- De-waxing denture teeth
- 70-75 PSI Pressure range
- 1 year warranty on parts and labor

Order No. E68-001. .\$1,290.00



Poly-Die[™] and Poly-Base[™] Material

Poly-Die and Poly-Base polyurethane material are specifically designed for dental laboratories. These products can be used with any type of impression material and many different styles of models, from regular pin to tray and plastic plate systems.

Poly-Die is ideal for implants, orthodontic models and partial denture work and can be used in combination with regular stone when using dowel pins, providing superior cost savings. The die material features no contraction or expansion, eliminating the fit problems sometimes encountered with

typical die stone. Poly-Die is odorless.

Poly-Base is a less-expensive model base material that shares many of the physical properties of Poly-Die. It is not quite as dimensionally accurate, and as such, is ideal as a base material. Due to the inclusion of advanced resin technologies, Poly-Die and Poly-Base weigh 40% less than other materials, producing more models with less consump-

Order No.	Description	
A95-090	Poly-Die Kit, 1000g	\$110.00
A95-091	Poly-Base Kit, 1800g	\$162.00
A95-018	Separating Spray (also for Klic-Lok)	\$21.70

Sta-Put[™] Adhesive for Models

Sta-Put Gel is used to temporarily hold the restoration in place on the model while it is being worked on. This low-viscosity gel-like material is squeezed into the restoration or onto

Once placed, the porcelain is held in position while the dentist or ceramist adds stain, glaze, contacts, creates characterizations, etc. It burns out clean when heated.

Sta-Put Adhesive, 3-pack of 3ml syringes	
Order No. A89-400\$47.4	40

Vident Diamond Restoration Polish, Jar

Specially formulated for restoring the glaze finish on ground porcelain or composite surfaces. Adds a smooth mirror finish to all high-noble, noble and base metal alloys.

Order No.	C15-100	\$68.00
1	Jar Diamond Restoration Polish, 3g	

- Mandrel
- Brush Wheel
- Felt Wheels

Diamond Restoration Polish, Syringe, 3g

Syringe delivery system for convenient dispensing of diamond impregnated polishing paste. Ideal for polishing ceramics and composite materials. Kit contains one 3g syringe, straight handpiece wheel brush, and three felt wheels.

Order No. C15-101 \$64.0	00
--------------------------	----

Technical Help Desk: 888-249-1640 USA | 800-324-6224 Canada

Wax-Up Assistant from Smile Line

With the Architect Wax-up Assistant, technicians have the ability to create accurate diagnostic wax-ups using a simple, fast technique. After selecting the correct wax pattern (available in 3 sizes), position the pattern on the die using positioning wax, fill and then finish with modeling wax.

All patterns are available in white only.

Wax-up Assistant, Anterior Set

Order No. CS-2010W	\$545.00
Wax-up Assistant, Posterior Set	
Order No. CS-2020W	\$545.00
Wax-up Reorder Materials	
Smile Line Modeling-wax 60g Order No. CS-2050W	\$65.25

Reorder Materials

Unner Anteriors

Opper Amenors	
Order No.	Description
CS-A11SW	Smile Line 6 pc. Refill in A11S
CS-A11LW	Smile Line 6 pc. Refill in A11L
CS-A12SW	Smile Line 6 pc. Refill in A12S
CS-A12MW	Smile Line 6 pc. Refill in A12M
CS-A12LW	Smile Line 6 pc. Refill in A12L
CS-A13SW	Smile Line 6 pc. Refill in A13S
CS-A13MW	Smile Line 6 pc. Refill in A13M
CS-A13LW	Smile Line 6 pc. Refill in A13L
CS-A21SW	Smile Line 6 pc. Refill in A21S
CS-A21MW	Smile Line 6 pc. Refill in A21M
CS-A21LW	Smile Line 6 pc. Refill in A21L
CS-A22SW	Smile Line 6 pc. Refill in A22S
CS-A22MW	Smile Line 6 pc. Refill in A22M
CS-A22LW	Smile Line 6 pc. Refill in A22L
CS-A23SW	Smile Line 6 pc. Refill in A23S
CS-A23MW	Smile Line 6 pc. Refill in A23M
CS-A23LW	Smile Line 6 pc. Refill in A23L
CS-B11SW	Smile Line 6 pc. Refill in B11S
CS-B11MW	Smile Line 6 pc. Refill in B11M
CS-B11LW	Smile Line 6 pc. Refill in B11L
CS-B12SW	Smile Line 6 pc. Refill in B12S
CS-B12MW	Smile Line 6 pc. Refill in B12M
CS-B12LW	Smile Line 6 pc. Refill in B12L
CS-B13SW	Smile Line 6 pc. Refill in B13S
CS-B13MW	Smile Line 6 pc. Refill in B13M
CS-B13LW	Smile Line 6 pc. Refill in B13L
CS-B21SW	Smile Line 6 pc. Refill in B21S
CS-B21MW	Smile Line 6 pc. Refill in B21M
CS-B21LW	Smile Line 6 pc. Refill in B21L
CS-B22SW	Smile Line 6 pc. Refill in B22S
CS-B22MW	Smile Line 6 pc. Refill in B22M
CS-B22LW	Smile Line 6 pc. Refill in B22L
CS-B23SW	Smile Line 6 pc. Refill in B23S
CS-B23MW	Smile Line 6 pc. Refill in B23M
CS-B23LW	Smile Line 6 pc. Refill in B23L



54

. \$55.75

CS-C11SW CS-C11MW CS-C11LW CS-C12SW CS-C12MW CS-C12LW CS-C13SW CS-C13MW CS-C13LW CS-C21SW CS-C21LW CS-C22LW CS-C22SW CS-C22SW CS-C22MW CS-C22SW CS-C23WW CS-C23SW CS-C23MW CS-C23MW CS-C23MW	Smile Line 6 pc. Refill in C11S Smile Line 6 pc. Refill in C11M Smile Line 6 pc. Refill in C11L Smile Line 6 pc. Refill in C12S Smile Line 6 pc. Refill in C12M Smile Line 6 pc. Refill in C12M Smile Line 6 pc. Refill in C12L Smile Line 6 pc. Refill in C13S Smile Line 6 pc. Refill in C13M Smile Line 6 pc. Refill in C13L Smile Line 6 pc. Refill in C21S Smile Line 6 pc. Refill in C21M Smile Line 6 pc. Refill in C21L Smile Line 6 pc. Refill in C22L Smile Line 6 pc. Refill in C22L Smile Line 6 pc. Refill in C22L Smile Line 6 pc. Refill in C23S Smile Line 6 pc. Refill in C23M Smile Line 6 pc. Refill in C23M Smile Line 6 pc. Refill in C23L
	\$11.60
Order No. CS-I31SW	Description Smile Line 6 pc. Refill in I31S
CS-I31LW	Smile Line 6 pc. Refill in 131L
CS-I32SW	Smile Line 6 pc. Refill in I32S
CS-I32MW	Smile Line 6 pc. Refill in I32M
CS-I33SW	Smile Line 6 pc. Refill in I33S
CS-I33MW	Smile Line 6 pc. Refill in I33M
CS-I33LW	Smile Line 6 pc. Refill in I33L
CS-I41SW	Smile Line 6 pc. Refill in I41S
CS-I41MW	Smile Line 6 pc. Refill in I41M
CS-I41LW	Smile Line 6 pc. Refill in I41L
CS-I42SW	Smile Line 6 pc. Refill in I42S
CS-I42MW	Smile Line 6 pc. Refill in I42M
CS-I42LW	Smile Line 6 pc. Refill in I42L
CS-143SW CS-143LW	Smile Line 6 pc. Refill in I43S Smile Line 6 pc. Refill in I43L
G3-143LVV	Sillie Lilie o pc. Neilli III 143L
Upper Posteriors	\$11.60
Order No.	Description
CS-UP14SW	Smile Line 6 pc. Refill in UP14S
CS-UP14LW	Smile Line 6 pc. Refill in UP14L
CS-UP15SW	Smile Line 6 pc. Refill in UP15S
CS-UP15LW	Smile Line 6 pc. Refill in UP15L
CS-UP16SW	Smile Line 6 pc. Refill in UP16S
CS-UP16MW	Smile Line 6 pc. Refill in UP16M
CS-UP17SW CS-UP17LW	Smile Line 6 pc. Refill in UP17S Smile Line 6 pc. Refill in UP17L
CS-UP24SW	Smile Line 6 pc. Refill in UP24S
CS-UP24LW	Smile Line 6 pc. Refill in UP24L
CS-UP25SW	Smile Line 6 pc. Refill in UP25S
CS-UP25MW	Smile Line 6 pc. Refill in UP25M
CS-UP26SW	Smile Line 6 pc. Refill in UP26S
CS-UP26MW	Smile Line 6 pc. Refill in UP26M
CS-UP26LW	Smile Line 6 pc. Refill in UP26L
CS-UP27SW	Smile Line 6 pc. Refill in UP27S
CS-UP27MW	Smile Line 6 pc. Refill in UP27M

Lower Posteriors	\$11.60
Order No.	Description
CS-LP35MW	Smile Line 6 pc. Refill in LP35M
CS-LP35LW	Smile Line 6 pc. Refill in LP35L
CS-LP36SW	Smile Line 6 pc. Refill in LP36S
CS-LP36MW	Smile Line 6 pc. Refill in LP36M
CS-LP36LW	Smile Line 6 pc. Refill in LP36L
CS-LP44SW	Smile Line 6 pc. Refill in LP44S
CS-LP44MW	Smile Line 6 pc Refill in LP44M
CS-LP45SW	Smile Line 6 pc. Refill in LP45S
CS-LP45MW	Smile Line 6 pc. Refill in LP45M
CS-LP46SW	Smile Line 6 pc. Refill in LP46S
CS-LP46MW	Smile Line 6 pc. Refill in LP46M
CS-LP46LW	Smile Line 6 pc. Refill in LP46L
CS-LP47SW	Smile Line 6 pc. Refill in LP47S
CS-LP47MW	Smile Line 6 pc. Refill in LP47M
CS-LP47LW	Smile Line 6 pc. Refill in LP47L

Laservest Soldering Investment

Vident's Laservest Soldering Investment is a premium quality, high-heat investment material specially formulated for use in soldering high-noble, noble and base metal crown and bridge alloys.

- Does not give off destructive gases
- Fine particle size allows the investment to be mixed and poured easily
- Stable at high temperatures for both high and low fusing soldering operations
- Assured dimensional accuracy and fit

Ceramic Alloys

Vident offers a variety of superior quality ceramic alloys that are matched to the thermal expansion of VITA VM13, VMK 95, Omega and Omega 900 Porcelain Systems. All Vident alloys exhibit improved working characteristics, excellent casting and

finishing properties and consistent marginal integrity. Call for market prices and available discounts.

V13 Alloy

V13 is a high-strength ceramic alloy, ideal for single units through bridges, in both the anterior and posterior areas of the mouth. Made as a match for VITA VM13, it is also compatible with a wide range of porcelains with a CTE of up to 14.2. With 45% Au, 44.3% Pd and no silver or Cu, V13 is easy to work with, offers a high bond strength, and produces a high-gloss polish.

Palladius V+ Alloy

Palladius V+ is a micro-fine grain, silver-free PFM alloy with only 1.8% gold content. The alloy creates a silver-free environment that virtually eliminates the discoloration in porcelain. It is compatible with most porcelain systems, and is recommended for fabricating single crowns, short and long-span bridges and implant super-structures. 1.8% Au, 75.7% Pd, 8.0% In, 6.3% Ga.

Aurix Alloy

Vident's Aurix is a silver-free, micro-fine-grained alloy formulated to be compatible with most

dental porcelains. Its structure and superior physical properties facilitate the fabrication of all types of crowns, multi-units and bridge substructures. 66% Au, 25.5% Pd, 8.1% In, 0.3% Ga.

Solders and Fluxes for Vident's Noble and Base Metal Alloys

To complete our alloy systems, Vident offers a comprehensive line of pre- and post-solders, as well as fluxes in liquid and paste form.

Flux Paste for Palladius V+, 15g	
Order No. M20-129	\$22.50
Pre-solder for Aurix/Palladius V+	
Order No. M20-133	CALL
Post-solder for Aurix/Palladius V+	
Order No. M20-134	CALL
Pre-solder for Aurix AG	
Order No. M20-136	CALL
Post-solder for Aurix AG	
Order No. M20-137	CALL

Klic-Lok™ Model Tray System

The Klic-Lok Model Tray System introduces a new dimension in model-making. This unique single-pour re-usable system features a unique locking retentive feature, which retains the dies securely in a durable base tray with a "Klic". It is easy to use and offers attractive benefits, such as full arch and quadrant trays, and your choice of mounting methods, including magnetic mounting plates. You no longer have to compromise accuracy with conventional dowel-pin systems. You won't have to worry about breaking lower anterior dies when sawing your model. The Vident Klic-Lok Model System offers an outstanding professional-looking model with extreme accuracy, time savings, and ease-of-use.



Klic-Lok Model Tray System

Klic-Lok™ Tray Intro Kit

o. A95-001	\$153.00
Full arch trays	
Left spines for full arch and quad trays	
Right spines for full arch and quad trays	
Left quad trays	
Right quad trays	
Cutting fixtures	
Release fixture	
Strike rod	
O Left quad labels	
O Right quad labels	
Release spray	
	Full arch trays Left spines for full arch and quad trays Right spines for full arch and quad trays Left quad trays Right quad trays Cutting fixtures Release fixture Strike rod Left quad labels Right quad labels

Klic-Lok™ Tray Mini Intro Kit

1 Articulator sample pack

Order No. A	.95-007
20	Full arch trays
30	Left spines for full arch and quad trays
30	Right spines for full arch and quad trays
10	Left quad trays
10	Right quad trays
1	Release fixture
1	Strike rod
100	Left quad labels
100	Right quad labels

1 Articulator sample pack

Klic-Lok™ Full Arch Kit

Order No.	. A95-002	\$150.00
40	Full arch trays	
80	Left spines for full arch and quad trays	
80	Right spines for full arch and quad trays	
5	Cutting fixtures	
1	Release fixture	
1	Strike rod	
100	Left quad labels	

Klic-Lok™ Full Arch Mini Intro Kit

Articulator sample pack

Release spray

Order No.	. A95-004\$74.2	5
40	Full arch trays	
40	Left spines for full arch and quad trays	
40	Right spines for full arch and quad trays	
5	Cutting fixtures	
1	Release fixture	
1	Strike rod	
1	Articulator sample pack	

Klic-Lok™ Magnet Kit

	•	
Order No.	A95-042	\$44.30
10	Klic-Lok magnets	
3	Klic-Lok caps	



ACCESSORIES

Klic-Lok™ Magnet Plate Kit

Order No.	A95-041	\$14.5
10	Plates	
3	Sning covers	

Klic-Lok™ Reorder Items

0 (=)	
Cutting fixtures (5)	ф0.00
Order No. A95-009	\$6.20
Release fixture (1)	#0.40
Order No. A95-010	\$3.10
Spine covers (3)	# 40.00
Order No. A95-011	\$10.30
Strike rod (1)	Ф0.10
Order No. A95-012	\$3.10
Flex die cutter	ቀንቦ ንቦ
Order No. B30-101	\$39.20
Left quad labels (100pk)	ቀ ጋ7 10
Order No. A95-017	\$37.10
Right quad labels (100pk) Order No. A95-019	ቀ 27 10
Release spray (1)	ჶა/.10
Order No. A95-018	¢21 70
Full arch trays (10)	ΦΖ1.70
Order No. A95-025	\$7.30
Right spines for full arch/quad trays(10)	ψ7.50
Order No. A95-026	\$3.10
Left spines for full arch/quad trays (10)	φσ. το
Order No. A95-027	\$3.10
Right guad trays (10)	
Order No. A95-028	\$6.20
Left quad trays (10)	, -
Order No. A95-029	\$6.20
Magnet plates (25)	
Order No. A95-036	\$21.70
Magnet caps (100)	
Order No. A95-037	\$70.25
Magnets (10)	
Order No. A95-038	\$36.10

The Quick Fix Mandrel System For the quick release and attachment of discs

This new mandrel system does away with the tedious process of changing discs and polishers. The system consists of two different mandrels.

The Yellow mandrels are designed to be used with two-layered flexible diamond discs for shaping porcelain and larger diamond discs used to cut/section stone/polyurethane models. To mount them, you simply insert the pin into the mandrel and press until it engages with a click.

The Black mandrels are used with Vident's Golden Devil separating discs. These discs are truly superior in terms of durability and cutting efficiency and can be used on all metals. To change a disc, you just turn the pin 90° .

56

Quick Fix Kits/Accessories

QF - Hi-Tech Diamond Plaster Disc, 45mm	
Order No. B30-101	\$56.75
QF - Hi-Tech Diamond Plaster Disc, 35mm	
Order No. B30-102	\$50.50
QF - Mandrel, yellow (for Diamond Discs)	
Order No. B30-201	\$21.70
QF Diamond porcelain contour set,	
1 mandrel (yellow), 2 discs	
Order No. B30-003	\$50.50
QF Sepaflex Diamond contour Disc,	
22mm x 0.17mm, 3-pk	
Order No. B30-103	\$50.50
QF Plexoflex Diamond contour Disc,	
22mm x 0.22mm, 3-pk	
Order No. B30-104	\$50.50
QF Golden Devil Separating set,	
1 mandrel black, 50 discs	
Order No. B30-002	\$48.30
QF Golden Devil Separating Discs, (50)	
Order No. B30-105	\$29.40
QF - Mandrel, black (for Golden Devil Discs only)	
Order No. B30-202	\$21.70



CONTOURING

Sepaflex

Ceramics and composite

Grit size: fine Double-sided Diameter: 22mm Thickness: 0.17mm Order No. B30-103 Price: \$50.50



Plexoflex

Separating and contouring of ceramics and acrylics

Grit size: fine Double-sided Diameter: 22mm Thickness: 0.22mm Order No. B30-104 Price: \$50.50



SEPARATING DISCS

Golden Devil

Extra thin separating disc for metals

Diameter: 22mm Thickness: 0.24mm Order No. B30-105 Price: \$29.40



PLASTER DISCS

High Tech

Separating disc for gypsum and resins

Grit size: fine Double-sided Diameter: 35mm Thickness: 0.30mm Order No. B30-102 Price: \$50.50

Separating disc for plaster and synthetic resins

Grit size: fine Double-sided Diameter: 45mm Thickness: 0.30mm Order No. B30-101 Price: \$56.75



ArCon® Refills

Order No.	le Mounts A75-003 Disposable mounts	\$64.00
	le Mounts for Model Tray Systems A75-014 Disposable mounts for model tray systems	\$64.00
	=	\$52.75
	le Upper & Lower Incisal Holder (w/ Angle Base) A75-007 Disposable incisal pin holders Disposable incisal pin bases (angled)	\$52.75
Hinges, C Order No. 3 6 6	A75-002	\$47.40



Slider[™] Disposable Articulator

The inexpensive Slider disposable articulator simulates movement of the condyle by providing smooth, true protrusive and lateral movement. Compatible with most tray/model base systems.



Slider[™] Disposable Articulator Kits

Hinges (50) with Doweled Connector Plates (100)	
Order No. A95-055	20.60
Hinges (50) with Flat Connector Plates (100)	
Order No. A95-057	20.60

Vident Die Sealer Pen, 2 pk

A quick and economical way of providing a sealed surface for porcelain application.

Order No. A41-351.....\$45.00

VITA Modisol Pen

Provides clean and accurate separation of the build up from the plaster model with real application time savings.

Order No. A41-251.....\$58.75

Precision Porcelain Scoop Set



Modeling Liquids

ACCESSORIES

These universal modeling liquids have been developed to work with the majority of VITA ceramics. VITA Modeling Liquid improves the working and handling characteristics of your porcelains. Use VITA Modeling Liquid with VMK 95, VMK Master, Omega, Omega 900 and Vitadur Alpha. This fluid replaces Omega Modeling Liquid, Omega 900 Modeling Liquid, VMK 95 Modeling Liquid and Vitadur Alpha Modeling Liquid.

VITA VM Modeling Liquid	
50ml Order No. V68-909	\$22.70
250ml Order No. V69-909	\$78.50
Opaque Fluid	
50ml Order No. VX50-2320	\$21.70
250ml Order No. VX50-2321	\$78.50
Paste Opaque Liquid, 15ml	
Order No. V68-907	\$30.90
VITA Modeling Liquid RS, Red	
50ml Order No. V68-913	\$24.70
250ml Order No. V69-913	\$79.50
VITA Modeling Liquid	
50ml Order No. V68-750	
250ml Order No. V69-750	\$68.00
VITA Opaque Liquid	
50ml Order No. V68-751	
250ml Order No. V69-751	\$67.00
*VITA Special Modeling Liquid	
50ml Order No. V68-752	
250ml Order No. V69-752	\$71.25
*Modeling Liquid Extra	
50ml Order No. V68-753	
250ml Order No. V69-753	\$71.25
*Omega 900 Opaque Liquid	
50ml Order No. V68-906	
250ml Order No. V69-906	\$72.25
VITA Paste Opaque Liquid, 15ml	
Order No. V68-927	\$29.90
30M Modeling Liquid	
50ml Order No. V68-728	
250ml Order No. V69-728	\$76.25

*While Supplies Last

58

Quick-Peg™ Custom Firing Paste

Quick Peg is a paste that is squeezed into/onto the restoration, and placed onto the peg or firing tray. It dries quickly, holds the restoration in place during

firing and is removed easily. It can be used with all-ceramic or porcelain-fused-to-metal restorations. It can also be used for degassing and staining and glazing. It is especially helpful in eliminating the frustrating balancing act experienced placing the crown onto the peg. With over 60 applications per 10cc syringe, Quick-Peg is economical enough to use with every type of restoration.

Quick-Peg, 10cc Syringe 3-Pack	
Order No. A47-320	\$35.50

Firing Trays

Firing Stand A & B Kit	
Order No. A61-201	\$39.00
Single-Crown Firing Stand A (pk of 6)	
Order No. A61-202	\$25.00
Double-Crown Firing Stand B (pk of 6)	
Order No. A61-203	\$27.00
Post Crown Firing Stand C (pk of 6)	
Order No. A61-204	\$27.00
Large Firing Tray G with Pins (pk of 2 with 12 pins)	
Order No. A61-208	\$63.00
Fibrous Firing Pad (set of 2) Pillow	
Order No. A61-211	\$13.00
Firing Tray W (honeycomb)	
Order No. A61-216	\$41.20
Platinum Pins (for firing Tray W)	
Order No. A61-215	\$288.00
Holding Pins for Tray G (pk of 50)	
Order No. A87-145	\$39.00



Instruments

VITA Metal Instruments are made of stainless steel with uniform, easy-to-use handles. The blades and tips are actually soldered into the handles. The full length of each instrument is approximately 14 to 16 cm.

Metal Instruments

*Modeling Knife for Ceramics, E22 For modeling interproximal spaces	
Order No. A51-152	\$ FPO
*Modeling Conical Instrument, E-25.	
For creating finer fissures on masticatory surfaces	
Order No. A51-155	\$63.00
*Scraping Instrument, E-26.	
For vibrating, pressing and scraping	
Order No. A51-154	\$84.00
Ceramic Crown Forceps.	
For risk free handling of all-ceramic restorations	
Order No. A51-250	\$115.00

Additional Instruments

*Glass Mixing Rod	
Order No. A51-251	\$12.00
*Straight Point Foil Tweezers	
Order No. A51-253	\$43.00
*Platinum Foil Scissors	
Order No. A51-255	\$63.00
Furnace Forceps	
Order No. A51-190	\$61.00
*Curved Artery Forceps	
Order No. A51-256	\$85.00
Glass Spatula	
Order No. A51-261	\$12.00
In-Ceram Glass Spatula	
Order No. A51-263	\$15.00
VITA Sponge	
Order No. A51-290	N/A
IMPAC Forceps	
Order No. A51-195	\$5.00

In-Ceram Brushes

In-Ceram Brush IC4	
Order No. A51-211	\$20.00
In-Ceram Brush IC3	
Order No. A51-213	\$14.00

Vident Kolinsky Ceramist Brushes

These high quality build-up and stain brushes are Kolinsky sable that keeps their shape with proper use. Each of the brush sets (2 and 5 brushes) come in a magnetic case.

Order No.	Description	
A51-283	Set of 5 (No. 0, 2, 4, 6, 8)	\$90.00
A51-285	Vident Brush, No. O	\$13.00
A51-286	Vident Brush, No. 2	\$16.00
A51-287	Vident Brush, No. 4 Flat	\$18.00
A51-288	Vident Brush, No. 6	\$21.00
A51-289	Vident Brush, No. 8	\$29.00
Smile Line Brush Tips		

Order No.	Description	
A51-275	Art & Experience Tip #8, pk of 2	\$105.00
A51-276	Art & Experience Tip #6, pk of 2	\$82.00
A51-277	Art & Experience Tip #4, pk of 3	\$35.00
A51-278	Art & Experience Tip #1, pk of 5	\$40.00
A51-279	Art & Experience Tip #000, pk of 5	\$35.00
A51-610	Art & Experience Tip #000, pk of 5	\$35.00

Flexible Diamond Discs

Order No

Vident's Flexible Diamond Discs have diamonds embedded onto a flexible base which allows for exceptional grinding and smoothing, ensuring durability and cutting efficiency. Can be used either wet or dry.

Urder No.	Description	
B25-002	Flexible Disc Assortment Kit, 1 dz	\$141.00
B25-301	Flexible Disc—Extra Fine/Wht., 1 dz	\$141.00
B25-302	Flexible Disc—Fine/Yellow 1 dz	\$141.00
B25-303	Flexible Disc—Medium/Red, 1 dz	\$141.00
B25-304	Flexible Disc—Coarse/Black, 1 oz	\$141.00

Porcelain Polishing Wheels

Specially designed and formulated to achieve a smooth, uniform surface on the restoration. Excellent for polishing delicate margins on porcelain veneers, inlays/onlays and crowns.

- Do not leave excessive residues or exude offensive odors
- Reduce the risks of porcelain cracking
- Generate minimal heat

Order No.	Description	
B35-675	Round Pre-Polisher Gray, 24	\$33.00
B35-676	Round Pre-Polisher Gray, 100	\$110.00
B35-677	Knife Edge Pre-Polisher Gray, 24	\$33.00
B35-678	Knife Edge Pre-Polisher Gray, 100	\$115.00
B35-685	Round Final Polisher, Pink, 24	\$33.00
B35-686	Round Final Polisher, Pink, 100	\$111.00
B35-687	Knife Edge Final Polisher, Pink, 24	\$33.00
B35-688	Knife Edge Final Polisher, Pink, 100	\$115.00

Super Gray Polishing Wheels

Designed to polish porcelains to a glazed finish or metals to a mirror finish. These wheels are used as a final step after grinding to avoid repeating the glaze firing and risking damage to finished cases.

Order No	Description	
B35-665	Round, Box of 24	\$32.00
B35-666	Round, Box of 100	\$115.00
B35-663	Knife Edge, Box of 24	\$32.00
B35-664	Knife Edge, Box of 100	\$115.00

*While supplies last

VITA Akzent Palette, Plastic

The glazed porcelain VITA Ceramic Stain Palette contains 21 mixing wells for mixing and storing dental ceramic stains plus a tight-sealing lid for dust-free storage. (Palette and lid only)

Order No. A51-240......\$56.00

VITA Glass Mixing Slab

The VITA Glass Mixing Slab is 7 1/4" x 5 1/4" and approximately 1/4" thick. It features an etched VITA logo and provides a simple and effective surface for mixing VITA porcelain.

Order No. A51-257\$30.00

VITA Akzent Glaze Spray

Easy to apply, spray ceramic powder for simple, time-saving glazing of all ceramic restorations with a sintering temperature of 850°C.

Akzent Glaze Spray 125ml, 15g	
Order No. S28-790	\$122.00
Patterson Order No. 706-9628	



"Motivation" by Claude Sieber

"Motivation" by Claude Sieber is a stunning pictorial tour of teeth as designed by nature and Claude. Telling the two apart is one of the challenges presented to the reader. With incredible photography, special textured papers, many foldouts and precise reproduction of dozens of restorations, you are sure to be motivated to take your own restorations to a higher level. This oversize hardcover book features a slipcase and a special preservation binding to permit limitless viewing. The photography is museum quality. A softcover textbook is included in a rear pocket that provides details about the photos and techniques used. "Motivation" is a must for anyone interested in the esthetics of modern dental restorative materials.

"Motivation" by Claude Sieber
Order No. L-1244E.....\$180.00

"The Art of Passion" by Ed McLaren

"The Art of Passion: a photographic journey" showcases Dr. McLaren's years of study in the dental profession and photography. The level of understanding and respect for the optical behavior of teeth shows through in stunning spreads, each lovingly reproduced with breathtaking four-color photography. Dr. McLaren takes his book one step further by providing a step-by-step guide to taking and communicating shade information and covers his innovative ceramic layering "Skeleton Buildup Technique" in great detail. This elegant coffee-table sized book is a must-have for everyone who appreciates fine dentistry and fine photography.





Terms and Conditions

The following terms and conditions covering Vident products, equipment, and price lists are subject to change by Vident without prior notice.

Orders Orders may be placed by mail, telephone, facsimile, e-mail or directly with a Vident salesperson. The buyer should include name, address (with zip or postal code), and Vident account number with any order. To facilitate accurate order fulfillment, the buyer should be as specific as possible and include the Vident product number, the desired quantity, and any descriptive information such as size and color shade, wherever indicated. There is no minimum order. For orders under \$30 (\$35 Canadian), a processing fee of \$5 (\$6.00 Canadian) will be applied. Some items shown in Vident's catalog may be available by special order only and may be shipped separately from your original order. In all cases, efforts will be made to consolidate shipments. Vident reserves the right to reject any order in its sole discretion.

Price and Terms of Payment Prices for products are those in effect at the time of purchase. For current retail prices, call Vident. Each shipment will be accompanied by an itemized packing slip. An invoice will either be included with the shipment, or will be delivered separately. Except in those cases where C.O.D. payment is specified or a credit card is used, buyer agrees to pay the full amount due on each invoice within thirty (30) days of the invoice date. All or any portion of an invoice amount which has not been paid within thirty (30) days of the invoice date, will accrue service charges at the rate of one and three-quarters percent (1 3/4%) per month, but in no event in excess of the maximum legal rate, commencing thirty (30) days after the invoice date. Returned checks are subject to a \$25 (\$29 Canadian) charge. Payment by Visa, Master Card or American Express is available.

Credit Approval All orders, shipments and deliveries at all times shall be subject to the approval of buyer's credit by Vident's Finance Department. Vident reserves the right to decline to make deliveries whenever, for any reason, there is doubt as to buyer's financial responsibility, and Vident shall not, in such event, be liable for nonperformance of contract in whole or in part.

Leasing Contact Vident for leasing information.

Taxes The buyer will be responsible for all applicable sales, excise and GST taxes which will be added to the sale price on each invoice and are due and payable on the same terms as the sales price. For tax exemption, the buyer must provide suitable documentation of tax exempt status.

Shipments Usually, products will be sent via ground carrier unless otherwise specified by Vident. Equipment will usually be sent via air freight. The buyer will be responsible for the shipping, handling and insurance costs, which will be prepaid by Vident and added to the sale price on each invoice. These charges are due and payable on the same terms as the sales price. COD shipments will be made at Vident's discretion and will be subject to the carrier's COD surcharges. Shortages or shipment errors must be reported within seven (7) business days of invoice date to be eligible for credit or adjustment as appropriate.

Title and Risk of Loss Title of the product and/or equipment passes to the buyer at the time of shipment. Buyer assumes all risk of damage to or loss or destruction of the product and/or equipment shipped, and no loss, injury or destruction of the product and/or equipment shall release buyer from buyer's obligations to Vident.

Return Policy No returns will be accepted without prior authorization from Vident. A buyer first must obtain a return authorization number from the Vident Customer Service Department. The return authorization number must be included with the return.

Standard catalog products Except as indicated below, Standard catalog products may be returned within thirty (30) days of the date of shipment; provided, however, that return will only be granted if all the following conditions are met:

- **A**. The buyer obtains a return authorization number from the Vident Customer Service Department.
- **B**. The product and/or equipment is unused, in resalable condition, in original packaging, unmarked, with seals and labels intact.
- **C**. A copy of the invoice of sale, a written reason for the return and a copy of the return authorization accompanies the product and/or equipment.
- **D.** The product and/or equipment is returned with shipping costs and insurance prepaid, at the buver's expense to:

U. S. Customers

Vident, 3150 E. Birch Street, Brea, CA 92821

Canadian Customers

Vident Canada c/o TCD, 3505 Laird Road, Unit 7 Mississauga, ON L5L 5Y7

E. The return is made within thirty (30) days or sixty (60) days of the date of the original shipment (see below).

If the above conditions are met and if in Vident's reasonable judgment the return is due to a defective product and/or equipment or Vident's error, the product or equipment will be exchanged or a full credit or refund will be issued, at Vident's discretion. For any returns received between thirty-one and sixty (60) days after the date of shipment, a restocking charge equal to fifteen percent (15%) of the value of the return, or \$10.00 (\$12 Canadian) for returns of less than \$70.00 (\$80 Canadian) will be deducted from the credit or refund, provided that a credit or refund will only be granted if the conditions specified in A, B, C, D and E above are met. No returns will be accepted later than sixty (60) days from the shipment date.

Special Returns Refunds or credits for returns of high noble and noble alloys, will be based upon the market price on the date Vident receives the return, but in no event will the credit or refund exceed the original purchase price. These materials must be unused and in original sealed packages, with seals intact.

Denture Teeth Contact your VITA Denture Tooth DealerReturn of Full Sets Full sets may be returned for exchange for new teeth of like kind only. Full set exchanges will be on a card for card basis. Only currently available teeth in Vident's product line may be exchanged.

Full sets may be returned for exchange within thirty (30) days of the original order date, if all the following conditions are met:

- 1.The buyer obtains a return authorization number from the Vident Customer Service Department prior to return
- 2. The returned denture teeth are unused, in resalable condition.

Orders: 800-828-3839 USA | 800-263-4778 Canada | www.vident.com

- 3.A copy of the return authorization accompanies the denture teeth.
- 4.The product is returned with shipping costs and insurance prepaid, at the buyers expense, to: Vident, 3150 E. Birch Street, Brea, CA 92821 (U.S. Customers only)

For any returns received between thirty-one (31) and sixty (60) days after the date of shipment, a restocking charge equal to fifteen percent (15%) of the value of the return, or \$10.00 for returns of less than \$70.00 will be deducted from the exchange, provided that an exchange or refund will only be granted if the conditions specified in 1, 2, 3 and 4 above are met. No returns will be accepted later than sixty (60) days from the shipment date.

Return of Broken Sets and Full Sets Not In Saleable Condition Broken and full sets may be returned for exchange for new teeth of like kind only. The exchange amount for broken sets will be determined by multiplying the number of teeth remaining on each card by the exchange credit of 90% per tooth based on the effective current published price.

The following items are not returnable for exchange:

- Broken sets or full sets of teeth not in resalable condition
 - o Dirty and/or dusty teeth or wax cards
 - o Teeth melted in the wax on the wax cards
 - o Cracked or scratched wax cards
 - Stickers, non Vita/Vident bar code labels or other markings on the wax cards
- Loose teeth or teeth not on their original card
- Teeth that exhibit grinding marks of any kind
- Teeth returned in an unsanitary manner
- Teeth with permanent markings, i.e. ink pen, pencil, razor edge
- Discontinued moulds or shades

No Return of Certain Products Special order products and custom products (including custom packaging), equipment, products requiring refrigeration, products within 9 months of their expiration date, opened chemical products, shade guides in packaging with broken seals, video tapes, DVDs, CDs and books may not be returned.

Changes Because products and equipment are continually undergoing refinement in design and manufacturing methods, Vident reserves the right to improve, modify or discontinue products or equipment or change pricing at any time without incurring any obligation and without prior notice.

Errors Vident reserves the right to correct clerical or stenographic errors or omissions in its terms, conditions, price lists, or invoices.

LIMITED WARRANTY AND LIMITATION OF REMEDIES: VIDENT WARRANTS ALL PRODUCTS AND EQUIPMENT TO BE FREE FROM DEFECTS IN MATERIAL AND WORKMANSHIP UNDER NORMAL USE AND MAINTENANCE. THIS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, WHETHER EXPRESS OR IMPLIED, AND THERE IS NO WARRANTY OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. THE BUYER'S EXCLUSIVE REMEDY AND VIDENT'S SOLE OBLIGATION UNDER THIS WARRANTY IS LIMITED TO ONE OF THE FOLLOWING, AT VIDENT'S SOLE DISCRETION: REPLACEMENT OR REPAIR OF THE PRODUCT OR A REFUND OF THE PURCHASE PRICE. VIDENT SHALL NOT BE LIABLE FOR INCIDENTAL, CONSEQUENTIAL, SPECIAL OR PUNITIVE DAMAGES. THE BUYER'S CAUSE OF ACTION FOR BREACH OF ANY WARRANTY WILL BE BARRED UNLESS VIDENT RECEIVES WRITTEN NOTICE AT THE ADDRESS SET FORTH BELOW WITHIN THE APPLICABLE WARRANTY PERIOD.

To make a warranty claim, written notice of an alleged defect or other breach of warranty should be sent to Vident at 3150 E. Birch Street, Brea, California 92821 within the applicable warranty period.

Vident Education and Training

Course Registration and Cancellation Policy

Course Registration Course tuition fees are due at the time of registration. Payment may be made by credit card or check. Customers with open accounts in good standing at Vident may charge tuition fees to their accounts.

Course Cancellations Vident reserves the right to cancel courses at any time and will provide as much advance notice as is reasonably possible. Vident's obligation in such instances shall be limited to credit of tuition paid for a future Vident course or refund of paid tuition. Course registrations may be cancelled by the registrant at any time. Paid course tuition will be refunded or credited to a future Vident course based on the following schedule:

Days Prior to Course % Refund/Credit 14 or more 100% 7-13 50% 0-6 0%

Attorney's Fees In the event of a dispute between Vident and the buyer arising out of or relating to an order, Vident shall be entitled to recover from the buyer all reasonable fees, costs and expenses, including, without limitation, reasonable attorneys' fees and expenses.

Force Majeure Should Vident delay or fail to effect delivery or shipment by reason of causes beyond its control, at Vident's option the time of shipment may be extended, or a part or all of an order may be canceled without liability. Such causes shall include but not be limited to acts of God, acts of government, wars, revolutions, strikes, floods, droughts, fires, storms, snow, perils of the sea, unavoidable interruption of transportation, and unavailability of essential materials. In such cases Vident shall notify buyer in writing of the occurrence of such contingency as soon as possible.







© 2011 Vident, Inc. All rights reserved.

Tradenames shown are registered trademarks of VITA and/or Vident unless otherwise indicated.



64

D-Master Shade Guide	5
D-Master Reorder Items	
Accu-Liquid	
Accu-Wax	
Accu-Scan™ Universal Scan Stone	
Accu-Scan Lab Putty	
kkzent Glaze Spray	
kzent Stains	
lleachedguide 3D-Master	
looks	
rushes	
AD-Temp™	
AD-Temp™ DISCS for KaVo Everest	
AD-Waxx	
eramic Alloys	
EREC Etch	
EREC Powder	
Classical Shade Guide	
oloring Liquid	
Diamond Discs	
Diamond Restoration Polish	
lie Sealer Pen	
asyglaze Furnace	
asymatch™ Software	
asyshade Compactasyshade Quality Manager Software	
iring Trays	
urnaces	
lot-Stik Waxer	
MPAC Duplicating Material n-Ceram	
1-Ceram AL	
n-Ceram for inLab	
n-Ceram YZ DISCS for Open CAD/CAM	
n-Ceram YZ	
n-Ceram YZ Pre-Shaded	
n-Ceramat III T Touchscreen Furnace	
nstruments	
nterno Stains	
lic-Lok™ Model Tray System	
aservest Soldering Investment	
inearquide 3D-Master	
ingoform Posterior Denture Teeth	
/IFT Denture Teeth (Multifunctional teeth)	
Modeling Liquids	
Nodisol Pen	80
Imega 900 Porcelain	18
Imega 900 Reorder Materials	
hysiodens Denture Teeth	
M9 Pressable Material	
olishing Sets	
oly-Die™ and Poly-Base™ Material	
orcelain Polishing Wheels	
orcelain Scoops	
luick Fix Mandrel System	79
luick-Peg™ Custom Firing Paste	
hadeAssist Software	
hading Paste	44
imulate Color Indicator	6
imulate Preparation Material	63
llider™ Disposable Articulator	80

Solders, Fluxes and Base Metal Alloys	77
Sta-Put™ Adhesive for Models	75
Super Steam Cleaner	74
Titanium Porcelain	22
VACUMAT 6000M Ceramic Furnace	72
VITABLOCS 3D-Master Guide	5
VITABLOCS Mark II	52
VITABLOCS New User Assortment Kit	57
VITABLOCS TriLuxe forte	53
VITABLOCS TriLuxe	53
VITABLOCS for KaVo Everest	56
VITACOLL™ Bond	68
VITAFOL H Paste	68
Vitapan Denture Teeth	66
Vitasonic Mixer	42
VM7 Porcelain	37
VM9 Porcelain	30
VM9 Esthetic Kit for VITABLOCS	60
VM13 Porcelain	11
VMCC Cold Cure Resin	48
VMLC Light Cure Resin	46
VMKMaster Porcelain	24
Wax-Up Assistant	75
YZ Coloring Liquid	59
YZ Finishing Kit	57
ZYrcomat Furnace	74

Orders: 800-828-3839 USA | 800-263-4778 Canada | www.vident.com