

"CYTO" PEN

FEATURES

- Motorized Micro Needle Therapy
- 12 Needle Tip using 30 gauge Acupuncture Needles to reduce discomfort
- Digital Motor in Hand-piece
- Highest quality for medical professionals (Aluminum Housing)
- Optimal absorption of products
- Adjustable Needle length from .25mm to 2.00mm
- Adjustable 5 level speed control from 7000rpm to 8100rpm
- Compact and ergonomic design

BENEFITS

- Noise free
- Scar-Less Healing
- Little to No Downtime
- Minimal discomfort design
- Self contained no leak cartridges (No cross contamination)
- Safety needle



"The Stem Facial"

Maximizing Aesthetic
Improvements using a
Minimally Invasive
Therapeutic Modality
Featuring The Patented
CYTO Pen®



1.866.233.3412 | www.EmageMedical.com


©2013 EMAGE, LLC. ALL RIGHTS RESERVED. CYTO PEN IS A TRADEMARK OF EMAGE MEDICAL, LLC. SPECIFICATIONS, APPEARANCE AND TECHNICAL FEATURES SUBJECT TO CHANGE WITHOUT NOTICE.

EMAGE

"CYTO" Pen

Medical micro-needling is rapidly gaining popularity as an effective treatment option for reducing wrinkles, scars, stretch marks, hyperpigmentation, and improving overall skin tone and texture. Compared to chemical peels, dermabrasion and laser treatments, microneedling causes substantially less tissue damage with minimal down time. STEM Facial™ consists of creating scores of microscopic porations using the CYTO Pen® with patented disposable multi-needle attachments. The length of needle used is determined


according to the area of skin being treated and the purpose of treatment. Shorter needles are used to enhance penetration of the STEM Facial serum through the stratum corneum, the outermost "barrier" layer of the skin.



Use of longer needles for micro-needling creates hundreds of small injuries in the papillary dermis. The subsequent healing response produces additional collagen and elastin which improves the appearance of wrinkles, scars, and stretchmarks.

EMAGE Medical, a recognized leader in providing advanced technology to professional esthetic practices, now offers a new patented CYTO Pen® that provides unsurpassed control of treatment depth and speed. The CYTO Pen® uses disposable sterile needling cartridges and provides easy "click to adjust" control of both needle penetration and oscillation speed. Adjustment can be made

rapidly and precisely so that each area of skin is appropriately treated. The electrical and mechanical components of the device are sealed from possible contamination by blood or topically applied



products assuring patients are protected from cross contamination. Removal and replacement of disposable needling cartridges is fast with a secure lock feature that prevents inadvertant dislodgement during treatment.

The STEM Facial™

The STEM Facial™ combines medical micro-needling with AnteAGE Professional, a high-science topical containing stem cell derived serums that modulate inflammation, promotes healing and removes pigment and wrinkles. Over the course of six weeks, patients will receive three micro-needling treatments with a take-home two-part system.

Developed by physician scientists, the major active is produced through laboratory culture of the specific human stem cells proved to be responsible for command and control of healing throughout the body. In AnteAGE Professional, the concentration of the important biosignals produced by the cells is maximized. The stem cell biosignals are only the beginning of what sets AnteAGE

apart from other anti-aging skincare products. The sixteen actives in STEM Facial™ include peptides, antioxidants, anti-inflammatories, vitamins, and more, each scientifically engineered to work in conjunction with CYTO Pen® micro-needleling.

EMAGE Medical partnered with AnteAGE Professional by Cellese as part of the STEM Facial™ for a simple reason: documented cosmetic improvement.

