


Phonak Lyric™

Clear natural sound

Lyric advantages:

- Being positioned 4mm from the tympanic membrane (TM) **increases access to spectral speech cues**
- Sound is collected and directed to the TM via the **natural path to the ear**
- **Short processing time** thanks to analog technology



Studies show:



No significant difference between Lyric and the open ear KEMAR for normal hearing subjects¹



Less errors locating the sound source **preserves normal localization accuracy**²



Average percentage **localization scores are higher** than for a daily wear hearing aid³

¹ Lough, M., Whiston, H., Dillon, H., & Stone, M. (2020). Subjective benefits of a deeply inserted extended wear hearing aid compared to other in-canal sound-delivering hearing aids. Manuscript in preparation.

² Brungart D et al, (2016). Military applications of Extended wear hearing technology.

³ Standaert, L. (2017). Benchmark comparison (unpublished). Contact claims@phonak.com if interested in more information.