



**Elite Mobility**  
**Head Office**  
**19 Broad Street**  
**Staple Hill**  
**Bristol**  
**BS16 5LN**

**Tel: 0800 1695910**  
**info@elitemobility.co.uk**  
**www.elitemobility.co.uk**



**Sarum**



Curved stairlift

**Stannah**



I didn't need to worry about the bend in my staircase – my Sarum handles it perfectly.



## Make a step change

### If a Stannah stairlift could make life easier, why wait?

Life is better when you have the energy to do the things you want to do, like enjoying a visit from the grandchildren or entertaining friends. Our experience has taught us that the sooner a stairlift is installed, the sooner you will be able to start enjoying the things you love doing. Just think how a Stannah stairlift could make your day-to-day life easier.

#### A stairlift that's right for you and your home.

We know that not all homes are the same. If you think your stairs are unsuitable for a stairlift, we think you might be pleasantly surprised.

- Our stairlifts can accommodate curved staircases of different sizes and layouts.
- We can install our rails on either side of the staircase.
- Our slimline rails and fold-away chairs make the best use of available space, leaving enough room for others to use the staircase.
- All our rails are tailor-made, so they fit your staircase perfectly.
- Contrary to popular belief our rails are fitted to the staircase and not the wall – ensuring minimal mess and bother.
- Installation is quick and easy.
- Our stairlifts can navigate all kinds of corners and landings with a range of rail solutions.

**“I am very satisfied with the stairlift; it is easy to use, good to look at and is making life much easier for me!”**

**Mrs Collyer, 2011**



Outside rail



Inside rail



## A company you can trust

We are a British, family-run company and have been for five generations. The current managing directors are the great, great grandchildren of Joseph Stannah, who started the company in the 1860s. His engineering principles and strong family values run as true today as they always have.

Now, 140 years later, we have manufactured more than 500,000 stairlifts and installed products in over 40 countries. The half-million milestone was commemorated with a visit from HRH The Prince of Wales – a proud day for all those who work at Stannah.

## The Sarum

The Sarum is one of our most popular stairlifts. It is perfect for navigating landings, curves and even spiral staircases smoothly, whilst being comfortable and easy-to-use.

Its design means that it won't obstruct your hallway or stairs, and you can pick from a choice of upholstery colours so it suits the interior of your home.

Over 500,000  
stairlifts sold  
globally

## Key features

### 1. Footrest lever

To store the stairlift away, you can fold away the arms and seat and press the lever to raise the footrest. There's no awkward bending down, which is particularly useful at the top of the stairs.

### 2. Swivel seat

Using the levers enables you to swivel the seat easily, making it simpler to get off the chair. Alternatively, you can choose to have a powered swivel seat, which turns the chair to face the landing automatically.

### 3. Retractable seatbelt

Fasten the seatbelt for added security.

### 4. Safety edges

In the event that something obstructs the stairlift, sensors on the edge of the carriage and footrest will automatically bring the stairlift to a safe stop.



### 5. Directional controls

Right-handed or left-handed, we can put the controls on either arm to suit you.



### 6. Wall controls

Ideal for when you need to summon your stairlift from another floor, or perhaps send it up or down to someone else who wants to use it.



**“It would be remiss of me if I did not write thanking the workmen for their kindness, cleanliness and expertise in their work.”**

Mr Raleigh, 2011

## Using your Sarum

The simple-to-use Sarum stairlift has been designed to make your life easier and will help you to rediscover your independence. This step-by-step guide demonstrates some of its key features.



1. The wall-mounted controls at the top and bottom of the stairs make it easy to call the stairlift wherever you are.



2. Simply unfold the seat and the arms and use the lever to lower the footrest.



3. Once unfolded, make yourself comfortable in the seat. Fasten your seatbelt for extra security.



4. Gently apply pressure to the control in the direction you wish to travel and your Sarum will glide smoothly up or down the stairs.



5. When at the top of the stairs, the footrest will stop level with your landing. You can then swivel the chair to face the landing to alight safely.



6. For added safety, the Sarum can remain in the swivelled position to form a barrier at the top of the stairs.



7. When not in use, your Sarum folds neatly against the wall, giving you and others plenty of space.

## Your Sarum - tailored for you

### Standard Package



Twin-rail pearl colour

In addition to the standard package, you can choose from a range of other optional extras that will ensure your Sarum is just right for you. Our sales consultant will guide you through the various options available to you.



Retractable seatbelt



Manual swivel seat



Directional controls – push button

Left hand  Right hand



Footrest with manual lever

## Upholstery colours and optional extras

### Vinyl upholstery



Clava



Forest



Skye



Mallum



Roeburn

### Curved rail colours



Pearl (standard)



Nut Brown



Jet Black



Window Grey



Beige



Ochre Brown

### Rail start



Inclined



Droop



Hinged



Run-Off

### Rail finish



Overrun



Inclined

### Swivel



Powered

### Controls



Joystick

### Seatbelt



Diagonal lap harness



Five-point harness

### Seat



Reduced seat  
(heightened seat option also available)

### Dimensions

Min. width of unfolded stairlift:  
679mm / 26.7"  
Min. width of folded stairlift:  
405mm / 16"  
Max. weight:  
120kg / 19 stone