

Bronchitrac® "L"

SELECTIVE BRONCHIAL SUCTIONING SET



Left Main Stem Bronchial Tube Suctioning

Supplies for respiratory care departments are in **critical need**. Suctioning of secretions from the left endobronchial tree is often difficult in intubated patients.¹

The single-use, sterile **Bronchitrac "L" Selective Bronchial Suctioning Set** is designed to assist in placement of the catheter **directly into the left main stem bronchial tube**, to remove lung secretions often caused by pneumonia and/or acute respiratory distress syndrome (ARDS).



- ✓ **No need to twist the catheter or patient's body**
- ✓ **Can be used for deep-lobe suctioning**
- ✓ **Patented curvature of catheter ensures consistent entry**
- ✓ **Sterile, single-use set**
- ✓ **Latex-free**

41-7814

14 Fr Bronchial Suction Catheter,
Slide-cover Connector and
(2) Vinyl Gloves, Sterile (20/Boxes)

To order, please call Heartland Medical Specialties, Inc.
(P) 913-498-1318 (F) 913-541-1897
or email: orders@heartlandmedspec.com

1. Panacek EA, Selective left endobronchial suctioning in the intubated patient. Chest. 1989 Apr;95(4):885-7