


COMFIT™ Endotracheal Tube Holder

Significantly Decreases the Incidence of Self-Extubation*

- No tape on the ET tube; no tape on the patient
- May eliminate need for bite block
- Completely adjustable
- Wide cotton-lined neckband minimizes skin irritation, providing maximum patient comfort
- Minimal plastic loop around the ET tube allows access to the oral cavity
- Economical in two ways: low initial cost, no frequent changing
- Latex-free product
- Significantly decreases incidence of Self-Extubation*

Complies with ACLS Guidelines

ACLS Guidelines recommend the use of commercially available endotracheal tube holders as a substitute for tape.



***The tapeless way to secure
an endotracheal tube***

Ordering Information

COMFIT

Product No. 42-1820

Plastic double securing loop, 1½" x 19" neckband.

Packaging

Non-sterile, 25 per box, 100 (4 boxes) per case.

* Reduction of Self-Extubation Using Comfit Tracheal Tube Holder
J. McAndrew RRT, Respiratory Care, Nov. 94, Vol. 39 No. 11.

ET Tube Fixation Made Simple